
**SoloPath® Re-Collapsible Access System
INSTRUCTIONS FOR USE**

Description

Onset Medical Corporation's SoloPath® Re-Collapsible Access System consists of a flexible, reinforced polymer sheath with an external collapsible outer jacket and specially folded distal end (the Sheath) pre-mounted over a central balloon dilatation catheter (the Expander). The folded distal region of the Sheath is small in diameter, thus facilitating passage through the vessel. The SoloPath Assembly is inserted percutaneously into the femoral artery, over a guidewire, with the deflated Expander in place. Once positioned into the vessel, the Expander balloon, when inflated with liquid, exerts controlled radial force, enlarging the folded distal region of the Sheath and surrounding anatomy. The Expander balloon is deflated and the Expander is removed leaving a large central lumen extending from the proximal end to the distal end of the Sheath. The Sheath is designed as a guide for catheters and/or devices introduced into the femoral artery. Prior to removal, the outer jacket is activated with liquid under low pressure, collapsing the outer sheath diameter for ease of removal.

Figure 1 – SoloPath® Re-Collapsible With and Without Expander

Product Specifications

The SoloPath® Re-Collapsible Access System is available with three (3) internal diameters and two (2) usable lengths. Refer to Table 1.

Table 1 – SoloPath® Re-Collapsible Access System

Model Number	Expanded Sheath I.D.		Expanded Sheath O.D.		Expandable Length	Distal Tip Outside Diameter		Usable Length	Maximum Rated Inflation Pressure		Max. Collapsible Jacket Pressure	
	French	mm	French	mm		French	mm		cm	atm/bar	kPa	atm/bar
SR-1925	19	6.33	23	7.67	20	5.3	1.77	25	20	2026	6	608
SR-1935	19	6.33	23	7.67	30	5.3	1.77	35	20	2026	6	608
SR-2225	22	7.33	26	8.67	20	5.3	1.77	25	20	2026	6	608
SR-2235	22	7.33	26	8.67	30	5.3	1.77	35	20	2026	6	608
SR-2425	24	8.00	28.5	9.50	20	5.3	1.77	25	20	2026	6	608
SR-2435	24	8.00	28.5	9.50	30	5.3	1.77	35	20	2026	6	608

Warning: The SoloPath® Re-Collapsible Access System is supplied sterile and is intended for single use only. Do not reuse, resterilize or reprocess this device. Do not use this device if product or sterile packaging is damaged.

Device Construction

- Recommended for travel over a 0.038" (0.97mm) or smaller guidewire.
- A flexible Sheath composed of reinforced polymer is compressed, folded and wrapped around a central, removable balloon dilatation catheter (the Expander).
- The distal fairing tip and entire sheath usable length is hydrophilically coated.
- One radiopaque marker is located on the Sheath to indicate the end of the expandable portion of the Sheath.
- A collapsible outer jacket on the Sheath exterior with collapsation port located on the proximal hub assembly.
- The Expander balloon is non-compliant and size limited. The Expander has both an inflation lumen and a guidewire lumen.
- The Expander is configured to be used with standard inflation devices having an appropriate pressure rating, a male Luer lock coupler, and a capacity of 25-CC.
- The Sheath is terminated, at its proximal end, by a hub having a straight through hemostasis port capable of accepting catheters equal to the full sheath inner diameter.
- A sheath flush line, terminated with a 3-way stopcock, permits flushing of heparinized saline into the central port of the Sheath hub. The Expander incorporates a T connector, the through lumen port being terminated with a hemostasis valve for guidewire introduction, while the right angle port is for inflation of the Expander balloon. The stopcocks are used for priming the device prior to insertion. The flush lines are not intended for high pressure injection.

Indications for Use

The SoloPath® Re-Collapsible Access System is intended to be inserted percutaneously into the femoral artery, over a guidewire and once expanded, to provide a guide for catheters and/or devices introduced into the femoral artery.

Caution: A thorough understanding of the technical principles, clinical applications and risks associated with vascular access is required before attempting to use this device.

Contraindications

There are no known contraindications for the SoloPath® Re-Collapsible Access System.

Precautions

- Inspect all components before use.
- Carefully read the instructions for use (IFU) before using this device to help reduce the potential dangers associated with the access technique such as air emboli and/or perforation of the vasculature.
- Prior to inserting the device into the patient, insert the guidewire through the balloon Expander central lumen, checking and correcting for excessive resistance as the tip of the guidewire advances through the Sheath/Balloon Expander assembly.
- Flush, with heparinized saline, the Sheath, Expander, and any catheters to minimize the potential for air emboli and/or clot formation.
- Do not remove the balloon Expander or catheter from the Sheath rapidly or damage to the hemostasis valves may occur, causing blood leakage out through the valve.
- If resistance is met when advancing or withdrawing the guidewire or Sheath, determine the cause by fluoroscopy and correct before continuing with the procedure.

Potential Complications

Complications associated with the use of the SoloPath® Re-Collapsible Access System include but are not limited to, vessel perforation, vessel avulsion, aortic perforation, infection, post-procedural discomfort, air embolism, tissue trauma, hemorrhage, thromboembolism, thrombosis, injury to the vascular introduction site, and death.

Warnings

- Do not alter the device in any way.
- Do not reuse this device. Thorough cleaning of biological and foreign material from the device is not possible. Adverse patient reactions may result from reuse of this device.
- Only those physicians who are specially trained in vascular access procedures should use this device.
- Maintain continuous patient blood pressure monitoring throughout the procedure.
- Immediate availability of anterior-posterior and lateral fluoroscopy is required throughout the procedure.
- Always observe acceptable pressure tracings prior to advancing the dilator or any other component.
- Do not permit air to enter the Sheath. Remove components and make catheter exchanges slowly.
- Using the flushing ports, purge all air from the device prior to insertion.
- Fibrin may accumulate in or on the Sheath tip during the procedure.

Suggested Procedure

Preparation

A standard preparation for endovascular access should be completed per hospital protocol. Standard fluoroscopic equipment should be available for use during the procedure. Proper radiological protection should be provided for all attending personnel. It is highly recommended to verify that the size and condition of the vessel is appropriate for the size sheath chosen (reference product specification – Table 1).

Caution: Use of incorrect sheath size may result in a higher risk of vessel damage and/or avulsion.

Device Inspection

Using aseptic technique, remove the *SoloPath*® from its sterile pouch. Visually inspect the device to make sure there is no distortion or kinking in the shaft of the sheath or in the folded distal end and that a smooth taper exists between the distal end of the Sheath and the balloon Expander. The balloon Expander shaft is fixed into position within the folded Sheath body.

Caution: Do NOT attempt to separate the Sheath from the balloon Expander at any time prior to expansion of the Sheath.

Device Preparation

Use aseptic technique for all steps of the procedure. Open the three-way stopcock on the Flush Port of the Sheath and gently flush with sterile, heparinized saline. Do not flush aggressively or create pressure during the flushing process or the sheath may open prematurely. Close the stopcock. Close the stopcocks to prevent saline loss or air embolism during the procedure.

Percutaneous Access

A standard Seldinger preparation and access are completed into the vessel using, for example, an 18-gauge thin wall access hollow needle. A 0.038" (0.97mm) stiff guidewire with a floppy tip is advanced through the access needle and into the vessel under fluoroscopic guidance.

Device Insertion

Ensure that all ports have been flushed and primed with heparinized saline and are free of air. The *SoloPath*® with its integral balloon Expander is advanced as a single unit, under fluoroscopic control, over the 0.038" (0.97mm) guidewire until its distal radiopaque sheath marker is positioned at the desired location within the vasculature.

Caution: The *SoloPath*® is to be advanced over the guidewire as a single unit. Do not allow separation or misalignment of the central balloon Expander from the unexpanded Sheath prior to expansion within the vasculature.

Caution: Should resistance be encountered, cease advancing the *SoloPath*® until the cause of the resistance can be determined and corrected.

Caution: The balloon Expander is intended for inflation only after the unexpanded Sheath has been completely advanced to its final target location within the vasculature.

Caution: Do not attempt to refold an expanded sheath.

Sheath Expansion

A dilute 50% solution of radiopaque contrast media and sterile saline is prepared and approximately 20-CC are drawn up into a high-pressure inflation syringe equipped with an appropriate pressure gauge (not supplied). Care should be taken to remove all air from the syringe, the Expander, and associated tubing. Remove cap from the Expander Balloon Inflation Port and discard it. Attach the pressure line of the inflation device to the Expander Balloon Inflation Port.

Under fluoroscopic control, the dilute contrast media is injected to the Maximum Rated Inflation Pressure per Table 1. It is normal to observe a drop in pressure as the Expander balloon progressively expands the folded distal section of the Sheath. The rated inflation pressure should be maintained for a minimum of 60 seconds to expand any "waist" that may remain along the length of the expanded Sheath. During inflation and or at rated pressure, verify the sheath is fully expanded under fluoroscopic visualization.

Warning: Only use liquid to inflate the balloon Expander. NEVER use gas, such as air, to inflate the balloon Expander.

Deflated Expander Removal and Sheath Use

Apply suction, using the high-pressure syringe, to the inflation/deflation port on the balloon Expander hub in order to deflate the balloon Expander. Remove the deflated Expander from the expanded Sheath slowly to minimize blood loss. Introduce the appropriate therapeutic device or catheter through the central working channel of the expanded Sheath.

Caution: Do NOT attempt to separate the balloon Expander from the Sheath while the Expander is under pressure.

Caution: Care should be taken to maintain the position of the expanded Sheath within the vasculature during the procedure. Advancing an already expanded Sheath must be avoided as this movement may not be possible and any attempt will potentially damage the vessel.

Caution: If resistance is encountered during device insertion, stop immediately, remove device from the sheath and reinsert dilator. Repeat steps for sheath expansion in order to optimize sheath expanded diameter. If resistance is still encountered, stop and remove sheath from patient.

Warning: Never attempt to force anything through a partially expanded sheath. Attempting to force a device through a partially expanded sheath may result in sheath damage, breakage, or device hang up.

Sheath Removal

Remove any instrumentation from the sheath slowly through the hemostasis valve, being careful to minimize blood loss. A dilute 50% solution of radiopaque contrast media and sterile saline is prepared and approximately 20cc are drawn up into an inflation syringe equipped with a pressure gauge (not supplied). Care should be taken to remove all air from the syringe. Attach the inflation device to the collapsible port stopcock on the proximal sheath hub. Verify that the stopcock is in the open position between the collapse port and inflation syringe. Pull a slight negative pressure on the syringe to remove any excess air in the system and inject the diluted contrast solution to 6 atmospheres per Table 1. It is normal to observe a drop in pressure as the sheath progressively collapses. Apply negative pressure on the syringe to remove fluid from the jacket and reduce sheath profile. Prior to removing the sheath, apply slight pressure with the syringe by hand to gently fill the jacket and obstruct blood flow during the removal process. As an alternative method for sheath removal, a 20ml syringe filled with 5ml of 50% contrast solution may be attached to the collapsible port stopcock to apply slight pressure with the syringe by hand to gently fill the jacket and obstruct blood flow during the removal process. While maintaining hemostasis with slight manual compression at the puncture site, gently withdraw and remove the sheath from the vessel.

Warning: Only use liquid to collapse the sheath. NEVER use gas, such as air, to collapse the sheath.

Caution: Do not exceed maximum jacket inflation pressure.

How Supplied

Contents of unopened, undamaged package are sterile and non-pyrogenic. The device has no components made of natural rubber latex.

This device is sterilized with ethylene oxide gas. It is intended for single use only. Do not use if the package is opened or damaged.

Storage

Store in a cool, dry place. Rotate inventory so that products are used prior to the sterilization expiration date on package label.

Warranty

Warranty and Limitations: Onset warrants to Customer that its products are free from defect in design, workmanship and materials. For its sterile products Onset further warrants to Customer that such products will remain sterile to the expiration date specified on the applicable product's label, provided that the original packaging remains intact, and is not compromised in any manner. This warranty shall not apply to any products, which have been re-sterilized, repaired, altered or modified in any manner, or to any products, which have been improperly stored, installed, operated, used or maintained in any manner. In the event of any breach of its warranty set forth above, Onset's sole obligation shall be to repair or replace, at its sole option, any product that Onset determines was defective in workmanship or materials at the time of shipment if notice thereof from Customer is received by Onset within the "use before" or expiration date, as applicable, described on such product's label. Except as expressly provided above, Customer assumes all responsibility and all other liability, whether based upon warranty, contract, negligence, or otherwise, for injury or damages resulting from handling, possession, use or misuse of any Onset product and agrees to defend and indemnify Onset for same. Because Onset has no control over the operation, use, or maintenance of its products, the selection of patients or their condition, and does not warrant the performance of its products, THE WARRANTY SET FORTH ABOVE IS EXPRESSLY ACKNOWLEDGED BY CUSTOMER TO BE IN LIEU OF ANY OTHER EXPRESS WARRANTY AND OF ANY OTHER OBLIGATION ON THE PART OF ONSET. The remedies set forth in this Warranty and Limitations shall be the exclusive remedy available to any person. No officer, employee, representative or any other agent of Onset has any authority to change or otherwise modify in any respect any of the foregoing, or to act in any manner to assume or bind Onset to any additional liability or responsibility in connection to its products. Customer's submission of any orders for any Onset product(s) shall be deemed irrevocable acceptance of all of the terms and conditions of the foregoing Warranty and Limitations.

EXCEPT FOR THE EXPRESS WARRANTY DESCRIBED ABOVE AND AS OTHERWISE SET FORTH HEREIN, ONSET GRANTS NO OTHER WARRANTIES OF ANY NATURE, DIRECTLY OR INDIRECTLY, EXPRESS OR IMPLIED, BY STATUTE OR OTHERWISE, REGARDING THE PRODUCTS, THEIR FITNESS OR SUITABILITY FOR ANY PURPOSE, THEIR QUALITY, CONDITION, MERCHANTABILITY, OR OTHERWISE. ONSET'S LIABILITY UNDER THE WARRANTY SHALL BE LIMITED TO REPAIR OR REPLACEMENT OF PRODUCTS IN EXCHANGE FOR PRODUCTS RETURNED TO ONSET. IN NO EVENT SHALL ONSET BE LIABLE FOR ANY LOSS OF PROFITS, SPECIAL, CONSEQUENTIAL OR INCIDENTAL DAMAGES FOR BREACH OF LIABILITY.

U.S. and International Patents pending
© Copyright 2010 Onset Medical Corporation
SoloPath® is a registered trademark of Onset Medical Corporation

CONTENTS

MAXIMUM RATED INFLATION PRESSURE

BATCH CODE

USABLE LENGTH

CONSULT INSTRUCTIONS FOR USE

CATALOGUE NUMBER

EXPANDABLE LENGTH

DO NOT USE IF PACKAGE IS DAMAGED

USE BY DATE

DO NOT REUSE

STERILIZED USING ETHYLENE OXIDE

DO NOT RESTERILIZE

WARNING: Federal law (USA) restricts this device to sale by or on the order of physician


MANUFACTURER

Onset Medical Corporation
13900 Alton Parkway, Suite 120
Irvine, CA 92618 USA
Phone 949.716.1100
Fax 949.716.1101

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SoloPath®, Onset Medical Corporation şirketinin tescilli ticari markasıdır