en	ACTIVE SAFETY NEEDLE	2-3
fr	AIGUILLE DE SÉCURITÉ ACTIVE	4-5
de	AKTIVE SICHERHEITSKANÜLE	6-7
es	AGUJA DE SEGURIDAD ACTIVA	8-9
pt	AGULHA DE SEGURANÇA ATIVA	10-11
it	AGO DI SICUREZZA ATTIVO	12-13
nl	ACTIEVE VEILIGHEIDSNAALD	14-15
sv	AKTIV SÄKERHETSKANYL	16-17
da	AKTIV SIKKERHEDSKANYLE	18-19
no	AKTIV SIKKERHETSKANYLE	20-21
fi	AKTIIVINEN TURVANEULA	22-23
el	ΒΕΛΟΝΑ ΕΝΕΡΓΗΤΙΚΗΣ ΑΣΦΑΛΕΙΑΣ	24-25
ru	АКТИВНАЯ БЕЗОПАСНАЯ ИГЛА	26-27
pl	AKTYWNA IGLA BEZPIECZNA	28-29
hu	AKTÍV BIZTONSÁGI TÚ	30-31
cs	AKTIVNÍ BEZPEČNOSTNÍ JEHLA	32-33
sk	AKTIVNA VARNOSTNA IGLA	34-35
tr	AKTİF GÜVENLİK İĞNESİ	36-37
et	AKTIIVSE KAITSEGA NÕEL	38-39
lv	ADATA AR AKTĪVĀS DROŠĪBAS MEHĀNISMU	40-41
lt	AKTYVIOJI APSAUGINĖ ADATA	42-43
sl	AKTIVNA BEZPEČNOSTNÁ IHLA	44-45
sr	IGLA SA AKTIVNOM ZAŠTITOM	46-47
ro	AC CU SIGURANȚĂ ACTIVĂ	48-49
bg	АКТИВНА ЗАЩИТЕНА ИГЛА	50-51
uk	ГОЛКА З АКТИВНОЮ СИСТЕМОЮ БЕЗПЕКИ	52-53
hr	AKTIVNA SIGURNOSNA IGLA	54-55

SYMBOL EXPLANATION

Manufacturer	Do not re-sterilize	Fragile, handle with care
Use-by date	Non-pyrogenic	Contents
Sterilized using ethylene oxide	Rx only Caution: Federal law (USA) restricts this device to sale by or on the order of a physician	Keep away from sunlight
Single sterile barrier system	Date of manufacture	Keep dry
Consult instructions for use	Catalogue number	Temperature limit
Do not use if package is damaged and consult instructions for use	Batch code	
Do not re-use	Medical device	

PRODUCT DESCRIPTION

Non toxic

Needle Gauge	Needle Length	Needle Outer Diameter	Needle Length	Wall Thickness	Bevel Angle
25G	1"	0.5 mm	25 mm	Thin Wall (TW)	12°
25G	5/8"	0.5 mm	16 mm	Thin Wall (TW)	12°

INTENDED PURPOSE

The K-Pack Embrace™ Active Safety Needle, being a hypodermic needle with safety shield, is intended to inject fluids into, or withdraw fluids from, parts of the body below the surface of the skin. The safety shield shall be manually locked (activated), after use, to cover the needle to minimize the risk of accidental needle stick.

INDICATIONS

The K-Pack Embrace™ Active Safety Needle is for general application - for treatment (injection of fluids) or diagnosis (withdrawal of fluids).

CONTRA-INDICATIONS

No contra-indications.

PATIENT TARGET GROUP

Intended for general application.

INTENDED USERS

Healthcare professional or lay person.

CLINICAL BENEFIT

The K-Pack Embrace™ Active Safety Needle has an indirect clinical benefit (indirect medical effect) since it is used for injection or withdrawal of fluids and since it prevents needle stick injury.

WARNINGS

- Do not use if unit package (= cap and case) is damaged.
- Use immediately after opening the unit package (= cap and case).
- Do not re-cap.
- The needle is made of stainless steel containing nickel and cobalt. Cobalt is classified as CMR* 1B and is present in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
* CMR = Carcinogenic, mutagenic or toxic to reproduction (CLP Regulation EU 1272/2008)

PRECAUTIONS

- If the needle is bent or damaged, no attempt should be made to straighten the needle or use the product.
- After use, dispose of safely as medical waste in a sharps disposal container and/or according to health institution policies. The product is biohazardous and is physically hazardous due to its sharp edge.

REPORT OF INCIDENT

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

PRECAUTIONS FOR STORAGE

- Store between 2°C and 30°C.
- Keep dry.
- Keep away from sunlight.
- Fragile, handle with care.

DIRECTIONS FOR USE


Aseptic technique, proper skin preparation and continued protection of the site are essential. Observe Universal Precautions on ALL patients.

Precaution: Handle with care to avoid needle sticks. Keep hands behind the needle at all times during use and disposal.

COMPONENTS OVERVIEW

1. PREPARING FOR USE – OPENING

Inspect tamper-evident label.


Discard the product if the tamper-evident label is damaged or missing.

Twist and pull to break the tamper-evident label and to remove the cap.

Caution: Avoid touching the hub after removing the cap.


COMPONENTS OVERVIEW

2. PREPARING FOR USE – ATTACHING

Firmly attach the hub onto the (pre-)filled syringe.


Caution: Do not overtighten.

Luer slip: pushing and twisting the needle


Luer lock: twisting the needle


Pull the case straight off. The needle is protected by the safety shield and protector.


COMPONENTS OVERVIEW

3. PREPARING FOR USE – REMOVE PROTECTOR


Move the safety shield away from the needle and toward the syringe barrel to the angle shown.


Pull the protector straight off.

COMPONENTS OVERVIEW


4. INJECTION AND/OR WITHDRAWAL OF FLUIDS

Perform injection and/or withdrawal of fluids according to established technique.


Caution: Do not touch the safety shield during injection and/or withdrawal of fluids.

Note: The needle "bevel up" position is oriented to the position of the safety shield.


5. LOCK (ACTIVATE) SAFETY SHIELD

Use a **one-handed technique** to lock the safety shield using any of the three (3) activation methods illustrated.

Note: Activation is verified by an audible and/or tactile "click".

6. INSPECT SAFETY SHIELD ACTIVATION

Visually confirm that the safety shield is fully engaged.

Caution: Do not attempt to unlock (deactivate) the safety device by forcing the needle out of the safety shield.

7. DISPOSAL

After use, dispose needle and syringe safely as medical waste in a sharps disposal container and/or according to health institution policies.

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