en	NEEDLE	2-3
fr	AIGUILLE	4-5
de	KANÜLE	6-7
es	AGUJA	8-9
pt	AGULHA	10-11
it	AGO	12-13
nl	NAALD	14-15
sv	KANYL	16-17
da	KANYLE	18-19
no	KANYLE	20-21
fi	NEULA	22-23
el	ΒΕΛΟΝΑ	24-25
ru	ИГЛА	26-27
pl	IGŁA	28-29
hu	TŰ	30-31
cs	JEHLA	32-33
sk	IHLA	34-35
tr	İĞNE	36-37
et	NÕEL	38-39
lv	ADATA	40-41
lt	ADATA	42-43
sl	IGLA	44-45
sr	IGLA	46-47
ro	AC	48-49
bg	ИГЛА	50-51
uk	ГОЛКА	52-53
hr	IGLA	54-55

SYMBOL EXPLANATION

LOT	Batch code	Do not re-use	MD	Medical Device
REF	Catalogue number	Do not use if package is damaged and consult instructions for use		Non pyrogenic
Rx only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician	Fragile, handle with care		Single Sterile Barrier System
	Consult Instructions for Use	Keep away from sunlight		Sterilized using ethylene oxide
	Contents	Keep dry	2°C - 30°C	Temperature limit
	Date of Manufacture	Manufacturer		Use-by date
	Do not resterilize			

PRODUCT DESCRIPTION

Non toxic

Needle Gauge	Needle Length	Needle Outer Diameter	Needle Length	Wall Thickness	Bevel Angle
25G	5/8"	0.5 mm	16 mm	Thin Wall (TW)	12°
27G	1/2"	0.4 mm	12 mm	Regular Wall (RW)	12°
27G	1/2"	0.4 mm	12 mm	Thin Wall (TW)	12°

INTENDED PURPOSE

The K-Pack Surshield™ Needle is a hypodermic needle with a passive sharps protection feature that covers the cannula immediately and permanently after use; and is intended for use in combination with a (pre-)filled syringe for subcutaneous and intramuscular injection.

INDICATIONS

The K-Pack Surshield™ Needle is for general application - for treatment (subcutaneous or intramuscular injection).

CONTRA-INDICATIONS

No contra-indications.

PATIENT TARGET GROUP

Intended for general application. Selection of the needle size for intramuscular or subcutaneous application, based on patient's morphology, is at the discretion of the medical practitioner following standard injection procedures.

INTENDED USERS

Healthcare professional or lay person.

CLINICAL BENEFIT

The K-Pack Surshield™ Needle has an indirect clinical benefit (indirect medical effect) since it is used for a subcutaneous or intramuscular injection. For general application. The passive sharps protection feature has an indirect clinical benefit since it prevents needle stick injury.

WARNINGS

- Do not use if unit package is damaged.
- Use immediately after opening the unit package.
- The needle is made of stainless steel containing nickel and cobalt. Cobalt is classified as CMR* 1B and is present in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
- * CMR = Carcinogenic, mutagenic or toxic to reproduction (CLP Regulation EU 1272/2008)

PRECAUTIONS

- Do not use for intravenous injection.
- After use, dispose of safely as medical waste in a sharps disposal container and/or according to health institution policies. The product is biohazardous.

DIRECTIONS FOR USE

COMPONENTS OVERVIEW


1. PREPARING FOR USE – OPENING

Inspect tamper-evident label.

Discard the product if the label is damaged or missing.


Twist and pull to break the tamper-evident label and to remove the cap.

Caution: Avoid touching the hub after removing the cap.

COMPONENTS OVERVIEW


2. PREPARING FOR USE – ATTACHING


Firmly attach the hub onto the (pre-filled syringe by:
a: pushing the needle (luer slip).
b: twisting the needle (luer lock).

Caution: Do not overtighten.

Pull the case straight off.
The needle is protected by the shield.

Caution: Keep your hands away from the shield at all times!


3. PREPARING FOR USE – CHECK

Caution: Do not re-cap the needle to avoid accidental activation of the shielding mechanism.

Caution: Do not touch the shield at any time prior to or during the injection. Any pressure on the shield may cause the shielding mechanism to lock, making the needle unusable.

Check the visual indication.

Blue = Ready for use.

No Colour = Shield already locked; dispose needle.


4. PERFORMING THE INJECTION

Caution: Do not have the thumb on the syringe plunger during needle insertion, to prevent inadvertent loss of medication.

Without hesitation and in one continuous motion, insert the shielded needle at a right angle (pinch when required) and push until the shield is completely retracted.

Maintain constant pressure against the skin while depressing the plunger, until the prescribed dose has been completely delivered.

Lift the needle away from the skin. The visual indication will disappear confirming that the shielding mechanism is automatically locked.

After use, dispose needle and syringe safely as medical waste in a sharps disposal container and/or according to health institution policies.

PRECAUTIONS FOR STORAGE

Store between 2°C and 30°C. Keep dry. Keep away from sunlight. Fragile, handle with care.

REPORT OF INCIDENT

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

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