
Needle Gauge	Needle Length	Needle Outer Diameter	Needle Length	Wall Thickness	Bevel Angle
19G	1"	1.1 mm	25 mm	Thin Wall (TW)	12°
21G	5/8"	0.8 mm	16 mm	Regular Wall (RW)	12°
	1"		25 mm	Thin Wall (TW)	
23G	5/8"	0.6 mm	16 mm	Thin Wall (TW)	12°
	1"		25 mm		
25G	5/8"	0.5 mm	16 mm	Thin Wall (TW)	12°
	1"		25 mm		
26G	1/2"	0.45 mm	12 mm	Regular Wall (RW)	12°
27G	3/4"	0.4 mm	20 mm	Regular Wall (RW)	12°

LOT	Batch code	Do not use if package is damaged	Non pyrogenic
REF	Catalogue number	Fragile, handle with care	Single Sterile Barrier System
#	Contents	Keep away from sunlight	Sterilized using ethylene oxide
Date of Manufacture	Date of Manufacture	Keep dry	Temperature limit
Do not re-sterilize	Do not re-sterilize	Manufacturer	Use-by date
Do not re-use	Do not re-use	MD Medical Device	

- ▶ Store between 2°C and 30°C. Keep dry. Keep away from sunlight. Fragile, handle with care.
- ▶ Do not use if unit package is damaged.
- ▶ For single use only. Do not reuse. Do not re-sterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- ▶ After use, dispose of safely as medical waste in a sharps disposal container and/or according to health institution policies. The product is biohazardous and is physically hazardous due to its sharp edge.
- ▶ Not validated for intravitreal injection.

- ▶ The needle is made of stainless steel containing nickel and cobalt. Cobalt is classified as CMR* 1B and is present in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
- * CMR = Carcinogenic, mutagenic or toxic to reproduction (CLP Regulation EU 1272/2008)
- ▶ If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

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