

# $^{99m}\text{Tc}$ -MAA Overestimates the Absorbed Dose to the Lungs in Radioembolisation: A Quantitative Evaluation in Patients Treated with $^{166}\text{Ho}$ Microspheres

Elschot, M. et al. *Eur. J. Nucl. Med. Mol. Imaging* 2014

## Background

- Radiation pneumonitis is a rare but serious complication of embolisation of liver tumors
- Estimation of the mean absorbed dose to the lungs based on pre-treatment diagnostic Technetium  $^{99m}\text{Tc}$  macro aggregated albumin (Tc- $^{99m}\text{Tc}$ -MAA) imaging, has not been evaluated compared to pre-treatment diagnostic  $^{166}\text{Ho}$  (166-Ho)

## Objective


To compare the performance of pretreatment diagnostic Tc- $^{99m}\text{Tc}$ -MAA imaging and pretreatment diagnostic  $^{166}\text{Ho}$  microsphere imaging for lung absorbed dose estimating in  $^{166}\text{Ho}$  radioembolisation.

## Methods

- A prospective clinical study included 14 patients with chemorefractory, unresectable liver metastases treated with  $^{166}\text{Ho}$  radioembolisation
- Tc- $^{99m}\text{Tc}$ -MAA-based and  $^{166}\text{Ho}$  microsphere-based estimation of lung absorbed doses was performed on pre-treatment diagnostic planar scintigraphy and SPECT/CT images
- The clinical analysis was preceded by an anthropomorphic torso phantom study with simulated lung shunt fractions of 0 to 30 % to determine the accuracy of the image-based lung absorbed dose estimates after  $^{166}\text{Ho}$  radioembolisation

## Results

- Mean dose of  $^{166}\text{Ho}$  scout dose was 250MBq
- The median actual lung absorbed dose was 0.02 Gy, based on post-treatment imaging
- Lung absorbed doses estimated on the basis of pre-treatment diagnostic  $^{166}\text{Ho}$  microsphere imaging were better predictors of actual lung absorbed dose than Tc-99m-MAA based imaging
- Lung absorbed doses estimated using SPECT CT were more accurate than estimations based on planar scintigraphy



## CONCLUSION

This study shows that Tc-99m-MAA pre-treatment imaging overestimates the lung absorbed dose, expected in  $^{166}\text{Ho}$  radioembolisation.

### Key Takeaway

- The  $^{166}\text{Ho}$  scout dose has been shown to be **more accurate** than the commonly used surrogate Tc-99m-MAA at **predicting lung shunting**