en	SYRINGE WITH FIXED NEEDLE	2-3
fr	SERINGUE AVEC AIGUILLE SERTIE	4-5
de	SPRITZE MIT EINGESCHWEISSTER KANÜLE	6-7
es	JERINGA CON AGUJA PREMONTADA	8-9
pt	SERINGA COM AGULHA PREMONTADA	10-11
it	SIRINGA CON AGO FISSO	12-13
nl	SPIJT MET VASTE INJECTIENAALD	14-15
sv	INJEKTIONSSPRUTA MED INTEGRERAD KANYL	16-17
da	SPRØJTE MED INTEGRERET KANYLE	18-19
no	SPRØYTE MED INTEGRERT KANYLE	20-21
fi	RUISKU KIINTEÄLLÄ KANYYLILLÄ	22-23
el	ΣΥΡΙΓΓΑ ΜΕ ΜΗ ΑΠΟΣΠΩΜΕΝΗ ΒΕΛΟΝΑ	24-25
ru	ШПРИЦ СО ВСТРОЕННОЙ ИГЛОЙ	26-27
pl	STRZYKAWKA POŁĄCZONA Z IGŁĄ	28-29
hu	TŰVEL SZERELT FECSKENDŐ	30-31
cs	STŘÍKAČKA S PERMANENTNĚ FIXOVANOU JEHLOU	32-33
sk	STRIEKAČKA SO ZABUDOVANOU IHLOU	34-35
tr	SABIT İĞNELİ ŞIRINGA	36-37
et	SÚSTAL FIKSEERITUD NÕELAGA	38-39
lv	ŠĻIRCE AR FIKSĒTU ADATU	40-41
lt	ŠVIRKŠTAS SU NENUIMAMA ADATA	42-43
sl	BRIZGA S FIKSNO IGLO	44-45
sr	ŠPRIC SA FIKSIRANOM IGLOM	46-47
ro	SERINGA CU AC FIX	48-49
bg	СПРИЦОВКА С ВГРАДЕНА ИГЛА	50-51
uk	ШПРИЦ ІЗ ФІКСОВАНОЮ ГОЛКОЮ	52-53
hr	ŠPRICA S FIKSNOM IGLOM	54-55

Instructions for Use

ENGLISH

SYMBOL EXPLANATION

Batch code	Do not re-use	Medical Device
Catalogue number	Do not use if package is damaged and consult instructions for use	Non pyrogenic
Consult Instructions for Use	Fragile, handle with care	Single Sterile Barrier System
Contents	Keep away from sunlight	Sterilized using ethylene oxide
Date of Manufacture	Keep dry	Temperature limit
Do not re-sterilize	Manufacturer	Use-by date

PRODUCT DESCRIPTION

Non-Toxic

Needle Gauge	Needle Length	Needle Outer Diameter	Needle Length	Wall Thickness	Syringe Nominal Capacity/ Nominal Volume (v)	Dead Space Volume
29G	1/2"	0.33 mm	12 mm	Regular Wall (RW)	1 ml	0.01 ml
Tolerance on Graduated Capacity						
Syringe Nominal Capacity/ Nominal Volume (v)		When syringe filled to a volume less than half of the nominal capacity		When syringe filled to a volume equal or greater than half of the nominal capacity		
1 ml		± (1.50% of v + 2% of expelled volume)		2% of expelled volume		

INTENDED PURPOSE

TERUMO SYRINGE with Fixed Needle for 600 IU/ml FSH/LH (Follicle Stimulating Hormone/Luteinizing Hormone) is a hypodermic syringe intended for manual aspiration of Menopur Multidose 600 IU/ml FSH/LH and for the injection of this solution into parts of the body below the surface of the skin.

INDICATIONS

The device is indicated for subcutaneous injection of fertility hormones to treat women suffering from hormonal deficiencies.

CONTRAINDICATIONS

No contraindications.

PATIENT TARGET GROUP

Women suffering from hormonal deficiencies.

INTENDED USERS

Healthcare professional or lay person.

CLINICAL BENEFIT

The device has an indirect clinical benefit (indirect medical effect) since it is used for subcutaneous injection of fertility hormones to treat women suffering from hormonal deficiencies.

WARNINGS


- Do not use if unit package is damaged.
 - Use immediately after opening the unit package.
 - For single use only. Do not reuse. Do not re-sterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
 - The needle is made of stainless steel containing nickel and cobalt. Cobalt is classified as CMR* 1B and is present in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
- *CMR = Carcinogenic, mutagenic or toxic to reproduction (CLP Regulation EU 1272/2008)

PRECAUTIONS

- The syringe is intended for use immediately after filling, as it is not suitable for containing fluids over extended periods of time.
- If the needle is bent or damaged, no attempt should be made to straighten the needle or use the product.
- After use, dispose of safely as medical waste in a sharps disposal container and/or according to health institution policies. The product is biohazardous and is physically hazardous due to its sharp edge.

Instructions for Use**ENGLISH****DIRECTIONS FOR USE**

1. Peel open blister
2. Remove cap
3. Push plunger down
4. Fill syringe
5. Remove air bubbles
6. Inject
7. Do not re-cap
8. Dispose of safely

**PRECAUTIONS FOR STORAGE**

Store between 2°C and 30°C. Keep dry. Keep away from sunlight.
Fragile, handle with care.

REPORT OF INCIDENT

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

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