**TERUMO**

CE 0482  
2022-08 (Rev.0)

# Ultimaster Nagomi™

Sirolimus eluting coronary stent system

Ultimaster Nagomi Sirolimus eluting coronary stent system  
Système de stent coronaire à libération de sirolimus Ultimaster Nagomi  
Sirolimus freisetzendes Koronar-Stentsystem Ultimaster Nagomi  
Sistema de stent coronario liberador de sirolímús Ultimaster Nagomi  
Sistema de stent coronário com eluição de sirolímús Ultimaster Nagomi  
Sistema di stent coronarico a rilascio di farmaco sirolimus Ultimaster Nagomi  
Ultimaster Nagomi Sirolimus-eluerend coronair stentsysteem  
Ultimaster Nagomi Sirolimusläkemedelseluerande koronarstentsystem  
Ultimaster Nagomi Sirolimus-afgivende koronarstent  
Ultimaster Nagomi Sirolimusfrigivende koronarstentsystem  
Ultimaster Nagomi Sirolimus-sepelvaltimostenttijärjestelmä  
Σύστημα στεφανιαίας ενδοπρόθεσης με έκλυση φαρμάκου sirolimus Ultimaster Nagomi  
Сиролимус-выделяющий коронарный стент Ultimaster Nagomi  
System stentu wieńcowego Ultimaster Nagomi uwalniającego lek sirolimus  
Ultimaster Nagomi Sirolimus-kibocsátó koronárisztent-rendszer  
Koronární stentový systém Ultimaster Nagomi uvolňující léčivo sirolimus  
Ultimaster Nagomi Systém koronárneho stentu uvoľňujúci sirolimus  
Ultimaster Nagomi Sirolimus salınımlı koroner stent sistemi  
Ultimaster Nagomi Sirolimust elueeriv koronaarstendi süsteem  
Ultimaster Nagomi Sirolimu eluėjošo koronāro artēriju stentu sistēma  
Vaistą sirolimuzą išskiriančio koronarinio stento sistema „Ultimaster Nagomi“  
Sistem koronarne opornice Ultimaster Nagomi, ki sprošča sirolimus  
Ultimaster Nagomi Koronarni stent sistem koji oslobađa lek sirolimus  
Ultimaster Nagomi Sistem de stent coronarian cu eliberare de sirolimus  
Елуираща сиролимус коронарна стент система Ultimaster Nagomi  
Стент-система коронарна сиролімус-елютуюча Ultimaster Nagomi  
Ultimaster Nagomi Sustav koronarnog stenta za izlučivanje sirolimusa

DE34M402-01

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
## 1. DESCRIPTION

Ultimaster Nagomi Sirolimus eluting coronary stent system with rapid exchange balloon delivery system consists of a balloon expandable intra-coronary L605 cobalt chromium (CoCr) stent with abluminal drug eluting coating, that consists of a blend of sirolimus and poly (D,L-lactide-co-caprolactone), pre-mounted onto a high pressure, semi-compliant balloon delivery catheter.

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### 3. DEVICE COMPONENT DESCRIPTION

#### 3.1. STENT PLATFORM

The Ultimaster Nagomi Sirolimus eluting coronary stent system platform with open cell design, made of medical grade L605 CoCr alloy (table 2), is available in three stent models, small, medium and large (table 1):

Table 1

Stent model	Nominal Stent diameter (mm)	Nominal Stent length (mm)	Post-dilatation limits (I.D.)* (mm)
Small	2,0; 2,25; 2,5	9; 12; 15; 18; 21; 24; 28; 33; 38; 44; 50	3,5
Medium	2,75; 3,0	9; 12; 15; 18; 21; 24; 28; 33; 38; 44; 50	4,5
Large	3,5; 4,0; 4,5	9; 12; 15; 18; 21; 24; 28; 33; 38; 44; 50	6,25

\* Do not dilate the stent beyond these limits.

#### 3.2. STENT COATING LAYER

##### 3.2.1. Sirolimus – Active Pharmaceutical Ingredient

Sirolimus (figure 2 and tables 2-3) is a drug with potent anti-proliferative, anti-inflammatory and immunosuppressive effects. It acts by binding to the cytosolic receptor FK506-binding-protein-12 (FKBP-12). The complex that is formed between sirolimus and FKBP-12 inhibits the activation of mammalian target of rapamycin (mTOR), which in turn causes cell cycle arrest (progression from phase G1 to S). The typical targets of sirolimus are activated T lymphocytes, while other target cells are smooth muscle cells (SMC) and endothelial cells. Sirolimus inhibits the proliferation and the migration of SMCs and shows an antiproliferative effect on endothelial cells. Sirolimus also inhibits several phases of the restenosis cascade such as inflammation, neointimal hyperplasia formation, total protein and collagen synthesis. Sirolimus as ancillary medicinal substance is intended to reduce restenosis when using Ultimaster Nagomi stent for treatment of stenotic or occlusive lesions in coronary arteries.

##### 3.2.2. Degradable Polymers – Inactive Ingredients

Ultimaster Nagomi stent has a coating consisting of two layers: a primer layer and a drug matrix layer. The coating is applied abuminally, leaving the luminal side of the stent free from drug as such enhancing endothelial coverage. The primer layer and drug carrier polymer, poly(D,L-lactide-co-caprolactone) copolymer (figure 3 and tables 2-3) is expected to degrade within 3-4 months. Upon degradation, the polyactides release non-toxic lactic acid, which is eventually converted to water and carbon dioxide. Polycaprolactone degrades to ε-caprolactone and is completely excreted from the body through urine and feces.

##### 3.3. BALLOON DELIVERY CATHETER

The delivery catheter is a high pressure, semi-compliant balloon delivery catheter that has two radiopaque markers, which fluoroscopically mark the ends of the stent to facilitate proper stent placement. The active balloon length is closely sized to the length of the stent to prevent over-expansion of the tissue proximal or distal to the stent. At the proximal end of the system is a female luer lock connector hub. This hub connects to the balloon inflation lumen. The guide wire enters the distal tip of the catheter and exits 25 cm proximal to the tip. The surface of catheter is partially coated with hydrophilic polymer coating which generates lubricity when wet.

Table 2: Overall qualitative and quantitative information of the Ultimaster Nagomi stent platform, drug and coating

Material	Substance	Low (w/w%)	High (w/w%)	Stent mass (mg)
Cobalt Chromium alloy (ASTM F90)	Cobalt**	Balance *		See table 3
	Chromium	19,00	21,00	
	Nickel	9,00	11,00	
	Carbon	0,05	0,15	
	Manganese	1,00	2,00	
	Silicon	0	0,40	
	Phosphorus	0	0,040	
	Sulfur	0	0,030	
	Tungsten	14,00	16,00	
	Iron	0	3,00	
	Sirolimus	See table 3		
Poly(D,L-lactide-co-caprolactone) copolymer 90:10 (D,L-lactide: caprolactone)	See table 3			

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### 5. INTENDED PURPOSE / INDICATIONS

#### 5.1. Intended purpose

Ultimaster Nagomi Sirolimus eluting coronary stent system is a sterile, single use implantable drug (sirolimus) eluting stent mounted on a semi-compliant balloon delivery catheter, intended to improve myocardial blood flow in patients with stenotic or occlusive lesions in coronary arteries.

#### 5.2. Indications

Ultimaster Nagomi Sirolimus eluting coronary stent system is indicated for treatment of patients with stenotic or occlusive lesions in coronary arteries, including but not limited to patients with chronic coronary syndrome, acute coronary syndrome (STEMI, NSTEMI and unstable angina), diabetes mellitus, multivessel disease, bifurcation lesions, patients older than 65 years, male and female patients, patients with totally occluded lesions, long lesions, lesions residing in small coronary vessels, restenotic lesions including in-stent restenosis, ostial lesions, lesions in left main coronary artery and arterial or venous bypass graft.

The Ultimaster Nagomi stent system is suitable for both femoral and radial approach.

### 6. CONTRA-INDICATIONS

#### 6.1. Contra-indications

- Patients in whom anti-platelet and/or anti-coagulant therapy is contraindicated
- Patients with known allergy to L605 cobalt-chromium alloy and nickel
- Patients with known hypersensitivity to sirolimus or its structurally related compounds, to lactide polymers and caprolactone polymers
- Patients with known hypersensitivity to contrast agent that cannot be controlled prophylactically prior to Ultimaster Nagomi stent implantation
- Patients with extreme vessel tortuosity that may impair stent placement

### 7. PATIENT TARGET GROUP

Patients with stenotic or occlusive lesions in coronary arteries.

#### Special population

The safety and effectiveness of Ultimaster Nagomi Sirolimus eluting coronary stent system has not been established in women who are pregnant or who might be pregnant, in breastfeeding women and in pediatric patients.

The risks and benefits of sirolimus eluting stents should however be considered for each patient before implanting Ultimaster Nagomi stent. Physicians are responsible for assessing patient appropriateness for undergoing stent implantation prior to procedure.

### 8. INTENDED USER

Only physicians who have received appropriate training shall perform implantation of the stent.

### 9. CLINICAL BENEFIT

To improve myocardial blood flow by reducing stenosis or occlusion of arteries supplying blood and oxygen to myocardium to relieve or reduce severity of angina and sparing myocardial muscle and improving outcomes in patients with myocardial infarction.

### 10. COMPLICATIONS

Potential adverse events associated with coronary stent placement include but are not limited to:

- Abrupt vessel closure
- Acute myocardial infarction
- Allergic reaction to anti-coagulation and/or anti-thrombotic therapy, contrast material, or stent and/or delivery system materials or any other PCI mandatory medication
- Aneurysm
- Arrhythmias, including ventricular fibrillation and ventricular tachycardia
- Arteriovenous fistula
- Cardiac tamponade
- Cardiogenic shock
- Death
- Emboli, distal (air, tissue or thrombotic emboli)
- Emergent Coronary Artery Bypass Surgery
- Failure to deliver the stent to the intended site
- Fever
- Heart failure
- Hematoma
- Hemorrhage, requiring transfusion
- Hypotension / Hypertension
- Infection and pain at insertion site
- Ischemia, myocardial
- Myocardial infarction
- Nausea and vomiting
- No reflow
- Prolonged angina
- Pseudoaneurysm
- Renal failure
- Respiratory failure
- Restenosis of stented segment
- Rupture of native and bypass graft
- Stent compression
- Stent embolization
- Stent migration
- Stent thrombosis / occlusion
- Stroke / cerebrovascular accident
- Thrombosis (acute, subacute, or late)
- Total occlusion of coronary artery
- Unstable or stable angina pectoris
- Vessel dissection
- Vessel perforation
- Vessel spasm

Potential adverse events that may be associated with sirolimus drug and polymer coating. Sirolimus administration is limited to intra-coronary stent delivery. Therefore adverse events are not fully characterized but are considered to be consistent to those noted in sirolimus oral administration including:

- Abnormal liver function tests
- Anemia
- Arthralgias
- Changes in lipid metabolism which may include hypertriglyceridemia or hypercholesterolemia
- Diarrhea
- Hypersensitivity to the drug (sirolimus or its excipients) or to the polymer (or individual components) including anaphylactic / anaphylactoid type of reactions
- Hypokalemia
- Immune suppression, especially in patients with hepatic insufficiency or who are taking medications that inhibit CYP3A4 or P-glycoprotein
- Infections
- Interstitial lung disease
- Leukopenia
- Lymphoma and other malignancies
- Myalgia
- Thrombocytopenia

Because of the low systemic exposure to sirolimus after stent implantation, it is very unlikely that any of the adverse events (apart from hypersensitivity reaction) associated with oral administration of sirolimus, will occur.

22

### 11. WARNINGS AND PRECAUTIONS

#### 11.1. Warnings

- Judicious selection of patients is necessary since Percutaneous Coronary Intervention with the use of stents carries the risk of stent thrombosis, vascular complications and/or bleeding events. Hence patients should be properly premedicated and maintained on clinically adequate post-procedural antiplatelet therapy (aspirin and thienopyridine, or appropriate antiplatelet agents.)
- Any advancement after introduction of the delivery catheter into the vessel should be done under high resolution fluoroscopy. When resistance is felt during manipulation, determine the cause of the resistance before proceeding.
- Proper judgment is necessary to select lesion for direct stenting since insufficiently prepared lesion may lead to stent dislodgement.
- Ensure that the aluminum pack and blister pouch have not been damaged or opened as this may compromise the device stability and the sterile barrier.
- This device is to be used under fluoroscopy, all protective measures regarding radiation are to be respected.

#### 11.2. Precautions

##### 11.2.1. Stent Handling - Precautions

- FOR SINGLE USE ONLY. DO NOT REUSE. DO NOT RESTERILIZE. DO NOT REPROCESS. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- Do not use a device that has reached or exceeded its expiry date.
- Ensure that the blister pouch has not been damaged or opened as this may compromise the sterile barrier.
- Use immediately after opening the blister pouch.
- The entire operation should be carried out aseptically.
- Do not use if the stent is exposed to abnormal rubbing or contact with objects other than the guiding catheter or opened hemostatic valve prior to implantation.
- Do not rub or scrape the stent coating.
- Do not displace or remove stent on or from its delivery system as it may damage the stent and/or lead to stent embolization. Stent system is intended to perform as a system.
- Stent should not be used in conjunction with other delivery systems.
- Delivery system should not be used in conjunction with other stents.
- Special care must be taken not to handle or in any way disrupt the stent on the balloon. This is most important during catheter removal from packaging, catheter removal from holder, removal of protector sheath from stent, catheter placement over guide wire and advancement through rotating hemostatic valve adapter and guiding catheter hub.
- Do not "roll" the mounted stent with your fingers as this action may loosen the stent from the delivery balloon.
- Do not expose system to organic solvent. Use only the appropriate balloon inflation media. Do not use air or any gaseous medium to inflate the balloon as this may cause uneven expansion and difficulty in deployment of the stent.
- Exposing the stent to fluids before implantation is not recommended. Exposure to fluids prior to implantation may result in premature release of drug.
- Do not attempt to straighten the proximal shaft (hypotube) as it may cause the catheter to break if it is accidentally bent.
- \*\*This device contains Cobalt (CAS N°7440-48-4, EC N°231-158-0), classified as CMR† 1B, in a concentration above 0.1 % weight per weight. Current scientific evidence supports that medical devices manufactured from alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

†CMR: carcinogenic, mutagenic and toxic to reproduction (CLP regulation EU 1272/2008)

##### 11.2.2. Stent Placement - Precautions

- Do not introduce negative pressure, or pre-inflate delivery system prior to stent deployment other than as directed.
- Always select an appropriate size of the stent as an undersized stent may result in inadequate expansion of the lesion while an oversized stent may lead to inadequate expansion of the stent or damage to the vessel wall.
- Always verify whether the stent is well apposed against vessel wall because incomplete stent apposition may lead to stent thrombosis.
- When treating multiple lesions in the same vessel, stent the distal lesion prior to stenting the proximal lesion. Stenting in this order avoids crossing the proximal stent with the distal stent and reduces the chances for dislodgement.
- Do not expand the stent if it is not properly positioned in the vessel. (See Stent System Removal –Precautions)
- Placement of a stent has the potential to compromise side branch patency.
- Do not exceed rated burst pressure as indicated on the device label. Use of pressure higher than specified may result in balloon rupture with possible intimal damage and dissection.
- Advance the stent delivery system slowly and carefully as excessive force to the delivery stent system can potentially result in stent dislodgement or damage to the stent and/or delivery system.
- An unexpanded stent may be retracted into the guiding catheter one time only. Subsequent movement in and out through the distal end of the guiding catheter should not be performed as the stent may be damaged or dislocated. In case of stent dislodgement, stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the coronary vasculature and/or the vascular access site. Complication may include bleeding, hematoma or pseudoaneurysm.

##### 11.2.3. Stent / System Removal - Precautions

Stent introduction into the coronary artery is limited to one time only as dislodgement may occur. Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system pre-stent implantation, carefully attempt to pull the stent delivery system back through the guiding catheter. If resistance is felt in doing so, or if resistance is felt during the removal of the stent delivery system post-stent deployment, the delivery system and guiding catheter must be removed as a single unit.

When removing the delivery system and guiding catheter as a single unit:

- Do not attempt to retract an unexpanded stent into the guiding catheter while engaged in the coronary arteries. Stent damage or dislodgement may occur.
- Position the proximal balloon marker just distal to the tip of the guiding catheter.
- Advance the guide wire into the coronary anatomy as far distally as safely possible.
- Tighten the rotating hemostatic valve to secure the delivery system to the guiding catheter, then remove the guiding catheter and delivery system as a single unit.
- Failure to follow these steps and/or applying excessive force to the delivery system can potentially result in stent dislodgement or damage to the stent and/or delivery system.
- It is necessary to retain guide wire position for subsequent artery/lesion access, leave the guide wire in place and remove all other system components.

##### 11.2.4. Post Implant - Precautions

- Care must be exercised when crossing a newly deployed stent with a coronary guide wire, IVUS catheter, OCT/OFDI catheter, balloon or other stent delivery system to avoid disruption of the stent geometry.
- Patients should be maintained on clinically adequate post-procedural antiplatelet therapy (aspirin, thienopyridine or other appropriate antiplatelet agents) according to the current guidelines. In case of need, dual antiplatelet therapy can be discontinued earlier, but not before one month.
- Carefully complete the implant card enclosed with the device and hand over to the patient. Patients should be instructed to carry the implant card at all times.
- The stent is a permanent implant and is not intended to be removed.

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**11.2.5. Magnetic Resonance Imaging (MRI)**

Non-clinical testing has demonstrated the expanded Ultimaster Nagomi stent is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 tesla and 3 tesla, with
- Maximum spatial field gradient of 57 T/m
- Maximum force product of 102 T<sup>2</sup>/m
- Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the expanded Ultimaster Nagomi stent is expected to produce a maximum temperature rise of less than 5.9°C (1.9 W/kg, 1.5 tesla) RF-related temperature increase with a background temperature increase of <math>-0.9^{\circ}\text{C}</math> (1.9 W/kg, 1.5 tesla) 5.7°C (1.4 W/kg, 3 tesla) RF-related temperature increase with a background temperature increase of <math>-1.3^{\circ}\text{C}</math> (1.4 W/kg, 3 tesla) after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 10.1 mm from the expanded Ultimaster Nagomi stent when imaged with a gradient echo pulse sequence and a 3 tesla MR system.

**11.2.6. Drug interaction**

Drugs that act through the same binding protein (FKBP) may interfere with the efficacy of sirolimus.

Sirolimus is metabolized by CYP3A4. Strong inhibitors of CYP3A4 (e.g. ketoconazol) might cause increased sirolimus exposure to levels associated with systemic effects, especially if multiple stents are deployed. Systemic exposure of sirolimus should also be taken into consideration if the patient is treated concomitantly with systemic immunosuppressive therapy. Based on the results of the human pharmacokinetic study (PK) the systemic effect of sirolimus after single stent implantation is considered negligible<sup>1</sup>.

<sup>1</sup> Stojkovic et al. (2015) *Fundam Clin Pharmacol.* 29:95-105

**11.2.7. Carcinogenicity, Genotoxicity**

Carcinogenicity studies for sirolimus conducted in mouse and rat showed increased incidences of lymphomas (male and female mouse), hepatocellular adenoma and carcinoma (male mouse) and granulocytic leukaemia (female mouse) due to its chronic immunosuppression. However, because of the low/short systemic exposure to sirolimus after stent implantation, the potential for carcinogenicity is not considered to be of any concern to the safe use of the Ultimaster Nagomi stent. A 4 – and 26 week implantation study conducted in male and female rabbits did not show any local and systemic changes including preneoplastic changes. Sirolimus was not mutagenic in the *in vitro* bacterial reverse mutation assays, the Chinese hamster ovary cell chromosomal aberration assay, the mouse lymphoma cell forward mutation assay, or the *in vivo* mouse micronucleus assay. A genotoxicity study (bacterial reverse mutation test) showed that the Ultimaster Nagomi stent is not genotoxic.

**12. PRECAUTIONS FOR STORAGE**

Keep dry, keep away from sunlight

STORE THE DEVICE BETWEEN 1 - 30°C IN THE ALUMINUM PACK.

The device is packed under oxygen free conditions.

Aluminum pack includes an oxygen absorber and a desiccant. Discard them without opening.

After opening the aluminum pack, use the device within 12 hours.

Do not store the device in the blister pouch.

Humidity levels are kept low inside the blister pouch by the presence of a desiccant.

**13. REPORT OF INCIDENT**

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

**14. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE**

For the summary of safety and clinical performance (SSCP), please visit <https://ec.europa.eu/tools/EUDAMED> (BASIC UDI-DI: 5413206UMNCE).

**15. HOW SUPPLIED / DISPOSAL**

STERILE AND NON PYROGENIC in undamaged and unopened blister pouch. This device is sterilized by e-beam.

CONTENTS: One Terumo balloon expandable, Sirolimus eluting coronary stent mounted on a rapid exchange delivery system. One flushing needle.

DISPOSAL: After use, dispose of safely as medical waste according to health institution policies. The delivery catheter is biohazardous and the flushing needle is physically hazardous due to its sharp edge.

**16. DIRECTIONS FOR USE**

**16.1. Inspection Prior to Use**

- Carefully inspect the stent delivery system package for damage to the sterile barrier. Prior to using the stent system, carefully remove the system from the package and inspect for bends, kinks, and other damage.

**16.2. Materials Required / Quantity Material**

- Guiding catheter with 1.42 mm (0.056") of minimum inner diameter is appropriate for stents with  $\phi$ 2.0 to 4.0 mm.
- Guiding catheter with 1.80 mm (0.071") of minimum inner diameter is appropriate for stents with  $\phi$ 4.5 mm.
- 2 – 3 syringes (10 – 20 ml)
- 1,000 u/500 ml Heparinized Normal Saline (HepNS)
- 0.36 mm (0.014") x 175 cm (minimum length) guide wire
- Rotating hemostatic valve with appropriate minimum inner diameter [2.44 mm (0.096")]
- Diluted contrast medium 1:1 with heparinized normal saline (HepNS)
- Inflation device
- Pre-deployment dilatation catheter
- Three-way stopcock
- Torque device
- Guide wire introducer
- Appropriate arterial sheath
- Appropriate anticoagulation and antiplatelet drugs

**16.3. Preparation**

Guide wire Lumen Flush

Step Action

1. Carefully remove the stent system from its holder, then remove the stent protector sheath from over the stent.

**CAUTION** Carefully slide off the protector sheath from the stent by fixing the sheath at the distal end between the thumb and finger while gently pulling sheath and attached stylet.

2. Verify that the stent is centered on the balloon and located between the radiopaque balloon markers.

**CAUTION** Do not use if any defects are noted.

3. Flush the guide wire lumen with HepNS using the flushing needle supplied with Ultimaster Nagomi stent system. Insert the flushing needle into the tip of the catheter and flush until solution exits the guide wire port.

**CAUTION** Avoid manipulating the stent while flushing the guide wire lumen as this may dislocate the stent on the balloon.

**16.4. Delivery Procedure**

Step Action

1. Prepare vascular access site according to standard practice.
2. Predilate lesion with PTCA catheter.
3. Remove the PTCA catheter.
4. Open rotating hemostatic valve on the guiding catheter as widely as possible.
5. Backload delivery system onto proximal portion of guide wire while maintaining guide wire position across target lesion.

**CAUTION** Confirm that the guide wire OD does not exceed 0.36 mm (0.014"). If a large size guide wire has been used, exchange the wire in the standard manner.

6. Advance stent delivery system over guide wire to target lesion. Utilize radiopaque balloon markers to position stent across lesion; perform angiography to confirm stent position.

**CAUTION** Take care not to damage the delivery catheter and stent when advancing the delivery catheter over the guide wire.

**16.5. Deployment Procedure**

Step Action

1. Before deployment reconfirm the correct position of the stent relative to the target lesion via the catheter markers.
2. Attach the inflation device to the delivery catheter hub and apply negative pressure to purge the balloon of air.
3. Under fluoroscopic visualization, inflate the balloon to at least nominal pressure and maintain for 15-30 seconds to deploy the stent but do not exceed the labeled rated burst pressure (see label on packaging or enclosed compliance sheet).
4. Optimal expansion requires the stent to be in full contact with the artery wall, and with the stent internal diameter matching the size of the reference vessel diameter.
5. Stent wall contact should be verified through routine angiography or intravascular imaging.
6. Deflate the balloon by pulling a vacuum with the inflation device. Make sure the balloon is fully deflated before any attempted movement of the catheter.
7. Confirm adequate stent expansion by angiographic injection through the guiding catheter.

**16.6. Removal Procedure**

Step Action

1. Ensure balloon is fully deflated.
2. Fully open rotating hemostatic valve.
3. While maintaining guide wire position, withdraw stent delivery system.  
**Note:** Should unusual resistance be felt at any time during either lesion access or removal of the stent delivery system before stent implantation, the entire system should be removed. See stent/system removal – Precautions section for specific stent delivery system removal instructions. The catheter clip is recommended to roll the delivery catheter for easier disposal.
4. Tighten rotating hemostatic valve.
5. Repeat angiography to assess stented area. If an adequate expansion has not been obtained, exchange back to the rapid exchange delivery catheter or exchange to another balloon catheter of appropriate balloon diameter to achieve proper stent apposition to the vessel wall. The Ultimaster Nagomi stent should not be expanded to a diameter beyond the post-dilatation limits indicated in table 1.
6. Final stent diameter should match reference vessel diameter.

**16.7. INSTRUCTION FOR SIMULTANEOUS USE OF TWO DEVICES IN GUIDING CATHETER(KISSING BALLOON TECHNIQUE)**

6Fr Compatibility –Any combination of one Ultimaster Nagomi Sirolimus eluting coronary stent system (I.D. 2.0mm – 4.5mm) and one PTCA balloon catheter (for example Accuforce (2.00-5.00 mm) can be used simultaneously within a 6Fr (I.D.:1.8mm) guiding catheter. The technique can be performed as per the instructions listed below:

1. Insert Ultimaster Nagomi Sirolimus eluting coronary stent system using the instructions provided.
2. Insert one balloon catheter, track to the target site and inflate the balloon.
3. Removing the catheters: Remove one catheter and its associated guide wire completely prior to removing the other catheter and its associated guide wire.

**CAUTION** Care should be taken when introducing, torquing and removing one or both devices to avoid entanglement.

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All brand names are trademarks or registered trademarks of TERUMO CORPORATION and their respective owners.
**AU Sponsor** **TERUMO AUSTRALIA PTY LTD**
**Ashitaka Factory of Terumo Corporation**