Ultimaster Nagomi Sirolimus eluting coronary stent system and Magnetic Resonance Imaging (MRI) Examination

ENGLISH (English).....	2
FRANÇAIS (French)	3
DEUTSCH (German)	4
ESPAÑOL (Spanish).....	5
PORTUGUÊS (Portuguese)	6
ITALIANO (Italian)	7
NEDERLANDS (Dutch).....	8
SVENSKA (Swedish).....	9
DANSK (Danish).....	10
NORSK (Norwegian).....	11
SUOMI (Finnish)	12
ΕΛΛΗΝΙΚΑ (Greek)	13
РУССКИЙ (Russian)	14
POLSKI (Polish).....	15
MAGYAR (Hungarian)	16
ČESKÝ (Czech).....	17
SLOVENSKÝ (Slovak).....	18
TÜRKÇE (Turkish).....	19
EESTI (Estonian)	20
LATVISKI (Latvian).....	21
LIETUVIŠKAI (Lithuanian)	22
SLOVENŠČINA (Slovenian)	23
SRPSKI (Serbian).....	24
ROMÂNĂ (Romanian)	25
БЪЛГАРСКИ (Bulgarian)	26
УКРАЇНСЬКА (Ukrainian)	27
HRVATSKI (Croatian).....	28

ENGLISH (English)

Ultimaster Nagomi Sirolimus eluting coronary stent system and Magnetic Resonance Imaging (MRI) Examination

To whom it may concern:

We, Terumo Europe N.V., Interleuvenlaan 40, 3001 Leuven, Belgium, being the manufacturer of the **Ultimaster Nagomi Sirolimus eluting coronary stent system**, hereby state that the Ultimaster Nagomi stents are MR conditional.

Non-clinical testing has demonstrated the expanded Ultimaster Nagomi is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 tesla and 3 tesla, with
- Maximum spatial gradient field of 57 T/m
- Maximum force product of 102 T²/m
- Theoretically estimated maximum whole body averaged (WBA) specific absorption rate (SAR) of 2 W/kg (Normal Operating Mode)

Under the scan conditions defined above, the expanded Ultimaster Nagomi stent is expected to produce a maximum temperature rise of less than

- 5.9°C (1.9 W/kg, 1.5 tesla) RF-related temperature increase with a background temperature increase of $\approx 0.9^\circ\text{C}$ (1.9 W/kg, 1.5 tesla) after 15 minutes of continuous scanning
- 5.7°C (1.4 W/kg, 3 tesla) RF-related temperature increase with a background temperature increase of $\approx 1.3^\circ\text{C}$ (1.4 W/kg, 3 tesla) after 15 minutes of continuous scanning

In non-clinical testing, the image artifact caused by the device extends approximately 10.1 mm from the expanded Ultimaster Nagomi stent when imaged with a gradient echo pulse sequence and a 3 tesla MR system.

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