**MD 0197**

en	SAFETY HYPODERMIC NEEDLE	2-3
fr	AIGUILLE HYPODERMIQUE DE SÉCURITÉ	4-5
de	HYPODERMISCHE SICHERHEITSKANÜLE	6-7
es	AGUJA HIPODÉRMICA DE SEGURIDAD	8-9
pt	AGULHA HIPODÉRMICA DE SEGURANÇA	10-11
it	AGO IPODERMICO DI SICUREZZA	12-13
nl	HYPODERMISCHE VEILIGHEIDSNAALD	14-15
sv	SÄKERHETSKANYL	16-17
da	SIKKERHEDSKANYLE	18-19
no	HYPODERMISK SIKKERHETSKANYLE	20-21
fi	HYPODERMINEN TURVANEULA	22-23
el	ΥΠΟΔΕΡΜΙΚΗ ΒΕΛΟΝΑ ΑΣΦΑΛΕΙΑΣ	24-25
ru	БЕЗОПАСНАЯ ПОДКОЖНАЯ ИГЛА	26-27
pl	BEZPIECZNA IGŁA DO INIEKCJI	28-29
hu	BIZTONSÁGI EGYSZER HASZNÁLATOS TŰ	30-31
cs	PODKOŽNÍ JEHLA S BEZPEČNOSTNÍM CHRÁNIČEM	32-33
sk	BEZPEČNOSTNÁ HYPODERMICKÁ IHĽA	34-35
tr	GÜVENLİ HİPODERMİK İĞNE	36-37
et	HÜPODERMILINE TURVANÖEL	38-39
lv	HIPODERMISKĀ DROŠĪBAS ADATA	40-41
lt	SAUGI POODINĖ ADATA	42-43
sl	VARNOSTNA IGLA ZA PODKOŽNO INJICIRANJE	44-45
sr	SIGURNOSNA HIPODERMIČKA IGLA	46-47
ro	AC HIPODERMIC DE SIGURANȚĂ	48-49
bg	БЕЗОПАСНА ХИПОДЕРМИЧНА ИГЛА	50-51
uk	БЕЗПЕЧНА ГОЛКА ДЛЯ ПІДШКІРНИХ ІН'ЄКЦІЙ	52-53
hr	SIGURNOSNA POTKOŽNA IGLA	54-55

Instructions for Use

ENGLISH

SYMBOL EXPLANATION

	Batch code		Do not re-use		Medical Device
	Catalogue number		Do not use if package is damaged and consult instructions for use		Non pyrogenic
Rx only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician		Fragile, handle with care		Single Sterile Barrier System
	Consult Instructions for Use		Keep away from sunlight		Sterilized using ethylene oxide
	Contents		Keep dry		Temperature limit
	Date of Manufacture		Manufacturer		Use-by date
	Do not resterilize				

PRODUCT DESCRIPTION

Non-Toxic

Needle Gauge	Needle Length	Needle Outer Diameter	Needle Length	Wall Thickness	Bevel Angle
19G	1"	1.1 mm	25 mm	Thin Wall (TW)	12°
19G	1 1/2"	1.1 mm	38 mm	Thin Wall (TW)	12°
21G	1 1/2"	0.8 mm	38 mm	Thin Wall (TW)	12°
25G	5/8"	0.5 mm	16 mm	Thin Wall (TW)	12°

2

INTENDED PURPOSE

The SurGuard™3 Safety Hypodermic Needle is intended for use in the aspiration and injection of fluids for medical purposes. The Safety Hypodermic Needle is compatible for use with standard luer slip and luer lock syringes. Additionally, after withdrawal of the needle from the body, the attached needle safety shield shall be manually activated to cover the needle immediately after use to minimize risk of accidental needle stick.

CLINICAL BENEFIT

The SurGuard™3 Safety Hypodermic Needle has an indirect clinical benefit (indirect medical effect) since it is used for injection or aspiration of fluids. For general application. The safety shield has an indirect clinical benefit since it prevents needle stick injury.

WARNINGS

- Do not use if unit package is damaged.
- Use immediately after opening the unit package.
- Do not re-cap.
- The needle is made of stainless steel containing nickel and cobalt. Cobalt is classified as CMR* 1B and is present in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

*CMR = Carcinogenic, mutagenic or toxic to reproduction (CLP Regulation EU 1272/2008)

INDICATIONS

The SurGuard™3 Safety Hypodermic Needle is for general application - for treatment (injection of fluids) or diagnosis (aspiration).

CONTRA-INDICATIONS

No contra-indications.

PATIENT TARGET GROUP

Intended for general application.

INTENDED USERS

Healthcare professional or lay person.

PRECAUTIONS

- If the needle is bent or damaged, no attempt should be made to straighten the needle or use the product.
- After use, dispose of safely as medical waste in a sharps disposal container and/or according to health institution policies. The product is biohazardous and is physically hazardous due to its sharp edge.

Instructions for Use

ENGLISH


DIRECTIONS FOR USE

Aseptic technique, proper skin preparation and continued protection of the site are essential. Observe Universal Precautions on ALL patients.

PRECAUTION Handle with care to avoid needle sticks. Keep hands behind the needle at all times during use and disposal.


STEP 1.

Tighten the syringe to the needle using aseptic technique. **Grip the base of the needle**, not the safety shield, push and turn the syringe clockwise.

**STEP 2.**

Move the safety shield away from the needle and toward the syringe barrel **to the angle shown**, prior to removing the needle cap. Perform injection procedure according to established technique.

(Note: The needle "bevel up" position is oriented to the position of the safety shield.)


PRECAUTIONS FOR STORAGE


Store between 2°C and 30°C. Keep dry. Keep away from sunlight. Fragile, handle with care.

STEP 3.

After the injection procedure, **use a one-handed technique to activate the safety mechanism** using any of the three (3) methods illustrated below. (Activation is verified by an audible and/or tactile "click", and can be visually confirmed.)

**STEP 4.**

Visually confirm that the safety shield is fully engaged. Do not attempt to deactivate the safety device by forcing the needle out of the safety shield.

**STEP 5.**

After use, dispose of safely as medical waste in a sharps disposal container and/or according to health institution policies.

REPORT OF INCIDENT

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

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