en	INJECTION FILTER NEEDLE	2-3
fr	AIGUILLE-FILTRE D'INJECTION	4-5
de	INJEKTIONSKANÜLE MIT FILTER	6-7
es	AGUJA DE FILTRO PARA INYECCIÓN	8-9
pt	AGULHA DE INJEÇÃO COM FILTRO	10-11
it	AGO CON FILTRO PER INIEZIONE	12-13
nl	INJECTIEFILTERNAALD	14-15
sv	FILTERKANYL FÖR INJEKTION	16-17
da	INJEKTIONSFILTERKANYLE	18-19
no	INJEKSJONSFILTERNÅL	20-21
fi	INJEKTIOSUODATINNEULA	22-23
el	ΒΕΛΟΝΑ ΦΙΛΤΡΟΥ ΕΓΧΥΣΗΣ	24-25
ru	ИГЛА ИНЪЕКЦИОННАЯ С ФИЛЬТРОМ	26-27
pl	IGŁA FILTRUJĄCA DO INIEKCJI	28-29
hu	INJEKCIÓS SZŰRŐTŰ	30-31
cs	INJEKČNÍ FILTRAČNÍ JEHLA	32-33
sk	INJEKČNÁ IHLA S FILTROM	34-35
tr	ENJEKSİYON FİLTRE İĞNESİ	36-37
et	FILTRIGA SÜSTENÕEL	38-39
lv	INJEKCIJAS ADATA AR FILTRU	40-41
lt	INJEKCIŅĖ FILTRAVIMO ADATA	42-43
sl	INJEKCIJSKA IGLA S FILTROM	44-45
sr	INJEKCIJONA IGLA SA FILTEROM	46-47
ro	AC PENTRU FILTRUL DE INECȚIE	48-49
bg	ИГЛА ЗА ИНЖЕКТИРАНЕ С ФИЛТЪР	50-51
uk	ІН'ЄКЦІЙНА ГОЛКА З ФІЛЬТРОМ	52-53
hr	IGLA S FILTROM ZA UBRIZGAVANJE	54-55

Instructions for Use

ENGLISH

SYMBOL EXPLANATION

LOT	Batch code		Do not re-use	MD	Medical Device
REF	Catalogue number		Do not use if package is damaged and consult instructions for use		Non pyrogenic
	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician		Fragile, handle with care		Single Sterile Barrier System
	Consult Instructions for Use		Keep away from sunlight		Sterilized using ethylene oxide
	Contents		Keep dry		Temperature limit
	Date of Manufacture		Liquid filter with a pore size of 5 µm		Use-by date
	Do not resterilize		Manufacturer		

PRODUCT DESCRIPTION

Non-Toxic

Needle Gauge	Needle Length	Needle Outer Diameter	Needle Length	Wall Thickness	Bevel Angle	Filter pore size
30G	1/2"	0.3 mm	12 mm	Extra Thin Wall (ETW)	12°	5 µm

INTENDED PURPOSE

The Terumo Injection Filter Needle is a sterile medical device for single use intended to inject fluids.

INDICATIONS

The device is indicated for general application - for treatment - injection of fluids into parts of the body below the surface of the skin and into the vitreous.

CONTRA-INDICATIONS

No contra-indications.

PATIENT TARGET GROUP

Intended for general application.

INTENDED USERS

Healthcare professional.

CLINICAL BENEFIT

The device has an indirect clinical benefit (indirect medical effect) since it is used for general application - for treatment - injection of fluids into parts of the body below the surface of the skin and into the vitreous.

COMPLICATIONS

The use of the Terumo Injection Filter Needle for injection of fluids into the vitreous may be accompanied but not limited to the following:

- Endophthalmitis
- Damage to the retina
- Visual disturbance

WARNINGS

- Do not use if unit package is damaged.
- Use immediately after opening the unit package.
- For single use only. Do not re-use. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- The needle is made of stainless steel containing nickel and cobalt and lubricated with silicone oil. Cobalt is classified as CMR* 1B and is present in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

*CMR = Carcinogenic, mutagenic or toxic to reproduction (CLP Regulation EU 1272/2008)

PRECAUTIONS

- Do not use for aspiration or withdrawal of fluids.
- In case of injection of fluids into the vitreous, the Terumo Injection Filter Needle shall be used only in combination with (pre-)filled syringes validated for intravitreal use.
- If the needle is bent or damaged, no attempt should be made to straighten the needle or use the product.
- After use, dispose of safely as medical waste in a sharps disposal container and/or according to health institution policies. The product is biohazardous and is physically hazardous due to its sharp edge.

DIRECTIONS FOR USE


Aseptic technique, proper preparation and continued protection of the injection site are essential. Observe Universal Precautions on ALL patients.

PRECAUTION Handle with care to avoid needle sticks. Keep hands behind the needle at all times during use and disposal.

1. Preparing for use - opening

Inspect the blister. Discard the product if the blister and/or the needle is damaged. Peel the blister open.

CAUTION: Avoid touching the hub after opening the blister.

2. Preparing for use - attaching


Firmly attach the hub onto the (pre-)filled syringe by pushing and twisting the needle.

CAUTION: Do not overtighten.

3. Preparing for use - removing needle cap


Remove the needle cap (needle protector) by pulling it straight off.

Perform priming.


4. Perform injection of fluids according to established technique/protocol

WARNING: Do not re-cap.

5. Disposal

After use, dispose of needle and syringe safely as medical waste in a sharps disposal container and/or according to health institution policies.

PRECAUTIONS FOR STORAGE

Store between 2°C and 30°C. Keep dry. Keep away from sunlight. Fragile, handle with care.

REPORT OF INCIDENT

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

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