
TERUMO

SURFLO™

**WINGED INFUSION SET WITH FILTER
AND NEEDLE PROTECTION**

CE
MD 0197

en	WINGED INFUSION SET WITH FILTER AND NEEDLE PROTECTION	2-3
fr	MICROPERFUSEUR À AILETTES AVEC FILTRE ET PROTECTEUR D'AIGUILLE	4-5
de	PERFUSIONSBESTECK MIT FLÜGELKANÜLE MIT FILTER UND NADELSCHUTZ	6-7
es	PALOMILLA DE INFUSIÓN CON FILTRO Y PROTECCIÓN DE LA AGUJA	8-9
pt	DISPOSITIVO COM ASAS PARA INFUSÃO COM FILTRO E PROTETOR DE AGULHA	10-11
it	AGO A FARFALLA CON FILTRO E COPRIAGO DI PROTEZIONE	12-13
nl	INFUSIESET MET VLEUGELNAALD, FILTER EN NAALDPROTECTOR	14-15
sv	INFUSIONSSET MED VINGAR, FILTER OCH KANYLSKYDD	16-17
da	INFUSIONSSÆT MED VINGE, FILTER OG KANYLEBESKYTTELSE	18-19
no	VINGET INFUSJONSSETT MED FILTER OG NÅLEBESKYTTELSE	20-21
fi	SIIPINEULALLA VARUSTETTU INFUUSIOSETTI SUODATTIMELLA JA NEULASUOJUKESELLA	22-23
el	ΣΕΤ ΕΓΧΥΣΗΣ ΜΕ ΠΕΤΑΛΟΥΔΑ ΜΕ ΦΙΛΤΡΟ ΚΑΙ ΠΡΟΣΤΑΤΕΥΤΙΚΟ ΒΕΛΟΝΑΣ	24-25
ru	УСТРОЙСТВО ДЛЯ ИНФУЗИИ С ИГЛОЙ «БАБОЧКА», ФИЛЬТРОМ И ПРОТЕКТОРОМ ИГЛЫ	26-27
pl	MOTYLKOWY ZESTAW INFUZYJNY Z FILTREM I ZABEZPIECZENIEM IGŁY	28-29
hu	SZÁRNYAS INFÚZIÓS KÉSZLET, SZŰRŐVEL ÉS TÚVÉDELEMMEL	30-31
cs	INFUZNÍ SOUPRAVA S FILTREM A JEHLOU S KŘÍDLY A CHRÁNIČEM	32-33
sk	KRÍDIELKOVÝ INFÚZNY SET S FILTROM A CHRÁNIČOM IHLY	34-35
tr	FİLTRELİ VE İĞNE KORUMALI KANATLI İNFÜZYON SETİ	36-37
et	LIBLIKÕEL INFUSIOONIKOMPLEKT, FILTRI JA NÕELAKAITSEGA	38-39
lv	TAURIŅVEIDA INFŪZIJAS KOMPLEKTS AR FILTRU UN ADATAS AIZSARDZĪBU	40-41
lt	PETELIŠKĖS FORMOS INFUZIJOS RINKINYS SU FILTRU IR ADATOS APSAUGA	42-43
sl	INFUZIJSKI SET S KRILCI TER FILTROM IN ZAŠČITO IGLE	44-45
sr	INFUZIONI SET SA KRILCIMA, FILTEROM I ZAŠTITOM ZA IGLU	46-47
ro	MICROPERFUZOR CU FILTRU, ARIPIOARE ŞI PROTECTOR DE AC	48-49
bg	ИНФУЗИОНЕН НАБОР С КРИЛЦА, С ФИЛТЪР И ПРЕДПАЗИТЕЛ НА ИГЛАТА	50-51
uk	ІНФУЗІЙНА СИСТЕМА-МЕТЕЛІК З ФІЛЬТРОМ І ЗАХИСТОМ ГОЛКИ	52-53
hr	KOMPLET ZA INFUZIJU S KRILCIMA, FILTROM I ZAŠTITOM IGLE	54-55

Instructions for Use


ENGLISH

SYMBOL EXPLANATION

LOT	Batch code	Do not re-use	MD	Medical Device
REF	Catalogue number	Do not use if package is damaged and consult instructions for use	Non pyrogenic	
Rx only	Caution: Federal law (USA) restricts this device to sale by or on the order of a physician	Fragile, handle with care	Single Sterile Barrier System	
	Consult Instructions for Use	Keep away from sunlight	STERILIZED	Sterilized using ethylene oxide
	Contents	Keep dry	2°C ^{30°C}	Temperature limit
	Date of Manufacture	Liquid filter with a pore size of 20 µm	Use-by date	
	Do not re-sterilize	Manufacturer		

PRODUCT DESCRIPTION

Non toxic

Needle Gauge (a)	23G 25G	Needle Length (d)	19 mm	Tube Length (l)	13.8"
Needle Length (b)	3/4"	Wall Thickness	Thin Wall (TW)	Dead Space Volume (V)	≤ 0.20 ml
Needle Outer Diameter (c)	0.6 mm 0.5 mm	Tube Length (l)	35 cm	Filter Pore Size	20 µm

INTENDED PURPOSE

The Surflo™ Winged Infusion Set with Filter and Needle Protection (Surshield™) is intended to access the peripheral vascular system, for single-dose or short-term intravenous administration of fluids using a syringe or another compatible/appropriate device. Additionally, after withdrawal of the needle from the patient's vein, the safety shield (needle protection (Surshield™)) shall be manually activated to cover the needle to minimize risk of accidental needle stick.

INDICATIONS

The Surflo™ Winged Infusion Set with Filter and Needle Protection (Surshield™) is for general application - for treatment (administration of fluids).

CONTRA-INDICATIONS

No contra-indications.

PATIENT TARGET GROUP

Intended for general application.

INTENDED USERS

Healthcare professional or lay person.

CLINICAL BENEFIT

The Surflo™ Winged Infusion Set with Filter and Needle Protection (Surshield™) has an indirect clinical benefit (indirect

medical effect) since it is used for administration of fluids. For general application. The safety shield (needle protection (Surshield™)) has an indirect clinical benefit since it prevents needle stick injury.

WARNINGS

- Do not use if unit package is damaged.
- Use immediately after opening the unit package.
- The needle is made of stainless steel containing nickel and cobalt. Cobalt is classified as CMR* 1B and is present in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
* CMR = Carcinogenic, mutagenic or toxic to reproduction (CLP Regulation EU 1272/2008)

PRECAUTIONS

- Do not use for blood collection.
- After use, secure the needle into the safety shield according to the directions for use.
- If the needle is bent or damaged, no attempt should be made to straighten the needle or use the product.
- After use, dispose of safely as medical waste in a sharps disposal container and/or according to health institution policies. The product is biohazardous and is physically hazardous due to its sharp edge.

DIRECTIONS FOR USE

Aseptic technique, proper skin preparation and continued protection of the site are essential. Observe Universal Precautions on ALL patients.

PRECAUTION Keep hands behind the needle at all times during use and disposal.

SET UP FOR ADMINISTRATION

1. Open the package.
2. Connect a syringe or another compatible/appropriate device and secure.
3. Proceed with priming and confirm that administration fluid comes from the needle.

VENIPUNCTURE AND ADMINISTRATION

4. Flip the safety shield back away from the needle towards the tubing. Grasp wings.
5. Remove the needle protector.
Caution: Care should be taken not to touch the needle.
6. Perform venipuncture and confirm proper positioning in the vessel (*Fig. 1*).
7. Secure the position of winged infusion set per facility protocol.

AFTER USE

8. Remove tape, if present, from wings.
9. Flip the safety shield forward towards the needle (*Fig. 2*). Grasp the wing with your thumb on top of the shield and your index finger under the wing. Completely remove the needle from the puncture site and apply digital pressure to the site using a sterile gauze pad (*Fig. 3*).
10. While maintaining the same grip on the device (per step 9 above), activate the safety measure by performing either of the following methods:

► Recommended (one-handed) method:

(Reference OSHA 1910.1030 (d)(2)(vii)(B))







Slide your thumb forward over the shield while pinching your finger and thumb together (or press the safety shield against a hard surface such as a bedside table) until an audible click is heard and the safety feature is visually confirmed to be activated (*Fig. 4a*).

OR
► Alternate (two-handed) method:

While keeping both hands behind the needle, grasp one of the wing tabs with your free hand to stabilize the device. (If necessary, reposition your original grip per step 9 being careful not to press against an exposed area of the needle.) Then, slide your thumb forward over the shield while pinching your finger and thumb together until an audible click is heard and the safety feature is visually confirmed to be activated (*Fig. 4b*).

Caution: Needle must be carefully and firmly placed in the lock position of the safety shield.

11. Do not attempt to deactivate the safety device by separating the needle from the safety shield.
12. Dispose of safely as medical waste according to health institution policies (*Fig. 5*).

PRECAUTIONS FOR STORAGE

Store between 2°C and 30°C. Keep dry. Keep away from sunlight. Fragile, handle with care.

REPORT OF INCIDENTS

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

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