

Patient Information Leaflet- Australia

Ultimaster Nagomi Sirolimus eluting coronary stent system

You recently received or will receive a sirolimus eluting coronary stent implanted in a coronary artery of your heart. After the procedure your physician will give you an implant card. Please carry this patient implant card at all times and show it to any medical personnel who may be treating you.

The following information about your stent is important for you to know:

1. When is Ultimaster Nagomi stent used?

Ultimaster Nagomi stent is used to improve blood supply to the heart muscle, alleviating your symptoms, such as chest pain and discomfort or shortness of breath.

2. <u>Product description</u>

The Ultimaster Nagomi stent is a mesh-like tube made of a metal alloy, coated with a drug and a polymer. The stent is placed and expanded against the wall of the coronary artery to increase blood supply to the heart muscle. The polymer allows the drug to be released in a controlled way into the wall of the artery to reduce the chance of a reblockage.

Stent	Material	Substance	Quantity	
component				
Stent	Cobalt	Cobalt*	46.38-56.95 (w/w%)	
platform	Chromium alloy			
-		Chromium	19.00-21.00 (w/w%)	
		Niekol	0.00.11.00(w/w)	
		мске	9.00-11.00 (w/w%)	
l		Carbon	0.05-0.15(w/w%)	
		Manganese	1.00-2.00 (w/w%)	
l		Silicon	0-0.40 (w/w%)	
1			0-0.04 (w/w%)	
		Phosphorus		
		Sulphur	0-0.03 (w/w%)	
		Tungsten	14.00-16.00 (w/w%)	
1		Iron	0-3.00 (w/w%)	
Drug	Sirolimus		35-197 μ g (stents with diameter 2.0; 2.25 or	
l			2.5 mm, depending on stent length)	
			36-198 μ g (stents with diameter 2.75 or 3.0	
l			mm, depending on stent length)	
1			39-197 μg (stents with diameter 3.5; 4.0 or 4.5	
1			mm, depending on stent length)	

Table 1: Information on the materials and substances used in the Ultimaster Nagomi stent:



Polymer	Poly(D,L-lactide-co- caprolactone) copolymer	 48-268 μg (stents with diameter 2.0; 2.25 or 2.5 mm, depending on stent length) 49-269 μg (stents with diameter 2.75 or 3.0 mm, depending on stent length) 57-288 μg (stents with diameter 3.5; 4.0 or 4.5 mm, depending on stent length)

* Caution: This device contains Cobalt (CAS N°7440-48-4), classified as CMR 1B, in a concentration above 0.1% weight per weight. Current scientific evidence supports that medical devices manufactured from alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects. (CMR 1B: carcinogenic and toxic to reproduction (CLP regulation EU 1272/2008)

Caution:

Inform your physician in case of:

- possible allergy to Cobalt-Chromium alloy and nickel
- possible allergy to sirolimus or related compounds, to lactide polymers and caprolactone polymers
- pregnancy/breastfeeding

3. Potential adverse events

Adverse events are listed in alphabetical order. Please consult your physician for further explanation on the potential adverse events.

Potential adverse events associated with the procedure of <u>coronary stent placement</u> include but are not limited to:

Abrupt vessel closure	Infection and pain at insertion site	
Acute myocardial infarction	Ischemia, myocardial	
Allergic reaction to anti-coagulation	Myocardial infarction	
and/or anti-thrombotic therapy, contrast	Nausea and vomiting	
material, or stent and/or delivery system	Prolonged angina	
materials or any other PCI mandatory	No reflow	
medication	Pseudoaneurysm	
Aneurysm	Renal failure	
Arrhythmias, including ventricular	Restenosis of stented segment	
fibrillation and ventricular tachycardia	Respiratory failure	
Arteriovenous fistula	Rupture of native and bypass graft	
Cardiac tamponade	Stent compression	
Cardiogenic shock	Stent embolization	
Death	Stent migration	
Emboli, distal (air, tissue or thrombotic	Stent thrombosis / occlusion	
emboli)		
Emergent Coronary Artery Bypass Surgery	Stroke / cerebrovascular accident	
Failure to deliver the stent to the intended	Thrombosis (acute, subacute, or late)	
site		



Fever	Total occlusion of coronary artery	
Heart failure	Unstable or stable angina pectoris	
Hematoma	Vessel dissection	
Hemorrhage, requiring transfusion	Vessel perforation	
Hypotension / Hypertension	Vessel spasm	

Potential adverse events that may be <u>associated with sirolimus drug</u> are listed below and are consistent for oral intake of sirolimus. In case of stent implantation however, sirolimus is only locally delivered into the vessel wall and amounts of the drug released into the blood are low. Therefore, it is very unlikely that any of the below adverse events (apart from an allergic reaction) associated with oral intake of sirolimus will occur:

Abnormal liver function tests				
Anemia				
Arthralgias				
Changes in lipid metabolism which may include hypertriglyceridemia or				
hypercholesterolemia				
Diarrhea				
Hypersensitivity including anaphylactic / anaphylactoid type of reactions				
Hypokalemia				
Immune suppression, especially in patients with hepatic insufficiency or who are taking				
medications that inhibit CYP3A4 or P-glycoprotein				
Infections				
Interstitial lung disease				
Leukopenia				
Lymphoma and other malignancies				
Myalgia				
Thrombocytopenia				

If you experience any of the above listed adverse events, please report it to your treating physician.

4. Symptoms that could indicate that the device is malfunctioning

If you experience any symptom other than what your physician told you before the procedure, immediately seek medical attention. Symptoms of potential device malfunction after procedure may include sudden shortness of breath and chest pain.

5. <u>Medication and examinations after stent placement</u>

After stent placement it is extremely important that you follow the specific medication regimen your physician prescribed. In case you need to discontinue your medication prematurely, please consult your physician.

Follow the physician's instructions for any further follow-up examinations after stent placement.



6. Lifetime of the implanted stent

This device is a permanent implant and is not intended to be removed.

7. MRI (Magnetic Resonance Imaging) information

The Ultimaster Nagomi stent has been proven to be MRI conditional. This means that you can safely undergo an MRI scan under certain conditions. When you inform your physician on the type of stent you have (e.g. by showing your implant card), your physician will be able to determine these safe scanning conditions.

8. Passing through metal detectors

The Ultimaster Nagomi stent will not set off metal detectors, you can safely pass through metal detectors during security checks (e.g. airport, stores).

9. Incident reporting

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or the Australian sponsor and to the Therapeutic Goods Administration

www.tga.gov.au

Manufacturer:

TERUMO EUROPE N.V. INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM www.terumo-europe.com

Distributor/Sponsor:

Terumo Australia Pty Ltd Macquarie Park, NSW, 2113, Australia T: 1800 837 866 Mangere, Auckland, 2022, New Zealand T: 0800 66 77 57 E: <u>cs_australia@terumo.co.jp</u>



Product Code	Nominal Expanded Stent	Nominal Expanded Stent
	Inner Diameter (mm)	Length (mm)
DE-RS2009ASM	2.00	9
DE-RS2012ASM	2.00	12
DE-RS2015ASM	2.00	15
DE-RS2018ASM	2.00	18
DE-RS2021ASM	2.00	21
DE-RS2024ASM	2.00	24
DE-RS2028ASM	2.00	28
DE-RS2033ASM	2.00	33
DE-RS2038ASM	2.00	38
DE-RS2044ASM	2.00	44
DE-RS2050ASM	2.00	50
DE-RS2209ASM	2.25	9
DE-RS2212ASM	2.25	12
DE-RS2215ASM	2.25	15
DE-RS2218ASM	2.25	18
DE-RS2221ASM	2.25	21
DE-RS2224ASM	2.25	24
DE-RS2228ASM	2.25	28
DE-RS2233ASM	2 25	33
DE-RS2238ASM	2.25	38
	2.25	30
	2.25	
DE-RSZZSUASIVI	2.25	50
DE-RS2509ASM	2.50	9
DE-RS2512ASM	2.50	12
DE-RS2515ASM	2.50	15
DE-RS2518ASM	2.50	18
DE-RS2521ASM	2.50	21
DE-RS2524ASM	2.50	24
DE-RS2528ASM	2.50	28
DE-RS2533ASM	2.50	33
DE-RS2538ASM	2.50	38
DE-RS2544ASM	2.50	44
DE-RS2550ASM	2.50	50
DE-RS2709ASM	2.75	9
DE-RS2712ASM	2.75	12
DE-RS2715ASM	2.75	15
DE-RS2718ASM	2.75	18
DE-RS2721ASM	2.75	21
DE-RS2724ASM	2.75	24
DE-RS2728ASM	2.75	28
DE-RS2733ASM	2.75	33
DE-RS2738ASM	2.75	38
DE-RS2744ASM	2.75	44
DE-RS2750ASM	2.75	50
DE-RS3009ASM	3.00	9
DE-RS3012ASM	3.00	12
DE-RS3015ASM	3.00	15
DE-RS3018ASM	3.00	18
DE-RS3021ASM	3.00	21



Product Code	Nominal Expanded Stent	Nominal Expanded Stent
	Inner Diameter (mm)	Length (mm)
DE-RS3024ASM	3.00	24
DE-RS3028ASM	3.00	28
DE-RS3033ASM	3.00	33
DE-RS3038ASM	3.00	38
DE-RS3044ASM	3.00	44
DE-RS3050ASM	3.00	50
DE-RS3509ASM	3.50	9
DE-RS3512ASM	3.50	12
DE-RS3515ASM	3.50	15
DE-RS3518ASM	3.50	18
DE-RS3521ASM	3.50	21
DE-RS3524ASM	3.50	24
DE-RS3528ASM	3.50	28
DE-RS3533ASM	3.50	33
DE-RS3538ASM	3.50	38
DE-RS3544ASM	3.50	44
DE-RS3550ASM	3.50	50
DE-RS4009ASM	4.00	9
DE-RS4012ASM	4.00	12
DE-RS4015ASM	4.00	15
DE-RS4018ASM	4.00	18
DE-RS4021ASM	4.00	21
DE-RS4024ASM	4.00	24
DE-RS4028ASM	4.00	28
DE-RS4033ASM	4.00	33
DE-RS4038ASM	4.00	38
DE-RS4044ASM	4.00	44
DE-RS4050ASM	4.00	50
DE-RS4509ASM	4.50	9
DE-RS4512ASM	4.50	12
DE-RS4515ASM	4.50	15
DE-RS4518ASM	4.50	18
DE-RS4521ASM	4.50	21
DE-RS4524ASM	4.50	24
DE-RS4528ASM	4.50	28
DE-RS4533ASM	4.50	33
DE-RS4538ASM	4.50	38
DE-RS4544ASM	4.50	44
DE-RS4550ASM	4.50	50