# QuiremSpheres™

## Holmium-166 Microspheres

**REF** QS-V001

<b>EN</b> - Instructions for use	3
<b>CZ</b> - Návod k použití	9
<b>DA</b> - Brugervejledning	15
<b>DE</b> - Gebrauchsanweisung	22
<b>ES</b> - Instrucciones de uso	29
<b>FI</b> - Käyttöohjeet	36
<b>FR</b> - Mode d'emploi	42
<b>HR</b> - Upute za uporabu	49
<b>HU</b> - Használati utasítás	55
<b>IT</b> - Istruzioni per l'uso	62
<b>KA</b> - გამოყენების ინსტრუქცია	69
<b>LT</b> - Naudojimo instrukcijos	77
<b>NL</b> - Gebruiksaanwijzing	83
<b>NO</b> - Bruksanvisning	90
<b>PL</b> - Instrukcja użytkowania	96
<b>PT</b> - Instruções de utilização	103
<b>RO</b> - Instrucțiuni de utilizare	110
<b>RU</b> - Инструкция по применению	117
<b>SK</b> - Návod na použitie	124
<b>SL</b> - Navodila za uporabo	131
<b>SR</b> - Uputstvo za upotrebu	137
<b>SV</b> - Bruksanvisning	144
<b>TR</b> - Kullanım talimatları	150

LC-80043 [09]

## 1. DEVICE DESCRIPTION

QuiremSpheres™ consists of poly-L-lactic acid (PLLA) microspheres that contain holmium-166 in a 2 mL suspension buffer. The microspheres have a mean diameter of 25-35 micrometres. Holmium-166 is a high-energy beta emitting isotope for therapeutic use. The maximum energy of the beta particles is 1.85 MeV (50.0%) and 1.77 MeV (48.7%). The maximum range of the emitted beta particles in tissue is 8.7 mm with a mean of 2.5 mm. In addition, holmium-166 emits primary gamma photons with an energy of 81 KeV (6.7%). The half-life of holmium-166 is 26.8 hours, which means more than 90% of the radiation is delivered within the first 4 days following the administration procedure. At planned moment of treatment, the activity per microsphere is 200-400 Bq. QuiremSpheres™ is a permanent implant.

QuiremSpheres™ is administered into the hepatic artery via a micro-catheter. QuiremSpheres™ will distribute non-uniformly in the liver. This is primarily due to differences in the hepatic arterial flow to the tumour(s) and non-tumour liver tissue, the ratio of vascularity between tumour and non-tumour liver tissue, and tumour size. QuiremSpheres™ can be visualised in-vivo with MRI and/or SPECT.

## 2. INTENDED PURPOSE

QuiremSpheres™ is intended for implantation into hepatic tumours by delivery via the hepatic artery for the treatment of patients with unresectable liver tumours.

## 3. INTENDED USER

The intended user of QuiremSpheres™ includes nuclear medicine and interventional radiology professionals.

Only intended users and other clinical personnel who have completed the QuiremSpheres™ training programme are authorized to order, handle and/or implant QuiremSpheres™ Holmium-166 Microspheres.

## 4. INDICATIONS FOR USE

QuiremSpheres™ is indicated for the treatment of patients with unresectable liver tumours.

## 5. CLINICAL BENEFIT

Based on the available clinical data of QuiremSpheres™, the following benefits are identified:

- Disease control of unresectable liver tumours, resulting in progression-free and overall survival in line with state-of-the-art (patient benefit).
- Well-tolerated treatment in the intended patient population (patient benefit).
- A treatment with limited duration of side-effects (patient benefit).
- Symptomatic relief in patients with symptoms stemming from the tumour (patient benefit).
- Post-treatment verification through SPECT and MRI (user benefit).

## 6. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For the summary of safety and clinical performance (SSCP), please visit <https://ec.europa.eu/tools/eudamed> (BASIC UDI-DI: 8719266QS-V001ZV).

## 7. CONTRAINDICATIONS

QuiremSpheres™ is contraindicated in:

- patients who have a predicted lung dose exceeding 30 Gy in a single treatment.
- patients who have significantly abnormal synthetic and excretory liver function tests (LFTs).
- patients who have had previous external beam radiation therapy to the liver.
- patients who have complete main portal vein thrombosis.
- patients who have uncorrectable extrahepatic deposition of a surrogate marker. Deposition of activity in the falciform ligament, portal lymph nodes and gallbladder is accepted.
- patients who have inadequate kidney function.
- patients who are pregnant.
- children (<18 years).

## 8. ADVERSE EVENTS

Based on the finished clinical studies, the following adverse event profile for QuiremSpheres™ was established:

- Grade 1, 2 events that were observed after SIRT with QuiremSpheres™ in >10% of patients are: nausea, abdominal pain, fatigue, back pain, vomiting, edema limbs, paroxysmal atrial tachycardia, dyspnea, allergic (contrast) reaction, ascites, fever, dizziness.
- Grade 3, 4, 5 events that were observed after SIRT with QuiremSpheres™ in >1-10% of the patients are: hepatic failure, nausea, fatigue, ascites, upper gastrointestinal hemorrhage, thoracic pain, esophageal varices hemorrhage and gastric ulceration. Grade 3,4,5 events observed in >10% of patients are: abdominal pain.

- Grade 1, 2 laboratory events that were observed after SIRT with QuiremSpheres™ in >10% of patients are: y-glutamyl transferase increase, lymphocytopenia, aspartate aminotransferase increase, alkaline phosphatase increase, hyperglycaemia, total bilirubin increase, alanine aminotransferase increase, hyponatremia, anaemia / hemoglobin, thrombocytopenia, albumin decrease, lactate dehydrogenase increase, leucopenia, ammonia increase, INR increased, creatinine increase, hypokalaemia, erythrocytes, prolonged APTT.
- Grade 3, 4 laboratory events that were observed after SIRT with QuiremSpheres™ in >1-10% of the patients are: aspartate aminotransferase increase, alkaline phosphatase increase, hyperglycaemia, total bilirubin increase, alanine aminotransferase increase, hyponatraemia, anaemia / hemoglobin, thrombocytopenia, albumin decrease, lactate dehydrogenase increase.
- Grade 3, 4 laboratory events that were observed after SIRT with QuiremSpheres™ in >10% of the patients are: y-glutamyl transferase increase, lymphocytopenia.

#### Potential Serious Adverse Events Due to High Radiation:

- **Acute pancreatitis** ---- causes immediate severe abdominal pain. Verify by SPECT imaging of the abdomen and test for serum amylase.
- **Acute gastritis** ---- causes abdominal pain. Verify by standard methods to diagnose gastric ulceration.
- **Acute cholecystitis** ---- causes significant upper abdominal pain and may require cholecystectomy for resolution. Verify by appropriate imaging studies.
- **Other ulceration of the gastrointestinal (GI) tract** ---- causes abdominal pain. Verify by standard methods to diagnose ulceration.
- **Radiation hepatitis (REILD)** ---- causes unexplained progressive deterioration of liver function. Verify by transcutaneous core biopsy of the liver.
- **Radiation pneumonitis** ---- causes excessive non-productive cough. Verify by x-ray or HR pulmonary CT evidence of pneumonitis.

#### 9. WARNINGS

No warnings were identified.

#### 10. PRECAUTIONS

The following precautions are drawn attention to:

- Excessive radiation to the normal liver parenchyma may result in radioembolisation-induced liver disease (REILD), also known as radiation hepatitis.
- This product is radioactive. Local regulations must be followed when handling this device.

#### 11. CHEMICAL COMPOSITION

The microspheres contain poly (L-lactic acid), acetylacetone and holmium chloride. The aqueous suspension buffer contains sodium dihydrogen phosphate, di-sodium hydrogen phosphate, poloxamer 188 and water for injection.

#### 12. PREPARATION

##### 12.1. PATIENT SELECTION AND PRE-TREATMENT TESTS

Only patients with unresectable liver tumours may be considered for treatment with QuiremSpheres™. Patient selection for treatment with QuiremSpheres™ requires a medical opinion that treatment of tumour(s) within the liver will result in patient benefit and that the tumour(s) are unresectable.

It is recommended that the following patient tests are performed before treatment:

- A hepatic angiogram to establish arterial anatomy of the liver.
- Serologic tests of liver function to determine the extent of liver function damage.
- A kidney function test to determine the extent of kidney function.
- A nuclear medicine scan of an administered surrogate marker to predict the lung dose, intrahepatic dose and any extrahepatic deposition from the QuiremSpheres™ treatment.
  - As a surrogate marker, QuiremScout™ Holmium-166 microspheres or Technetium-99m MAA can be injected via a catheter into the hepatic artery.
  - To optimise the predictive value for this test, the catheter tip during the injection of either surrogate marker is recommended to be placed at the exact same anatomical position from which QuiremSpheres™ will be administered.
  - Do not place the catheter tip too close to a bifurcation or the arterial wall as this may introduce vortices that limit the predictive value of the surrogate marker.

*! Note ! The patient must recover from any major surgical operations before being treated with QuiremSpheres™.*

*! Note ! The potential effect of systemic-anti cancer agents on liver vascularisation and/or a potential radio-sensitising effect has not been evaluated in combination with QuiremSpheres™.*

*! Note ! Safety and effectiveness of this device in nursing mothers has not been established.*

## 12.2. DOSE PLANNING

As a general principle, the amount of QuiremSpheres™ radioactivity  $A$  [MBq] to achieve a radiation absorbed dose  $D$  [Gy] in a tissue with mass  $M$  [kg] is determined by the following formula:

$$A[\text{MBq}] = D[\text{Gy}] \times M[\text{kg}] \times 63 [\text{MBq/J}] \quad [1]$$

where  $M$  [kg] is determined by the volume  $V$  [L] and the density  $\rho$  [kg/L] of the tissue in the following way:

$$M[\text{kg}] = V[\text{L}] \times \rho [\text{kg/L}] \quad [2]$$

The density  $\rho$  [kg/L] that can be used for determining  $M$  [kg] of liver tissue (both normal and tumorous) is 1.06 kg/L.

The radioactivity to be administered for a patient treatment can be calculated based on the liver volume to be treated (Section 5.2.2.) or based on a dose simulation using a surrogate marker (Section 5.2.3.).

### 12.2.1. RECOMMENDATIONS FOR DOSE PLANNING

The following applies for dose planning:

- The recommended average whole-liver planned absorbed dose is 60 Gy.
- When only one lobe or (sub)segment of the liver is treated, the planned average radiation absorbed dose to that treatment volume may exceed 60 Gy. The average absorbed dose to the whole liver should still be evaluated and the recommended average absorbed dose of 60 Gy still applies.

**! Precaution !** *Excessive radiation to the normal liver parenchyma may result in radioembolisation-induced liver disease (REILD), also known as radiation hepatitis.*

**! Note !** *Excessive shunting to the lungs may lead to radiation pneumonitis.*

**! Note !** *Inadvertent delivery of QuiremSpheres™ to the gallbladder, gastrointestinal tract or pancreas may cause acute gastritis, ulceration of the GI tract, acute pancreatitis or acute cholecystitis.*

### 12.2.2. ACTIVITY CALCULATION BASED ON LIVER VOLUME

This method uses solely the liver volume in which QuiremSpheres™ will be injected as input for calculating the required activity needed to achieve a certain dose on the liver tissue according to formula [1] and [2]. The volume  $V$  is measured based on medical imaging data and is used to determine the mass  $M$  of liver tissue.

Due to non-uniform distribution of QuiremSpheres™ to the tumour and non-tumour liver tissue, it can be expected that a proportionally higher radiation absorbed dose will be delivered to the tumour tissue. Nevertheless, it is recommended that the intrahepatic distribution is studied as part of the work-up using a surrogate marker.

### 12.2.3. ACTIVITY CALCULATION BASED ON DOSE SIMULATION USING A SURROGATE MARKER

This method uses a surrogate marker to simulate the expected radiation absorbed dose for a given treatment activity in tumour tissue and normal liver tissue. As a surrogate marker QuiremScout™ Holmium-166 Microspheres or Technetium-99m MAA can be used. The distribution of either surrogate marker is assessed by SPECT(-CT) imaging based on which the fraction  $F$  of the total activity in a volume of interest (VOI) is determined according to:

$$F = \frac{\text{Counts}_{\text{VOI}}}{\text{Counts}_{\text{Total Volume}}} \quad [3]$$

Where *counts* refer to SPECT counts measured within the corresponding volumes and *Total Volume* refers to the total liver volume to which QuiremSpheres™ activity is planned to be delivered. If a lung shunt is detected, the *Total Volume* also includes the lung volume. The activity required to achieve an absorbed dose  $D$  [Gy] in a VOI is subsequently calculated based on formula [4] in the following way:

$$A = \frac{D_{\text{VOI}}[\text{Gy}] \times M_{\text{VOI}}[\text{kg}] \times 63 [\text{MBq/J}]}{F} \quad [4]$$

With  $M_{\text{VOI}}$  calculated according to formula [2]. The volume  $V$  is measured based on medical imaging data and is used to determine the mass  $M$  of tissue in the VOI for both normal liver tissue and tumour tissue.

## 12.3. LUNG DOSE PREDICTION USING A SURROGATE MARKER

A surrogate marker can be used to simulate the expected radiation absorbed dose in the lungs for a given treatment activity. As a surrogate marker QuiremScout™ Holmium-166 Microspheres or Technetium-99m MAA can be used. The distribution of either surrogate marker can be assessed by SPECT(-CT) imaging or planar scintigraphy based on which the percentage of lung shunting  $L$  [%] is calculated using the following formula:

$$L [\%] = \frac{\text{counts of total lung}}{\text{counts of total lung} + \text{counts of liver}} \times 100 \quad [5]$$

Where *counts* refer to SPECT-CT or planar scintigraphy counts measured on images within the corresponding regions of interest. The radiation absorbed dose to the lung can be predicted based on the measured lung shunt  $L$  [%], the amount of holmium-166 activity  $A$  [MBq] to be administered and the mass of the lungs  $M_{\text{lung}}$  [kg] using the following formula:

$$D_{\text{lung}} [\text{Gy}] = \frac{A[\text{MBq}] \times L/100}{M_{\text{lung}}[\text{kg}] \times 63 [\text{MBq/J}]} \quad [6]$$

Lung tissue mass  $M_{lung}$  [kg] can be calculated using a patient-specific lung tissue density and lung volume based on a CT scan or can be estimated as 1 kg.

### 13. INSTRUCTIONS FOR USE

#### 13.1. PRESENTATION AND STORAGE

Preparation can either be performed in a hospital facility for handling radioactive materials (e.g. a nuclear medicine laboratory) or in a hospital facility (e.g. an angiography room) approved for the QuiremSpheres™ and/or QuiremScout™ administration procedure.

***! Precaution ! The product is radioactive. Local regulations must be followed when handling this device.***

QuiremSpheres™ microspheres are supplied in suspension with a patient-specific activity in up to three capped V-vials. Details of the calibration date, activity at calibration and expiry information are indicated on the package labels. QuiremSpheres™ is supplied sterile and is intended for single-patient, single-use. Each V-vial comes in a lead container, packed in a type A transport box. The vial and its content should be stored inside its transportation container at room temperature (15-25°C, 59-77°F) until use.

***! Note !*** Upon receiving QuiremSpheres™, check for damage of the packaging. Do not use the product in case the V-vial or package is damaged.

***! Note !*** QuiremSpheres™ is for single patient use only.

***! Note !*** Do not use the product after the expiry date and time, as indicated on the product label.

***! Note !*** Before using QuiremSpheres™, verify the activity of the received vial(s) by using a dose calibrator calibrated for measuring holmium-166.

#### 13.2. INTERACTING DEVICES & ACCESSOIRES

Use the QuiremSpheres™ Delivery Set (QS-D001) for the transfer of QuiremSpheres™ from the patient V-vial to the patient catheter.

Use the QuiremSpheres™ Customer Kit (QS-C001) as means for protecting against unintended radiation during the QuiremSpheres™ administration procedure.

It is advised to use a micro-catheter with an inner diameter of at least 0.65 mm. If a catheter of smaller inner diameter is used, a lower flow rate is expected during the administration process. A lower flow rate:

- will lead to a faster rise in the fluid level in the V-vial, which takes longer to normalise. Therefore, allow sufficient time between administration cycles to allow the fluid level in the V-vial to normalise.
- may lead to an increased retention of microspheres in the QuiremSpheres™ Delivery Set and micro-catheter.

***! Note !*** Failure to allow sufficient time between administration cycles to normalise the fluid level in the V-vial may result in excessive pressure in the V-vial, ultimately leading to septal leakage and loss of pressure.

#### 13.3. RADIOLOGICAL PLACEMENT OF CATHETER

For transarterial implantation of QuiremSpheres™ a catheter is inserted either via the femoral or the radial artery under x-ray guidance. This should only be performed by a trained interventional radiologist.

Using standard techniques, position the microcatheter into the hepatic artery (common, right, left or one of the other branches to the liver). The interventional radiologist must be familiar with the frequent arterial variations in the blood supply to the liver and from the liver to the surrounding organs. Every attempt should be made to deliver the microspheres into the hepatic arteries in such a way that radiation is administered to the target area in the liver only. If the tumours are limited to one lobe or (sub)segment of the liver, the catheter can be inserted selectively into the artery supplying the target lobe or (sub)segment, thereby sparing healthy liver tissue.

#### 13.4. ADMINISTRATION OF QUIREMSPHERES™

Administration of QuiremSpheres™ shall be performed in a hospital facility (e.g. an angiography room) approved for the QuiremSpheres™ Holmium-166 microspheres administration procedure.

For an extended and illustrated description of the preparation and execution of the QuiremSpheres™ administration procedure, please refer to the Manual Administration Procedure (LC-80072) and the instructions for use of the QuiremSpheres™ Delivery Set (LC-80076). These instructions should be read and understood in their entirety prior to use.

The following recommendations apply for the administration of QuiremSpheres™

- Delivery of QuiremSpheres™ into the hepatic artery must be performed slowly (0.1 ml per push per second) using saline (0.9% NaCl). Rapid delivery may cause reflux.
- During the delivery procedure, the catheter should be flushed with saline (0.9% NaCl) at regular intervals to prevent blockage.
- Routinely check the catheter tip to ensure it remains in the planned position throughout the administration procedure.
- Check on stasis of blood flow and potential microspheres reflux by intermittent administration of contrast agent to prevent inadvertent delivery of microspheres.

### 13.5. DISPOSAL OF RADIOACTIVE WASTE

The QuiremSpheres™ V-vial, QuiremSpheres™ Delivery Set, (micro-)catheters and other single-use disposables will contain small residual quantities of microspheres and require monitoring for radioactivity. These items are to be disposed of according to local procedures. This may involve storage to decay prior to disposal through the usual facility waste system. All gowns and surgical gear must be monitored at the end of each procedure. Contaminated items should be bagged, labelled and returned to the medical physics department or other designated area for decay until safe for laundering or disposal.

### 13.6. POST-TREATMENT TESTS

It is recommended that a post-treatment SPECT/CT scan of the upper abdomen is performed to review the intrahepatic distribution and extrahepatic deposition of the administered microspheres.

### 14. MAGNETIC RESONANCE IMAGING (MRI) SAFETY

Non-clinical testing has demonstrated that the Holmium-166 Microspheres of QuiremSpheres™ are MR Conditional. A patient with this device may be safely scanned under the conditions shown in Table 1. Failure to follow these conditions may result in harm to the patient. Additional MRI Safety information may be requested via [info.quirem@terumo-europe.com](mailto:info.quirem@terumo-europe.com).

**Table 1 – MR scanning conditions**

**MR Conditional**

Parameter	Condition of use / Information
Static Magnetic Field Strength (B <sub>0</sub> )	1.5 T, 3 T
Static Magnetic Field (B <sub>0</sub> ) Orientation	Horizontal, Cylindrical Bore
Maximum Spatial Field Gradient (SFG)	9 T/m (900 gauss/cm)
RF Polarization	Circular Polarized (CP)
RF Transmit Coil	Any Transmit RF coil may be used.
RF Receive Coil	Any Receive RF coil may be used.
MR System (RF) Operating Modes or Constraints	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	3.2 W/kg (Normal Operating Mode)
B1 <sup>+</sup> <sub>rms</sub>	No maximum B1 <sup>+</sup> <sub>rms</sub>
Scan Duration	2 W/kg whole-body average SAR for 60 minutes of continuous RF (a sequence or back-to-back series/scan without breaks) followed by a cooling period of 15 minutes.
MR Image Artifact	The presence of this implant may produce image artifacts. Some manipulation of scan parameters may be needed to compensate for the artifacts.

### 15. REPORTING INCIDENTS

If an incident related to QuiremSpheres™ occurs, report this incident to Quirem Medical B.V.: [info.quirem@terumo-europe.com](mailto:info.quirem@terumo-europe.com).

If a *serious* incident related to QuiremSpheres™ occurs, report this incident to Quirem Medical B.V and the Competent Authority of the Member State in which the user or patient is established.

### 16. RADIATION SAFETY

The preparation and administration procedure must be regarded as being a potentially serious radiation hazard to the clinical staff. Regulatory and local radiation handling guidelines should be adhered to in terms of set-up, microsphere administration, waste disposal and post-implantation care.

Radiation hygiene principles should be taken into consideration at all times. In short, this means that dose exposure to clinical staff, nursing staff and unintended dose exposure to the patient should be 'as low as reasonably achievable' (ALARA) by considering the following aspects:

- **TIME** – Minimise the time of exposure
- **DISTANCE** – Increase the distance between the radiation source and body/ body extremities as much as possible
- **SHIELDING** – Take appropriate shielding measures

Avoiding contact with urine and bile of patients having undergone SIRT (in collection bags and wound dressings) can minimise the risk of radiation exposure for individuals, particularly during the first 24 hours following implantation.

In case of (suspected) radioactive contamination of staff, equipment or treatment room:

- Determine the extent of the radioactive contamination by measuring the gamma photon emission of holmium-166 with an appropriate hand-held radiation detector.
- Follow applicable local guidelines on cleaning or quarantining the contaminated surface(s).

## APPENDIX I – MEASUREMENT OF POST-IMPLANT EXPOSURE

Dose rates were measured at 1.0 m distance from the abdomen of patients (n=3) implanted with holmium-166 microspheres activity, associated with a recommended average whole-liver absorbed dose of 60 Gy.

The results of the measurements are presented in Table 2.

**Table 2 – Dose rates at 1.0 m distance from treated patients at t = 0, 6, 24 and 48 h after treatment.**

	t = 0 h	t = 6 h	t = 24 h	t = 48 h
Dose rate $\dot{D}$ ( $\mu\text{Sv/h}$ )	31 (18-53)	27 (16-46)	17 (10-29)	9 (5-15)

## APPENDIX II – CORRECTION FOR DECAY


The physical half-life of holmium-166 is 26.8 hours. To calculate the activity of QuiremSpheres™ at any moment after the reference time multiply the reference activity with the appropriate decay factor shown in Table 3.

The reference time must be in the user's local time before correcting for decay.

**Table 3 – Decay factors of QuiremSpheres™**

Hours	Decay Factor	Hours	Decay Factor
0.5	0.987	9	0.792
1	0.974	10	0.772
2	0.950	11	0.752
3	0.925	12	0.733
4	0.902	24	0.538
5	0.879	36	0.394
6	0.856	48 (day 2)	0.289
7	0.834	72 (day 3)	0.155
8	0.813	96 (day 4)	0.083

## APPENDIX III – EXPLANATION OF SYMBOLS ON PACKAGING AND LABELING

	Manufacturer		Temperature limitation		Use before date
	Date of manufacture		Caution, contains radioactive material		This side up
	Serial Number		Product reference code		Keep away from sunlight
	Sterile barrier indication		Sterile using aseptic processing techniques		Keep dry
	Medical Device indication		Unique Device Identifier		Do not use if the package has been damaged or opened
	European Conformity		European time of calibration UTC		Stacking limit is three boxes
	Only for single patient use		Peel-off labels on both sides		Fragile content
	Consult instructions for use		MR Conditional		

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