*QuiremScout*TM

Holmium-166 Microspheres

REF QS-S001

EN - Instructions for use	3
CZ - Návod k použití	8
DA - Brugervejledning	13
DE - Gebrauchsanweisung	18
ES - Instrucciones de uso	23
FI - Käyttöohjeet	28
FR - Mode d'emploi	33
HR - Upute za uporabu	38
HU - Használati utasítás	43
IT - Istruzioni per l'uso	48
KA - გამოყენების ინსტრუქცია	53
LT - Naudojimo instrukcijos	58
NL - Gebruiksaanwijzing	63
NO - Bruksanvisning	68
PL - Instrukcja użytkowania	73
PT - Instruções de utilização	78
RO - Instrucțiuni de utilizare	83
RU - Инструкция по применению	88
SK - Návod na použitie	93
SL - Navodila za uporabo	98
SR - Uputstvo za upotrebu	103
SV - Bruksanvisning	108
TR - Kullanım talimatları	113

Qirem Medical B.V.
Zutphenseweg 55
7418 AH Deventer
The Netherlands
www.qirem.com

1. DEVICE DESCRIPTION

QiremScout™ consists of a small dose of biocompatible poly-L-lactic acid (PLLA) microspheres that contain holmium-166 in a 2 mL suspension buffer. The standard dose contains approximately 3 million microspheres with a maximum total activity of 250 MBq at planned moment of administration. The microspheres have a mean diameter of 25-35 micrometres. Holmium-166 is a beta emitting isotope for therapeutic use. The maximum energy of the beta particles is 1.85 MeV (50.0%) and 1.77MeV (48.7%). The maximum range of the emitted beta particles in tissue is 8.7 mm with a mean of 2.5 mm. In addition, holmium-166 emits primary gamma photons with an energy of 81 KeV (6.7%). The half-life of holmium-166 is 26.8 hours, which means that more than 90% of the radiation is delivered within the first 4 days following the administration procedure. QiremScout™ is a permanent implant.

QiremScout™ is administered into the hepatic artery via a micro-catheter. QiremScout™ will distribute non-uniformly in the liver. This is primarily due to differences in the hepatic arterial flow to the tumour(s) and non-tumour liver tissue, the ratio of vascularity between tumour and non-tumour liver tissue, and tumour size. The microspheres can be visualised in-vivo with SPECT.

2. INTENDED PURPOSE

QiremScout™ is intended for evaluating lung shunt, extrahepatic deposition and intrahepatic distribution of microspheres injected in the hepatic artery in patients eligible for SIRT.

3. INTENDED USER

The intended user of QiremScout™ includes nuclear medicine and interventional radiology professionals.

Only intended users and other clinical personnel who have completed the QiremSpheres™ training programme are authorised to order, handle and/or implant QiremScout™ Holmium-166 Microspheres.

4. INDICATIONS FOR USE

QiremScout™ is indicated for patients that are considered for SIRT of liver tumours.

5. CLINICAL BENEFIT

Based on the Clinical Evaluation of QiremScout™, the following benefits are identified:

- Evaluation of intra- and extrahepatic depositions of intra-arterially injected microspheres, including lung shunting (patient benefit).
- Support informed decision-making about the expected safety and efficacy of the QiremSpheres™ treatment (patient benefit).

6. SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

For the summary of safety and clinical performance (SSCP), please visit <https://ec.europa.eu/tools/eudamed> (BASIC UDI-DI: 8719266QS-S001Z8).

7. CONTRAINDICATIONS

There are no contra-indications for QiremScout™.

8. ADVERSE EVENTS

Based on the finished clinical studies, the following adverse event profile for QiremScout™ was established:

- Grade 1, 2 clinical adverse events that were observed after SIRT with QiremScout™ in >10% of patients are: back pain.
- Grade 3, 4 clinical adverse events that were observed after SIRT with QiremScout™ in >1-10% of patients are: back pain, abdominal pain.
- No grade 5 adverse events have been reported. No laboratory events have been reported.

9. WARNINGS

No hazard has been identified that requires a warning.

10. PRECAUTIONS

This product is radioactive. Local regulations must be followed when handling this device.

11. CHEMICAL COMPOSITION

The microspheres contain poly (L-lactic acid), acetylacetone and holmium chloride. The aqueous suspension buffer contains sodium dihydrogen phosphate, di-sodium hydrogen phosphate, poloxamer 188 and water for injection.

12. PREPARATION

It is recommended that the following tests are performed prior to administration of QuiremScout™:

- A hepatic angiogram to establish arterial anatomy of the liver.
- Serologic tests of liver function to determine the extent of liver function damage.

! Note ! Safety and effectiveness of this device in pregnant women, nursing mothers or children (<18 years) has not been established.

12.1. ACTIVITY CALCULATION

A standard QuiremScout™ dose order consists of two vials with the total activity of 250 MBq split in ratio of 1/3 (80 MBq) and 2/3 (170 MBq).

QuiremScout™ dose order can also be customised by ordering a single vial or multiple (up to three) vials to be divided amongst multiple injection positions per targeted liver volume.

In case of customisation, calculate the fraction per vial (Fr_i) by dividing the targeted liver volume for each vial (V_i) by the total targeted liver volume (V_{total}), i.e.:

$$Fr_i = V_i/V_{total}$$

Consequently, all fractions summed should equal 1, i.e.:

$$\sum Fr_i = 1.$$

Multiply the different fractions times the maximum total dose of 250 MBq.

13. INSTRUCTIONS FOR USE

13.1. PRESENTATION AND STORAGE

Preparation can either be performed in a hospital facility for handling radioactive materials (e.g. a nuclear medicine laboratory) or in a hospital facility (e.g. an angiosuite) approved for the QuiremSpheres™ and/or QuiremScout™ administration procedure.

! Precaution ! The product is radioactive. Local regulations must be followed when handling this device.

QuiremScout™ is supplied in suspension as a standard dose in a single capped V-vial or is divided over multiple (up to three) V-vials. Details of the calibration date, activity at calibration and expiry information are indicated on the package labels. The device is sterile and has been processed using aseptic techniques. Each V-vial(s) comes in a lead container, packed in a type A transport box. The vial and its content should be stored inside its transportation container at room temperature (15-25°C, 59-77°F) until use.

! Note ! Upon receiving QuiremScout™, check for damage of the packaging. Do not use the product in case the V-vial or package is damaged.

! Note ! QuiremScout™ is for single patient use only.

! Note ! Do not use the product after the expiry date and time, as indicated on the product label.

Before using QuiremScout™, verify the activity of the received vials by using a calibrated well-chamber or dosimeter.

13.2. INTERACTING DEVICES & ACCESSOIRES

It is mandatory to use the QuiremSpheres™ Delivery Set (QS-D001) for the transfer of QuiremScout™ from the patient V-vial to the patient catheter.

It is mandatory to use the QuiremSpheres™ Customer Kit (QS-C001) as means for protecting against unintended radiation during the QuiremScout™ administration procedure.

It is advised to use a micro-catheter with an inner diameter of at least 0.65 mm. If a catheter of smaller inner diameter is used, a lower flow rate is expected during the administration process (Chapter 12.4). A lower flow rate:

- will lead to a greater fluid level rise in the V-vial that needs a longer time to normalise. Therefore, allow sufficient time for the fluid level in the V-vial to normalise before switching between administration cycles.
- may lead to higher retainment of microspheres in the QuiremSpheres™ Delivery Set and micro-catheter.

! Note ! Failure to allow sufficient time between administration cycles to normalise the fluid level in the V-vial may result in excessive pressure in the V-vial, ultimately leading to septal leakage and loss of pressure.

13.3. RADIOLOGICAL PLACEMENT OF CATHETER

For transarterial implantation of QuiremScout™ a catheter is inserted either via the femoral or the radial artery under x-ray guidance. This should only be performed by a trained interventional radiologist.

Using standard techniques, position the micro-catheter into the hepatic artery (common, right, left or one of the other branches to the liver). The interventional radiologist must be familiar with the frequent arterial variations in the blood supply to the liver and from the liver to the surrounding organs. Every attempt should be made to deliver the microspheres into the hepatic arteries in such a way that the microspheres are administered to the target area in the liver only. If the tumours are limited to one lobe or (sub)segment of the liver, the catheter can be inserted selectively into the artery supplying the target lobe or (sub)segment.

13.4. ADMINISTRATION OF QUIREMSCOUT™

Administration of QuiremScout™ shall be performed in a hospital facility (e.g. an angiography room) approved for the QuiremSpheres™ and QuiremScout™ administration procedure.

For an extended and illustrated description of the preparation and execution of the QuiremScout™ administration procedure, please refer to the Manual Administration Procedure (LC-80072) and the instructions for use of the QuiremSpheres™ Delivery Set (LC-80076). These instructions should be read and understood in their entirety prior to use.

The following recommendations apply for the administration of QuiremScout™

- Delivery of QuiremScout™ into the hepatic artery must be performed slowly (0.1 ml per push per second) using saline (0.9% NaCl). Rapid delivery may cause reflux.
- During the delivery procedure, the catheter should be flushed with saline (0.9% NaCl) at regular intervals to prevent blockage.
- Routinely check the catheter tip to ensure it remains in the planned position throughout the administration procedure.

13.5. DISPOSAL OF RADIOACTIVE WASTE

The QuiremScout™ V-vial, QuiremSpheres™ delivery set, (micro-)catheters and other single-use disposables will contain small residual quantities of microspheres and require monitoring for radioactivity. These items are to be disposed of according to local procedures. This may involve storage to decay prior to disposal through the usual facility waste system. All gowns and surgical gear must be monitored at the end of each procedure. Contaminated items should be bagged, labelled and returned to the medical physics department or other designated area for decay until safe for laundering or disposal.

14. EVALUATION OF ADMINISTRATION

Assessing the in-vivo QuiremScout™ distribution to evaluate lung shunt, extra hepatic deposition and intra-hepatic distribution can be carried out by means of gamma scintigraphy or Single Photon Emission Computed Tomography (SPECT). To ensure sufficient photon detection, it is recommended that imaging is performed on the same day of injection.

15. MAGNETIC RESONANCE IMAGING (MRI) SAFETY

Non-clinical testing has demonstrated that the Holmium-166 Microspheres of QuiremScout™ are MR Conditional. A patient with this device may be safely scanned under the conditions shown in Table 1. Failure to follow these conditions may result in harm to the patient. Additional MRI Safety information may be requested via info.quirem@terumo-europe.com.

Table 1 – MR scanning conditions

MR Conditional

Parameter	Condition of use / Information
Static Magnetic Field Strength (B ₀)	1.5 T, 3 T
Static Magnetic Field (B ₀) Orientation	Horizontal, Cylindrical Bore
Maximum Spatial Field Gradient (SFG)	9 T/m (900 gauss/cm)
RF Polarization	Circular Polarized (CP)
RF Transmit Coil	Any Transmit RF coil may be used.
RF Receive Coil	Any Receive RF coil may be used.
MR System (RF) Operating Modes or Constraints	Normal Operating Mode
Maximum Whole-Body SAR	2 W/kg (Normal Operating Mode)
Maximum Head SAR	3.2 W/kg (Normal Operating Mode)
B1 ⁺ _{rms}	No maximum B1 ⁺ _{rms}
Scan Duration	2 W/kg whole-body average SAR for 60 minutes of continuous RF (a sequence or back-to-back series/scan without breaks) followed by a cooling period of 15 minutes.
MR Image Artifact	The presence of this implant may produce image artifacts. Some manipulation of scan parameters may be needed to compensate for the artifacts.

16. REPORTING INCIDENTS

If an incident related to QuiremScout™ occurs, report this incident to Quirem Medical B.V.: info.quirem@terumo-europe.com.

If a *serious* incident related to QuiremScout™ occurs, report this incident to Quirem Medical B.V and the Competent Authority of the Member State in which the user or patient is established.

17. RADIATION SAFETY

The preparation and administration procedure must be regarded as being a potentially serious radiation hazard to the clinical staff. Regulatory and local radiation hygiene guidelines should be adhered to in terms of set-up, microsphere administration, waste disposal and post-implantation care.

Radiation hygiene principles should be taken into consideration at all times. In short, this means that dose exposure to clinical staff, nursing staff and unintended dose exposure to the patient should be 'as low as reasonably achievable' (ALARA) by considering the following aspects:

- **TIME** – Minimise the time of exposure.
- **DISTANCE** – Increase the distance between the radiation source and body/body extremities as much as possible.
- **SHIELDING** – Take appropriate shielding measures.

In case of (suspected) radioactive contamination of staff, equipment, or treatment room:

- Determine the extent of the radioactive contamination by measuring the gamma photon emission of holmium-166 with an appropriate hand-held radiation detector.
- Follow applicable local guidelines on cleaning or quarantining the contaminated surface(s).

APPENDIX I – CORRECTION FOR DECAY




The physical half-life of holmium-166 is 26.8 hours. To correct for physical decay of this radionuclide, the fractions that remain at selected intervals after the time of calibration are shown in **Table 2**.

The reference time must be in the user's local time before correcting for decay.

Table 2 – Physical decay table of holmium-166. Half-life: 26.8 hours

Hours	Fraction remaining	Hours	Fraction remaining
0.5	0.987	9	0.792
1	0.974	10	0.772
2	0.950	11	0.752
3	0.925	12	0.733
4	0.902	24	0.538
5	0.879	36	0.394
6	0.856	48 (day 2)	0.289
7	0.834	72 (day 3)	0.155
8	0.813	96 (day 4)	0.083

APPENDIX II – EXPLANATION OF SYMBOLS

	Manufacturer		Temperature limitation		Use before date
	Date of manufacture		Caution, contains radioactive material		This side up
	Serial Number		Product reference code		Keep away from sunlight
	Sterile barrier indication		Sterile using aseptic processing techniques		Keep dry
	Medical Device indication		Unique Device Identifier		Do not use if the package has been damaged or opened
	European Conformity		European time of calibration UTC		Stacking limit is three boxes
	Only for single patient use		Peel-off labels on both sides		Fragile content
	Consult instructions for use		MR Conditional		

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