QS-Q001

- EN - Instructions for Use
- DE - Gebrauchsanweisung
- ES - Instrucciones de uso
- FR - Mode d'emploi
- IT - Istruzioni per l'uso
- NL - Gebruiksaanwijzing
- PL - Instrukcja użycia
- PT - Instruções de utilização
- SV - Bruksanvisning
- TR - Kullanım Talimatları
- CS - Návod k použití

Quirem Medical B.V.
Zutphenseweg 55
7418 AH Deventer
The Netherlands
Website: www.quirem.com

1. INTRODUCTION

This Instructions For Use is a manual for using Q-Suite 1.2 in a safe and responsible manner. Q-Suite 1.2 is a software device developed and manufactured by Quirem Medical B.V. for post-treatment evaluation of holmium-166 selective internal radiation therapy (SIRT). It is assumed that the user has enough knowledge of operating personal computers and the operating system Microsoft (MS) Windows 7 and/or 8.1 and/or 10 to be able to use the relevant functionality of the system.

This manual describes the system, the intended use, safety related issues and the daily operation and interpretation of the systems interface and response. In this manual, you can find important information about how to operate the Q-Suite 1.2 software. The manual helps you with the operation and the maintenance of the device, in a safe and responsible manner. This document should be studied carefully before attempting to use the system.

Throughout this manual certain conventions are used. These are warnings, cautions and notes. They provide a means of prioritizing information to be brought to the attention of the user. They are given as follows:

! Warning ! A warning tells you that there is a risk of personal injury or death

! Caution ! A caution tells you that there is a risk of damage to the device and/or there is a risk of damage to other equipment.

! Note ! A note gives more information

2. DESCRIPTION

Q-Suite 1.2 is a medical image processing software designed to determine and evaluate the radiation absorbed dose following an injection of holmium-166 microspheres. The results give an estimate of the dose distribution within the liver tissue of the patient that was treated with the radioactive microspheres.

Q-Suite 1.2 determines 3-dimensional dose calculations based on multi-slice MR images or 3D SPECT-CT images that were acquired during or after administration of holmium-166 microspheres. When combining these dose maps with anatomical images, the user can measure the dose values within the tissues of interest and store these results into a report.

Input for the dose calculations is an image data set from an MRI examination, with multiple images that are acquired using certain pre-defined parameters, or an image data set from a SPECT examination, acquired and reconstructed using methods optimized for holmium-166 imaging. In addition to these image data sets Q-Suite 1.2 requires treatment related parameters for dose calculation like the specific activity and holmium content of the microspheres that were used for treatment or the time the microspheres were injected. Based on the inputs, Q-Suite 1.2 calculates the distribution of the radioactivity in the liver tissue on a voxel level and subsequently the voxel wise absorbed radiation dose.

For dose measurements within regions of interest the user can load additional MR or CT images, combine these with the calculated dose values, and create the volumes of interest in which the dose must be measured. The results can be stored in a report for further distribution or archiving purposes.

3. INTENDED USE

3.1 INTENDED USE

Q-Suite 1.2 is intended to support evaluation of medical treatment with radioactive holmium-166 microspheres. Medical indications for this treatment technique include primary and secondary cancers and treatment techniques include intra-arterial, intra-tumoral and intra-ductal administration of microspheres. The medical device software calculates the tissue absorbed radiation dose after injection of radioactive holmium-166 microspheres based on medical images.

When interpreted by a trained physician, the results generated by the software provide information that can aid in treatment evaluation.

3.2 INTENDED USER

Q-Suite 1.2 is intended to be used by medically trained staff in a clinical environment. The anticipated users are radiologists, radiotherapists, nuclear medicine physicians, radiation oncologists or other medical persons that are involved in radioembolization treatment.

3.3 INDICATIONS FOR USE

Medical indications for treatment with radioactive holmium microspheres include primary and secondary cancers and treatment techniques include intra-arterial, intra-tumoral and intra-ductal administration of microspheres.

3.4 INTENDED USE ENVIRONMENT

Q-Suite 1.2 is intended to be used in an environment with normal, non-stressful, office conditions.

4. COMPATIBLE DEVICES

Q-Suite 1.2 is a stand-alone software device. Q-Suite 1.2 reads DICOM image data generated by DICOM compliant devices

5. SAFETY AND SECURITY

5.1 PRECAUTIONS

Only qualified staff may use the system. The system manager at the clinic must ensure that individuals authorized to use Q-Suite 1.2 have been suitably trained in using the system by Quirem staff or other competent people.

It is important to ensure that the operator is aware that the quality of the output depends critically on the quality of the input data, and any irregularities or uncertainties about input data units, identification, or quality of any other nature shall be thoroughly investigated before the data is used. It is therefore highly recommended that a medical physics or imaging expert is involved throughout the introduction of Q-Suite 1.2 in the clinic.

The system manager should ensure that no custom settings in MS Windows© make information disappear or become illegible; for instance, when selecting color schemes for title bars, selecting font type and color, etc.

5.2 WARNINGS

Carefully read all the instructions before use. Observe all warnings and cautions throughout these instructions. This manual gives instructions for use of the system. It does not provide information/advice on administering treatment with microspheres or whatsoever.

Warnings are based on hazard analysis performed during the product life cycle. The following major concerns are taken into consideration:

- The system is only intended to support evaluation of treatment with Holmium-166-microspheres, like QuiremSpheres®, and should not be used for evaluation of treatments with Yttrium-90 microspheres or other SIRT treatments.

5.3 SAFETY AND SECURITY

Installation and upgrading of Q-Suite 1.2 shall be done by or under supervision of Quirem Medical B.V. or one of its official distributors. An upgrade of software and/or hardware of the machine Q-Suite 1.2 is installed on, requires that the user conducts a verification test, whether or not Quirem Medical has performed a complete system validation according to regulatory requirements.

Q-Suite 1.2 contains database functionality but it is not intended to be used as a long-term archive of patient data. The user is responsible for storage and backup of in and output data.

The use of Q-Suite 1.2 is protected through licenses.

It is strongly recommended to use anti-virus software for all systems Q-Suite 1.2 is installed on, whether or not connected to the Internet. Anti-virus software and firewalls can have negative impact on the functionality of Q-Suite 1.2 and should be configured correctly.

Data stored in the Q-Suite 1.2 database is not encrypted. It is recommended to restrict access to the computer system on which Q-Suite 1.2 is installed to authorized users only to prevent unintended access to database files. The use of Q-Suite 1.2 is protected by a login code to prevent unauthorized login. System settings are secured by an additional code. The hospital system manager must ensure that each individual has the correct access to Q-Suite 1.2. It is strongly recommended to give the login code only to people that are authorized to use Q-Suite 1.2.

6. INSTRUCTIONS FOR USE

6.1 INSTALLATION AND CONFIGURATION

6.1.1 Installation

Prior to installation, verify that the computer system where Q-Suite 1.2 will be installed on meets the requirements as specified in section 8.

- Install Q-Suite 1.2 using the Q-Suite 1.2 installation DVD. Installation can be started by running the file SuiteSetup.exe (run as administrator). The Q-Suite Setup window will be shown.
- Select 'Install' to start installation. All required components will automatically be installed.
- After successful installation, select 'Close' to exit the installation program.

After correct installation, Q-Suite 1.2 can be found at the following location: C:\Program Files\Quirem Medical\Q-Suite. Run Q-Suite.exe to start the program.

6.1.2 De-installation

Q-Suite 1.2 can be removed from the computer system in the following way.

- Uninstall Q-Suite 1.2 using the Microsoft Windows functionality *Uninstall or Change a Program*

6.1.3 Removing Patient Data

After uninstalling Q-Suite 1.2, patient data in the database can be **permanently** removed in the following way:

- Stop the service *MongoDB_Quirem* using the Microsoft Windows *Task Manager*
- Delete the folder *C:\MongoDB\data* from the local hard drive
- Delete the folder *C:\QuiremMedical\CalculationJobs* from the local hard drive

How to remove patient data without uninstalling Q-Suite 1.2 is described in section 6.10

6.2 PREPARATION OF DATA

All image data to be processed by Q-Suite 1.2 shall be stored in classic DICOM format on a local disk or mapped network drive.

6.2.1 PREPARATION FOR DOSE RECONSTRUCTION

Q-Suite 1.2 can reconstruct radiation absorbed dose values based on input SPECT data or MR data.

Acquisition recommendations MRI

MRI-based holmium-166 dose reconstruction is based on measuring the effect of the presence of the microspheres on MR signal decay. This requires a multi gradient echo (MGE) read-out of the free induction decay (FID) using at least 2 echoes. Data with less than 2 echoes will not be accepted for dose reconstruction by Q-Suite 1.2. Two MGE datasets should be acquired, one prior to and one after SIRT. The images should include the entire liver of the patient and there should be no overlap between image slices.

The accuracy of the output dose results of Q-Suite 1.2 will always depend on the image quality of the input MR images. To maximize the accuracy, acquisition settings should be optimized. Factors that influence the accuracy of the dose calculations include:

- the signal to noise ratio (SNR) of the images (which depends on the MR hardware, acquisition times, resolution, acceleration techniques etc)
- the number of gradient echoes used to sample the FID (minimum of 4 echoes recommended)
- the echo times of the gradient echoes (fast signal decay requires fast sampling)
- the presence of artifacts in the images (e.g. ghosting, fold-over, object induced)
- movement during acquisition (patient and/or organ)

For more background information on the subject of image acquisition see a paper published by van de Maat et al (Eur Radiol 2013;23:827-35).

! Note ! : The quality of the input MR images directly influences the accuracy of the calculated dose.

Acquisition recommendations SPECT

For SPECT-based dose reconstructions, a 3D reconstructed SPECT image is required with an acquisition energy window that includes the holmium-166 photopeak of 80.6 keV. The SPECT image can be accompanied by a multi-slice CT dataset in case of a SPECT-CT acquisition. The accuracy of the output dose results of Q-Suite 1.2 will always depend on the image quality of the input SPECT images. Therefore it is strongly recommended to use a state-of-the-art SPECT-CT scanner and an advanced reconstruction method that includes scatter correction, attenuation correction and modelling the collimator-detector response

! Note ! : The quality of the input SPECT images directly influences the accuracy of the calculated dose.

6.2.2 PREPARATION FOR DOSE EVALUATION

For the purpose of dose evaluation, the Q-Suite 1.2 generated dose reconstructions can be fused with MR or CT images for proper identification of tissues of interest (See also section 6.9). Guidelines for acquisition of these images are not part of this instruction.

6.3 GENERAL WORKING PRINCIPLES

Three main functional areas can be distinguished in Q-Suite 1.2:

- Patient data administration and selection (home screen)
- Dose reconstruction (workflow)
- Dose evaluation (workflow)

On starting Q-Suite 1.2, a home screen is displayed that shows all patient of which image data is stored in the Q-Suite 1.2 database. On selection of a patient, the user can choose to start a dose reconstruction workflow or a dose evaluation workflow for that patient, depending on the data that is available. If data required for a specific workflow is not present in the database, that workflow will not be available. After selection of a workflow, the selected patient becomes the active patient. There can be only one active patient at a time. By closing the active patient, the user returns to the home screen.

In the specific workflows the user is guided through the steps to create the required end results. See section 6.5 and 6.6 for more information on these workflows.

6.4 DATA IMPORT

Via the home screen, new patient image data can be imported to the Q-Suite 1.2 database. For importing data, the user can select a folder to search for data. DICOM data in that folder will be recognized by Q-Suite 1.2 and presented to the user for selection for import. Only Q-Suite 1.2 compatible data will be imported. The software will recognize if data selected for import is already present in the database and will not import that again.

6.5 DOSE RECONSTRUCTION

On starting a dose reconstruction for a selected patient, the user has three options: starting a new MR based dose reconstruction, starting a new SPECT based dose reconstruction or loading a previously stored dose reconstruction composition from the database. The available options depend on the data that is stored in the database.

Each dose reconstruction workflow consists of specific steps that depend on the data that will be processed.

6.5.1 MRI based dose reconstruction

On starting an MRI based dose reconstruction workflow, Q-Suite 1.2 presents a workflow that guides the user through all required steps to convert MRI data into absorbed dose. For each step an indication is shown whether the step is completed or not. Workflow steps that depend on activities in a previous step will only be available if that previous step is completed.

Image data selection pre-treatment

In this step an MGE dataset acquired before administration of Holmium-166 microspheres must be selected (see section 6.2.1 for recommendations on the type of data). Only data suitable for dose reconstruction can be selected. After selection of the dataset the images are shown for inspection.

Creating contour sets pre-treatment

MRI-based dosimetry requires segmentation of Volumes of Interest (VOI's) prior to processing of the images. Two VOI's must be created: 1) 'Liver', including the entire liver and 2) 'Noise Measure' including a region in the images that is suitable for image noise measurements. See section 6.8 for a description of contouring functionality.

The accuracy of the output dose results does depend on the contour data. Factors that influence the accuracy of the dose calculations include:

- Including non-liver tissue in the 'Liver' contour or excluding liver tissue from the 'Liver' contour.
- Defining a 'Noise Measure' area that is not suitable or too small for measuring the image noise.

! Note ! : The quality of the input contour data directly influences the accuracy of the calculated dose.

Image data selection post-treatment

In this step an MGE dataset acquired after administration of Holmium-166 microspheres must be selected (see 6.2.1 for recommendations on the type of data). Only data suitable for dose reconstruction can be selected.

Creating contour sets post-treatment

As for the pre-treatment data also for the post-treatment data, a VOI 'Liver' and a VOI 'Noise Measure' must be segmented (see "Creating contour sets pre-treatment").

Reconstruction details

In addition to the image input data, the following treatment related input data is required:

- 1) the holmium content of the microspheres that were used for the treatment in percentage of weight and
- 2) the specific activity at the moment of injection of the microspheres that were used for treatment. These values must be entered in the appropriate fields.

The following reconstruction settings are available for MR based dose reconstructions:

- 1) Dose calculation method
- 2) Noise threshold

For the dose calculation method, the user can choose between applying a Dose Point Kernel to convert radiation activity into absorbed dose or applying the method Local Dose Deposition. The latter method assumes that all radiation energy is absorbed in a single voxel that contains the radiation source.

For the noise threshold, the user can choose between a default SNR of 3 or a value that can be selected from a list. The noise threshold is used to exclude signal intensities from the calculation that are too low. The default value of 3 is taken from scientific literature (Eur Radiol 2013;23:827-35).

Check and Confirm

In this last step an overview is provided of the data that will be used to reconstruct dose. After inspection of this overview the user can start the reconstruction.

During processing the user is informed about the progress. After successful processing Q-Suite 1.2 automatically starts the Dose Evaluation workflow.

6.5.2 SPECT based dose reconstruction

On starting a SPECT based dose reconstruction workflow, Q-Suite 1.2 presents a workflow that guides the user through all required steps to convert SPECT(-CT) data into absorbed dose. For each step there will be an indication whether the step is completed or not. Workflow steps that do depend on activities in a previous step will only be available if that previous step is completed.

Image data selection

In the first step a SPECT dataset acquired after administration of Holmium-166 microspheres must be selected (see section 6.2.1 for recommendations on the type of data). Only data suitable for dose reconstruction can be selected. The data can be a SPECT only image set or a SPECT-CT set.

Reconstruction details

In addition to the image input data, the following treatment and SPECT system related input data is required:

- 1) the date and time of injection of the microspheres used for treatment and
- 2) a calibration factor that allows conversion of measured photon counts into radiation activity concentration.

The date/time of injection must be entered in the appropriate field.

For the calibration factor two options are available: using a pre-defined SPECT system specific factor or determine the factor based on the current patient data.

A pre-defined calibration factor can be selected from a list of factors that have previously been stored in the Q-Suite 1.2 database (see section 6.11.3). Recommendations on how to determine this factor can be found in scientific literature [Elschot M, et al. Quantitative Monte Carlo-based holmium-166 SPECT reconstruction. Med Phys 2013;40:112502]. Prior to processing data, at least one calibration factor must be defined in the Q-Suite 1.2 database.

A patient specific calibration factor can be determined by:

1. Creating a VOI for calibration (see section 6.8)
2. Providing the radiation activity that is assumed to be within that VOI at time of injection

Q-Suite 1.2 will determine the calibration factor based on the total SPECT counts in the VOI, the assumed radiation activity and the decay time between injection and acquisition. The VOI itself shall include the entire volume where the radiation activity is assumed to be. This can for example be the liver only, the liver plus lungs or the entire body of the patient within the field of view.

Check and Confirm

In this last step an overview is provided of the data that will be used to reconstruct dose. After inspection of this overview the user can start the reconstruction.

During processing the user is informed about the progress. After successful processing Q-Suite 1.2 automatically starts the Dose Evaluation workflow.

6.6 DOSE EVALUATION

The dose evaluation activity consists of three steps:

- 1) creating dose review sets by combining dose data with image data
- 2) reviewing the dose by visual inspection and measuring dose values in VOI's
- 3) reporting the dose measurements

6.6.1 REVIEW SETS

A dose review set is a combination of a 3D dose distribution, reconstructed by Q-Suite 1.2, and a 3D image dataset. The purpose of a dose review set is to measure the radiation absorbed dose in tissues of interest that can be identified on the images.

Dose review set properties

Dose review sets have the following properties:

- Creation date: Date and time the dose review set was created
- Name: A user editable name

- Dose Series: Description of the dose data set including the type (SPECT or MR based) and reconstruction details
- Image series: Description of the image data set
- Status: Indicates the relation between the dose and image data. The status can have the following values:
 - o Dose with source image (see *Auto-generated sets*); The spatial relation between the dose and images is defined by the data itself and cannot be edited.
 - o Fused by data; the spatial relation between the dose and the image set is defined by the data itself (datasets share their frame of reference) and cannot be edited.
 - o Fused by user; the spatial relation between the dose and the image set has been set and confirmed by the user.
 - o Not Fused; The spatial relation between the dose and image set is unknown and has to be set and confirmed by the user (via edit fusion). Review sets with status 'not fused' cannot be used for dose review.

Auto-generated sets

Q-Suite 1.2 automatically generates the following dose review sets after a successful dose reconstruction:

- Dose combined with the input SPECT for a SPECT-based dose reconstruction. These review sets have the status 'dose with source image'.
- Dose combined with the CT of the input SPECT-CT for a SPECT-based dose reconstruction. These review sets have the status 'fused by data'.
- Dose combined with the first echo the input MGE MRI for a MRI-based dose reconstruction. These review sets have the status 'dose with source image'.

The auto-generated sets can be immediately used for dose review.

User defined sets

New sets can be created via 'Create set'. A window opens in which a dose set and an image set can be selected from the database. Supported image modalities for the image sets are CT and MRI. After selection of data and closing the window a new dose review set is created and the user can change the name of the set.

Edit fusion

If the spatial relation between the dose and the image set is not defined, i.e. the data is not fused, the user can set the relationship by applying a rigid registration. By selecting the dose review set to be fused and clicking 'edit fusion' the fusion screen opens. (See section 6.9). After confirming the fusion, the status of the dose review set will change to 'fused by user'.

6.6.2 DOSE REVIEW

In the dose review step the dose can be visually inspected and dose metrics can be performed on tissues of interest.

Dose visualization

For a selected dose review set the dose is displayed as a heat map overlaid on its image series in three orthogonal orientations. Scaling of the heat map can be controlled by left mouse click and dragging the mouse up or down on the color bar.

When moving the mouse over the dose display, the dose value of the voxel the mouse is located is displayed near the mouse cursor.

For controlling the image appearance see section 6.7

Dose Volume calculations

The dose review screen contains a dose volume table showing dose metrics for VOI's that are defined for the image set. For each VOI the following data is shown:

- Volume of the VOI
- The mean dose in the VOI
- The percentage of the VOI that received a certain dose range. This dose range is divided in 6 pre-defined bins.

VOI's can be added or edited via 'Edit VOI' which will open a dedicated contouring screen (see section 6.8)

6.6.3 REPORTING

Auto-generated data

The last step in the dose evaluation activity is creating a report for a certain dose review set. All relevant details of input images, dose reconstruction details and dose review sets are automatically copied to the dose report together with the dose volume table.

Image screen captures

Screenshot of dose review sets as well as screenshots of input data for the dose reconstruction can be added to the report. Screenshots of dose review sets can be created via the dose review screen by clicking on the camera icon in the image viewers. Screenshots of the dose reconstruction input data can be created via the dose reconstruction workflow that contains a read-only version of the dose reconstruction once a dose review set is selected.

User data

Each section of the report contains a area for user input. The first section contains a specific area to record the name of the user that created the report, all other section contain areas for general comments.

Saving report

When finished, the dose report can be saved as a PDF file to a local disk via 'Save to File'.

6.7 IMAGE VIEWERS

Several screens of Q-Suite 1.2 contain areas in which image data is shown. These image viewers all have standardized behavior and controls.

Navigation

Slice navigation: Navigation through slices is controlled by the mouse scroll wheel
Zoom: Zooming of images is controlled via right mouse click and dragging the mouse
Panning: Panning of images is controlled via left mouse click and dragging the mouse

Appearance

Each image viewer contains functionality to change the brightness and contrast of the image. By moving the mouse over the 'brightness/contrast' icon in the right top corner, sliders appear to change the levels. Image viewers that show dose overlayed on a background image also contain functionality to change the opacity of the overlayed dose. A slider to control that appears when moving the mouse over the brightness/contrast icon

6.8 CREATING VOLUMES OF INTEREST

Volumes of interest (VOI's) can be created for several purposes. For all these purposes Q-Suite 1.2 has a dedicated contouring screen in which VOI's can be created and deleted and contours can be drawn for these VOI's.

Image Views

The contouring screen contains two image viewers, a big one that shows the images in their original orientation (in general transversal) and a small one that shows one of the two orthogonal views at a time with the possibility to switch between the two. Contouring is only possible in the original orientation although the contours are shown in the other views.

Creating and editing VOI's

The contouring screens contain a list of VOI's that belong to the image data set that has been loaded. It depends on the purpose of creating contours if the list contains already predefined VOI's and if VOI's can be added or deleted. For MRI and SPECT based dose reconstruction only predefined VOI's are present that cannot be deleted. For dose review, the user can add VOI's up to a maximum of 16.

For each VOI:

- a name can be defined (double click on the name)
- a color can be chosen from a predefined list
- contours can be shown or hidden

User defined VOI's can all be deleted via the 'waste bin' button.

Contouring

On selection of a VOI, a ball-contour tool appears at the mouse cursor when the mouse is positioned in the original image view. Contouring can be controlled in the following way:

- Create initial contour on a slice:
Click and hold left mouse button while moving the mouse
- Edit existing contour:
Click and hold left mouse button while pushing against a contour from the inside or the outside
- Add a contour besides existing contour in the same image plan
Switch mode to 'add/cut' mode and click outside the existing contour
- Create a cut out in a existing contour
Switch mode to 'add/cut' mode and click inside a existing contour
- Delete contour on a single image plane
Click 'remove'
- Change size of ball contour tool
Click right mouse button and drag mouse

Navigating through the images during contouring can be done by the mouse scroll wheel. To zoom or pan, switch from 'Contour' to 'Zoom/Pan' mode. Temporarily switching from 'Contour' to 'Zoom/Pan' can be done by pressing and holding the Ctr key.

6.9 EDIT IMAGE FUSION

Editing an image fusion means displacing two 3D image datasets with respect to each other to align the objects they contain. To do this both datasets need to be visualized in a useful way and the user should be able to displace one of the two datasets.

Datasets for fusion

Fusion of dose with image data in Q-Suite 1.2 is never done via the dose data itself but always via a corresponding (linked) image set. This corresponding image set is:

- The CT data of the input SPECT-CT in case of SPECT based dose reconstruction
- The SPECT data of input SPECT in case of a SPECT-based dose reconstruction
- The 1st echo of the MGE MRI data in case of an MRI-based dose reconstruction

To distinguish both sets, the image dataset to be fused with the dose is identified as 'Static background' and the image data corresponding with the dose is identified as 'Dynamic overlay'.

Visualization

The image sets are displayed in three orthogonal orientations. The following options are available to change the appearance of these sets:

Static Background:

Brightness: changing the brightness of the background image

Contrast: changing the contrast of the background image

Dynamic Overlay

Brightness: changing the brightness of the overlay image

Contrast: changing the contrast of the overlay image

Opacity: changing the opacity (transparency) of the overlay from 0 to 100%

Coverage: changing how far the overlay image covers the background image

Color Map: changing the color appearance of the overlay image (grayscale or rainbow)

Pattern: changing the overlay pattern (cover the background from top to bottom or left to right)

Displacement

The overlay image can be displaced with respect to the background either by translation in a plane or rotation in a plane. This translation and rotation can be done in all three orientations.

For translation select the 'Translate' mode, click the left mouse button and drag the image around.

For rotation, select the 'Rotate' mode, click the left mouse button and move the mouse down to rotate clockwise or up to rotate counterclockwise

6.10 DATABASE MAINTENANCE

All user imported image data, dose reconstructions and dose review sets are stored in the Q-Suite 1.2 database. Deleting data can be done in the following ways:

- In the home screen (Patient data administration and selection) all data of a selected patient can be deleted
- In 'Select dose reconstruction composition', incomplete (interrupted) compositions can be deleted
- In 'Create and fuse dose review sets' user created dose review sets can be deleted.

6.11 SETTINGS

To enter the settings page an additional login code must be provided. The 'Settings' page contains three screens in tabulated form: 'General', 'Calculations' and 'SPECT calibration'.

6.11.1 General

Default Import folder

A default folder can be set for importing image data in the Q-Suite 1.2 database. This folder is automatically used for import on starting a new Q-Suite 1.2 session. This folder can be set by clicking 'Change' and selecting a specific folder of choice in the Explorer window that appears.

Interface Language

The user can change the Q-Suite 1.2 user interface by selecting one of the available languages. Changing the language requires a restart of Q-Suite 1.2.

Database limit

A limit can be set to the maximum size of the Q-Suite 1.2 database. Setting a limit that is lower than the data size already stored in the database is prevented.

6.11.2 Calculations

Q-Suite 1.2 uses certain calculation parameters to calculate dose values based on MRI and SPECT images. Some of these values can be changed by the user.

! Note ! : changing the calculation parameters will have impact on the accuracy of the dose calculations!

MRI

For MRI the Holmium relaxivity can be adjusted. The user can choose between a predefined value, based on measurements described in scientific literature [van de Maat GH et al. Eur Radiol 2013;23:827-35] or a user defined value based on own relaxivity measurements. Recommendations to measure the holmium relaxivity can be found in the referenced paper.

6.11.3 SPECT calibration

The user can store multiple SPECT calibration factors in Q-Suite 1.2 which can be selected during preparation of a dose reconstruction. Factors can be added, edited and removed. Additional information that can be stored for each factor includes: used system, used collimator and used reconstruction method. Only the calibration factor is mandatory, the other fields are for information only.

7. TROUBLESHOOTING AND MAINTENANCE

7.1 TROUBLESHOOTING

For reporting errors in the software, please contact your local Q-Suite representative.

7.2 MAINTENANCE

The Q-Suite 1.2 software does not require any preventative maintenance, such as log files to be checked or to be cleared, database maintenance, etc.

7.3 SUPPORT

The lifetime of Q-Suite 1.2 is 5 years. During this time, Quirem Medical B.V. will periodically provide patches to fix non-critical issues on the newest release. Critical updates (safety bugs) will be provided for all releases within their lifetime. Quirem Medical B.V. may decide to release a critical update to more or all releases in use in the field. In case of a safety incident that relates to the use of Q-Suite 1.2, Quirem Medical B.V. will provide adequate follow-up actions during the lifetime as required by law.

8. TECHNICAL DATA

8.1 ACCURACY

! Note ! : The accuracy of the dose values calculated by Q-Suite 1.2 always depends on the quality and the content of the input images and the way the input data is prepared. Recommendations for image acquisition and preparation of data are provided in section 6.2. !

8.1.1 SPECT based calculations

Dose calculations using input SPECT images are performed based on a direct translation of voxel intensities (counts) into absorbed dose using system and decay time specific factors. This means that the SPECT voxel intensities should represent the actual activity concentration in that specific volume. All factors that lead to inaccurate reconstruction of the activity concentration in the final SPECT image will therefore introduce deviations in the dose calculations. Factors known for their influence of image quality include patient movement during acquisition, detector and collimator response and used reconstruction techniques.

Dose deviations introduced by Q-Suite 1.2 in the calculation process of transforming voxel intensity information into absorbed dose have been measured based on computer generated digital phantom data representing body, liver and tumor volume with specified intensities for each of these volumes correlating with a certain activity concentration. Measurements have been performed for a range of activity concentrations in the tumor volume.

A deviation of less than 1% has been determined for a dose range of 0 to 1000 Gy for this digital phantom data. Any of the above mentioned factors that influence the accuracy can lead to a deviation that is larger than 1%.

8.1.2 MRI based calculations

Dose calculations using input MR images include a sequence of processing steps to convert image intensities into absorbed dose. The effect that the presence of the paramagnetic holmium-166 microspheres induces changes in the MR parameter R_2^* is utilized to determine the local concentration of microspheres. Measuring this change requires pre-treatment and post-treatment MGE datasets as described in section 6.2. Subsequently the local concentrations of microspheres are translated into activity and dose values utilizing the known microsphere properties.

Factors that directly influence the accuracy of MRI based dose calculations include the presence of magnetic field distortions not caused by the microspheres, the signal-to-noise level in the images, the presence of image artifacts induced by for example patient/organ movement and used reconstruction techniques.

The implemented R_2^* based calculation model is sensitive to three other effects:

- deviations in voxel R_2^* values from the mean value of the entire liver volume (including all other tissue) in the pre-treatment image set will induce dose deviations.
- changes in voxel R_2^* values between the pre and post-treatment case induced by for example changes in tissue characteristics will induce dose deviations.
- very high concentrations of microspheres will lead to signal decay too fast to be able to measure, leading to underestimation of the dose.

The accuracy of Q-Suite 1.2 for MR based dose calculations has been validated on computer generated datasets simulating real image data in respect to expected image intensities but free of noise and artifacts. For these datasets the following observations have been made:

- in case the voxels of interest do have pre-treatment R_2^* values equal to the mean value of the entire volume (including liver and tumor) and there are no tissue induced changes, errors are limited to less than 1% for a dose range of 0 to 1000 Gy when appropriate echo times are used.
- In case the voxels of interest do have pre-treatment R_2^* values that are lower or higher than the mean value of the entire volume (including liver and tumor), the dose is respectively under- or overestimated where the actual deviation depends on the specific activity of the microspheres. Dose errors of 9 Gy were observed for R_2^* deviations of 10 s^{-1} in case of a specific activity of 6.7 MBq/mg but can increase in case of larger R_2^* deviations and higher specific activity

Any of the factors that were described before that influence the MR based dose calculations can lead to local deviations larger than the values mentioned in this section.

8.2 SYSTEM CONFIGURATIONS

Q-Suite 1.2 configurations

There is only a single configuration of Q-Suite 1.2. Customization of Q-Suite 1.2 is not supported.

Required and recommended computer systems

Q-Suite 1.2 requires a x64 based computer running Microsoft Windows 7, 8.1 or 10. In addition, the Microsoft .NET version 4.5.2 framework should be installed.

The following hardware is recommended to ensure a good user experience:

- Dual-core CPU, 2.0GHz or higher
- 16 GB RAM
- SSD with 20 GB free disk space

Data size and performance

The responsiveness of Q-Suite 1.2 and the time needed to process data depends on the size of the image data to be processed and the computer system configuration like processor speed and available RAM. Q-Suite 1.2 accepts large image datasets, but the user shall be aware that performing dose reconstructions and evaluation on high resolution/large volume datasets may cause slow or delayed system response.


8.3 UNITS

The following quantities and units are used in Q-Suite 1.2:

Quantity	Unit
Holmium content	% (weight/weight)
Specific activity	MBq/mg
Dose	Gy
Volume	mL
Activity	MBq
Magnetic field strength	T
Energy	keV
Time duration	H, s, ms,
Relaxivity	$\text{s}^{-1} \cdot \text{ml} \cdot \text{mg}^{-1} \cdot \text{T}^{-1}$
Distance	mm

9. SYMBOLS AND NAMES

The following symbols are used in Q-Suite 1.2

-  Date of manufacture
-  LOT Number
-  Manufacturer
-  Copyright
-  Brightness/contrast
-  Patient
-  Complete
-  Attention
-  Search

-  Consult instructions for use
-  eIFU, weblink refers to: www.quirem.com/IFU
-  weblink
-  CE-mark that show European Conformity
-  Screenshot
-  Delete
-  Close
-  Expand
-  Collapse

Q-Suite is a trademark of Quirem Medical

Quirem Medical B.V.
Zutphenseweg 55
7418 AH Deventer
Niederlande
Website: www.quirem.com

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