Q-Suite™

Imaging Software

EN	-	Instructions for Use	3
DE	-	Gebrauchsanweisung	26
ES	-	Instrucciones de uso	49
FR	-	Mode d'emploi	72
IT	-	Istruzioni per l'uso	95
NL	-	Gebruiksaanwijzing	118
PL	-	Instrukcja użycia	141
PT	-	Instruções de utilização	164
TR	-	Kullanım Talimatları	188
SV	-	Bruksanvisning	211
LT	-	Naudojimo instrukcijos	234
CS	-	Pokyny k použití	257
DA	-	Bruksanvisning	280
SK	-	Inštrukcie na používanie	303
SL	-	Navodila za uporabo	326
HU	-	Használati útmutató	349
HR	-	Upute za uporabu	372
RU	-	Инструкция по применению	395
FI	-	Käyttöohjeet	418
NO	-	Bruksanvisningen	441

REF QS-Q002

Quirem Medical B.V.
Zutphenseweg 55
7418 AH Deventer
The Netherlands
www.quirem.com

1. INTRODUCTION

This Instructions For Use is a manual for using Q-Suite™ 2.1 in a safe and responsible manner. Q-Suite™ 2.1 is a software device developed and manufactured by Quirem Medical B.V. for pre-treatment planning and post-treatment evaluation of holmium-166 selective internal radiation therapy (SIRT). It is assumed that the user has enough knowledge of operating personal computers and the operating system Microsoft (MS) Windows 7 and/or 8.1 and/or 10 to be able to use the relevant functionality of the system.

This manual describes the system, the intended use, safety related issues and the daily operation and interpretation of the systems interface and response. In this manual, you can find important information about how to operate the Q-Suite™ 2.1 software. The manual helps you with the operation and the maintenance of the device, in a safe and responsible manner. This document should be studied carefully before attempting to use the system.

Observe all warnings and cautions throughout these instructions. Throughout this manual certain conventions are used. These are warnings, cautions and notes. They provide a means of prioritizing information to be brought to the attention of the user. They are given as follows:

! Warning ! A warning tells you that there is a risk of serious personal injury or death

! Precaution ! A precaution tells you that there is a risk of non-serious injury.

! Note ! A note gives more information

2. DESCRIPTION

Q-Suite™ 2.1 is a medical image processing software designed to support planning and verification of holmium-166 SIRT treatment.

For planning, Q-Suite™ 2.1 determines the required holmium-166 radioactivity needed to deliver a specified radiation dose to the whole liver or a part of the liver and determines the expected dose on the lungs, liver and tumor tissue.

Input for determining the required radioactivity is an MR or CT image dataset on which the user has to identify the liver tissue volume. Based on this liver volume and the dose calculation method recommended by the manufacturer of currently the only available commercial product for holmium-166 SIRT (QuiremSpheres, Quirem Medical, the Netherlands) the required activity at moment of treatment is calculated.

Input for determining the expected dose on the lungs, liver and tumor tissue is a SPECT-CT dataset of a SIRT scout dose like 99mTc-MAA or holmium-166 microspheres (QuiremScout, Quirem Medical, the Netherlands) and the previously determined required activity for therapy. The recommended activity at moment of treatment can be adjusted by the user to evaluate the impact on the expected therapeutic mean dose to the tissue.

For verification, Q-Suite™ 2.1 reconstructs a 3-dimensional dose distribution based on 3D SPECT-CT images or multi-slice MR images that were acquired during or after administration of the therapeutic holmium-166 microspheres. When combining these dose reconstructions with anatomical MR or CT images, dose values within tissues of interest can be measured and reported. Input for a dose reconstruction is a Multi Gradient Echo MRI data set or a SPECT dataset acquired and reconstructed using methods optimized for holmium-166 imaging. Combined with treatment related parameters like the specific activity and holmium content of the therapeutic microspheres or the time the microspheres were injected, Q-Suite™ 2.1 calculates the distribution of the radioactivity in the liver tissue on a voxel level and subsequently the voxel wise absorbed radiation dose. For dosimetry within regions of interest, additional MR or CT images can be combined with the calculated dose values, volumes of interest can be drawn and dose metrics can be performed on these volumes. The results can be stored in a report for further distribution or archiving purposes..

3. INTENDED USE

3.1 INTENDED USE

Q-Suite™ 2.1 is intended to support planning and evaluation of SIRT treatment with holmium-166 microspheres.

3.2 INTENDED USER

Q-Suite™ 2.1 is intended to be used by medically trained staff in a clinical environment. The anticipated users are radiologists, nuclear medicine physicians, radiation oncologists or other medical persons that are involved in SIRT treatment.

3.3 INDICATIONS FOR USE

Q-Suite™ 2.1 is indicated for patients that are eligible for SIRT treatment of liver tumors

3.4 INTENDED USE ENVIRONMENT

Q-Suite™ 2.1 is intended to be used in an environment with normal, non-stressful, office conditions.

4. COMPATIBLE DEVICES

Q-Suite™ 2.1 is a stand-alone software device. Q-Suite™ can exchange data with DICOM compliant devices as described in the Q-Suite™ 2.1 DICOM conformance statement (LC-80094).

5. SAFETY AND SECURITY

5.1 WARNINGS

Warnings are based on hazard analysis performed during the product life cycle. No hazard has been identified that requires a warning.

5.2 PRECAUTIONS

Precautions are based on hazard analysis performed during the product life cycle. The following hazards have been identified that requires a precaution

- Q-Suite™ calculates the required activity to achieve a dose on a Target Liver VOI based on the volume enclosed by the user defined contour for that VOI. Inaccurate contouring of the Target Liver VOI can therefore lead to an aimed dose that is higher or lower than intended
- A planned dose of 60 Gy on the whole liver volume is recommended by the QuiremSpheres IFU. To avoid calculation of activities that will result in exceeding a dose of 60 Gy on the whole liver, the user needs to contour a Whole Liver VOI that is used by Q-Suite™ to limit the aimed dose that can be defined by the user for the specific Target Liver VOIs. Inaccurate contouring of the Whole Liver VOI can therefore lead to incorrect application of this limitation and as consequence lead to a total calculated activity that will result in a dose on the whole live that is higher than 60 Gy.
- For accurate lung dose prediction based on SPECT-CT, the SPECT and CT images need to be aligned in sapce. Incorrect aligned may lead to incorrect lung dose prediction.
- Q-Suite™ predict the lung dose based on the volumes enclosed by the user defined contours for the 'Liver' and 'Lung' VOIs. Inaccurate contouring of these VOIs can therefore lead to incorrect lung dose prediction
- If the SPECT image does not cover the entire 'Liver' and/or 'Lung' VOIs, the predicted lung dose value may not be representative for the lung VOIs
- Only qualified staff may use the system. The system manager at the clinic must ensure that individuals authorized to use Q-Suite™ 2.1 have been suitably trained in using the system by Quirem staff or other competent people.
- It is important to ensure that the operator is aware that the quality of the output depends critically on the quality of the input data, and any irregularities or uncertainties about input data units, identification, or quality of any other nature shall be thoroughly investigated before the data is used. It is therefore highly recommended that a medical physics or imaging expert is involved throughout the introduction of Q-Suite™ 2.1 in the clinic.

- Data stored in the Q-Suite™ 2.1 database is encrypted and the use of Q-Suite™ 2.1 is protected by a login code to prevent unauthorized login. System settings are secured by an additional code. The hospital system manager must ensure that each individual has the correct access to Q-Suite™ 2.1. It is strongly recommended to give the login code only to people that are authorized to use Q-Suite™ 2.1.
- Data exchange with DICOM compliant device is not encrypted. The hospital system manager must ensure protection of the network against unauthorized access.

5.3 SAFETY AND SECURITY

Installation and upgrading of Q-Suite™ 2.1 shall be done by or under supervision of Quirem Medical B.V. or one of its official distributors.

Q-Suite™ 2.1 contains database functionality but it is not intended to be used as a long-term archive of patient data. The user is responsible for storage and backup of in and output data.

The use of Q-Suite™ 2.1 is protected through licenses.

It is strongly recommended to use anti-virus software for all systems Q-Suite™ 2.1 is installed on, whether or not connected to the Internet. Anti-virus software and firewalls can have negative impact on the functionality of Q-Suite™ 2.1 and should be configured correctly. The system manager should ensure that no custom settings in MS Windows© make information disappear or become illegible; for instance, when selecting color schemes for title bars, selecting font type and color, etc

! Precaution ! Data stored in the Q-Suite™ 2.1 database is encrypted and the use of Q-Suite™ 2.1 is protected by a login code to prevent unauthorized login. System settings are secured by an additional code. The hospital system manager must ensure that each individual has the correct access to Q-Suite™ 2.1. It is strongly recommended to give the login code only to people that are authorized to use Q-Suite™ 2.1.

Raw image pixel data without any patient identifiable details are temporarily stored for calculation purposes. It is recommended to restrict access to the computer system on which Q-Suite™ 2.1 is installed to authorized users only.

! Precaution ! Data exchange with DICOM compliant device is not encrypted. The hospital system manager must ensure protection of the network against unauthorized access.

6. INSTRUCTIONS FOR USE

6.1 INSTALLATION AND CONFIGURATION

6.1.1 Installation

Prior to installation, verify that the computer system where Q-Suite™ 2.1 will be installed on meets the requirements as specified in section 8.

- Install Q-Suite™ 2.1 using the Q-Suite™ 2.1 installation DVD. Installation can be started by running the file SuiteSetup.exe (run as administrator). The Q-Suite™ Setup window will be shown.
- Select 'Install' to start installation. All required components will automatically be installed.
- After successful installation, select 'Close' to exit the installation program.

After correct installation, Q-Suite™ 2.1 can be found at the following location: C:\Program Files\Quirem Medical\Q-Suite™ 2.1. The software can be started from the start menu.

6.1.2 De-installation

Q-Suite™ 2.1 can be removed from the computer system in the following way.

- Uninstall Q-Suite™ 2.1 using the Microsoft Windows functionality *Uninstall or Change a Program (Windows 7 and 8.1)* or *Uninstall Apps & Features (Windows 10)*.

6.1.3 Removing Patient Data

After uninstalling Q-Suite™ 2.1, patient data can be **permanently** removed in the following way:

- Delete the folder C:\ProgramData\Quirem Medical\QSuite-2.1-...* from the local hard drive.
- Delete the folder C:\QuiremMedical\CalculationJobs from the local hard drive.

*ProgramData is a hidden folder on Windows systems. To navigate to it, the option 'show hidden items' in windows explorer should be selected.

How to remove patient data without uninstalling Q-Suite™ 2.1 is described in section 6.10.

6.2 PREPARATION OF IMAGE DATA

All image data to be processed by Q-Suite™ 2.1 shall be stored in classic DICOM format on a local disk, mapped network drive or Picture Archiving and Communication System (PACS).

6.2.1 PREPARATION FOR ACTIVITY PLANNING

Q-Suite™ 2.1 supports activity planning based on single or multi-slice CT or MR data. The image quality of this data must be such that the whole liver and targeted liver volumes can be properly delineated based on visual appearance of these volumes. Guidelines for acquisition of these images are not part of this instruction.

6.2.2 PREPARATION FOR LUNG DOSE PREDICTION

Q-Suite™ supports lung dose prediction based on planar or SPECT-CT images of QuiremScout, a holmium-166 based surrogate marker or ^{99m}Tc-MAA, a technetium-99m based surrogate marker. For this reason, Q-Suite™ only accepts image data that includes the holmium-166 photopeak (80.6 keV) or the ^{99m}Tc photopeak (140.5 keV) in the acquisition energy window. The accuracy of the lung dose prediction will always depend on the image quality of the input planar or SPECT-CT images. Therefore it is strongly recommended to use a state-of-the-art SPECT-CT scanner and an advanced reconstruction method that includes scatter correction, attenuation correction and modelling of the collimator-detector response. For accurate lung dose prediction the image data should cover the whole liver and lungs. Guidelines for acquisition of these images are not part of this instruction.

! Precaution ! For accurate lung dose prediction based on SPECT-CT, the SPECT and CT images need to be aligned in space. Incorrect aligned may lead to incorrect lung dose prediction.

! Note ! The quality of the input Planar or SPECT-CT images directly influences the accuracy of the lung dose prediction

6.2.3 PREPARATION FOR LIVER, LUNG AND TUMOR DOSE SIMULATION

Q-Suite™ supports dose simulation to liver, lung and tumor tissue based on SPECT or SPECT-CT images of QuiremScout, a holmium-166 based surrogate marker or ^{99m}Tc-MAA, a technetium-99m based surrogate marker. For this reason, Q-Suite™ only accepts image data that includes the holmium-166 photopeak (80.6 keV) or the ^{99m}Tc photopeak (140.5 keV) in the acquisition energy window.

The accuracy of the dose simulation will always depend on the image quality of the input SPECT or SPECT-CT images. Therefore it is strongly recommended to use a state-of-the-art SPECT-CT scanner and an advanced reconstruction method that includes scatter correction, attenuation correction and modelling of the collimator-detector response. For accurate dose simulation the image data should cover the whole liver and lungs. Guidelines for acquisition of these images are not part of this instruction.

! Note ! The quality of the input SPECT or SPECT-CT images directly influences the accuracy of the dose simulation

For dose simulation, the SPECT or SPECT-CT data of the surrogate marker can be fused with MR, CT, SPECT-CT or PET-CT images for proper identification of tissues of interest (See also section 6.9). Guidelines for acquisition of these images are not part of this instruction.

6.2.4 PREPARATION FOR DOSE RECONSTRUCTION

Q-Suite™ 2.1 can reconstruct radiation absorbed dose values based on input SPECT data or MR data.

Acquisition recommendations MRI

MRI-based holmium-166 dose reconstruction is based on measuring the effect of the presence of the microspheres on MR signal decay. This requires a multi gradient echo (MGE) read-out of the free induction decay (FID) using at least 2 echoes. Data with less than 2 echoes is not compatible with Q-Suite™ dose reconstruction. Two MGE datasets should be acquired, one prior to and one after microsphere delivery. The images should include the entire volume for which dose needs to be reconstructed and there should be no overlap or gaps between image slices.

The accuracy of the output dose results of Q-Suite™ 2.1 will always depend on the image quality of the input MR images. To maximize the accuracy, acquisition settings should be optimized. Factors that influence the accuracy of the dose calculations include:

- the signal to noise ratio (SNR) of the images (which depends on the MR hardware, acquisition times, resolution, acceleration techniques etc). The higher the SNR, the more accurate the dose reconstruction will be
- the number of gradient echoes used to sample the FID (minimum of 4 echoes recommended). The more echoes that are available as sample point, the more accurate the dose reconstruction will be
- the echo times of the gradient echoes (fast signal decay requires fast sampling)
- the presence of artifacts in the images (e.g. ghosting, fold-over, object induced), which in general lead to decreased dose reconstruction accuracy
- movement during acquisition (patient and/or organ), which in general leads to decreased dose reconstruction accuracy

For more background information on the subject of image acquisition see a paper published by van de Maat et al (Eur Radiol 2013;23:827-35).

! Note ! The quality of the input MR images directly influences the accuracy of the calculated dose.

Acquisition recommendations SPECT

For SPECT-based dose reconstructions, a 3D reconstructed SPECT image is required with an acquisition energy window that includes the holmium-166 photopeak of 80.6 keV. The SPECT image can be accompanied by a multi-slice CT dataset in case of a SPECT-CT acquisition. The accuracy of the output dose results of Q-Suite™ 2.1 will always depend on the image quality of the input SPECT images. Therefore it is strongly recommended to use a state-of-the-art SPECT-CT scanner and an advanced reconstruction method that includes scatter correction, attenuation correction and modelling the collimator-detector response.

! Note ! The quality of the input SPECT images directly influences the accuracy of the calculated dose.

! Note ! For accurate dose reconstruction based on SPECT-CT using Patient specific calibration, the SPECT and CT images need to be aligned in space. Incorrect aligned may lead to incorrect dose reconstruction.

6.2.5 PREPARATION FOR DOSE EVALUATION

For the purpose of dose evaluation, the Q-Suite™ 2.1 generated dose reconstructions can be fused with MR, CT, SPECT-CT or PET-CT images for proper identification of tissues of interest (See also section 6.9). Guidelines for acquisition of these images are not part of this instruction.

6.3 GENERAL WORKING PRINCIPLES

Three main functional areas can be distinguished in Q-Suite™ 2.1:

- Patient data administration and selection (home screen)
- Treatment planning including activity planning, dose prediction and dose simulation
- Treatment evaluation including dose reconstruction and dose evaluation

On starting Q-Suite™ 2.1, a home screen is displayed that shows all patients of which image data is stored in the Q-Suite™ 2.1 database. On selection of a patient, the user can choose to start the main activity treatment planning or the main activity treatment evaluation for that patient, depending on the data that is available. On entering one of these main activities for a selected patient, that patient becomes the active patient. There can be only one active patient at a time. By closing the activity for the active patient the user returns to the home screen.

6.4 DATA IMPORT

Via the home screen, new patient data can be imported to the Q-Suite™ 2.1 database. Data can be imported from a local folder, a PACS system or the Q-Suite™ Inbox. Only Q-Suite™ 2.1 compatible data as documented in the Q-Suite™ 2.1 DICOM Conformance Statement (LC-80094) will be imported. The software will recognize if data selected for import is already present in the database and will not import that data again.

6.4.1 import from local folder

For importing data from a local folder or a folder on a mapped network location, the user can select a folder to search for data. DICOM data in that folder will be recognized by Q-Suite™ 2.1 and presented to the user for selection for import.

6.4.2 IMPORT FROM PACS

For importing data from PACS, the user can select a PACS system that has been configured to search for data (see 6.11.4 for configuration of a PACS). The user can search the PACS based on patient ID or patient name and the patients that have been found are presented by the user for selection for import.

6.4.3 IMPORT FROM INBOX

Data that has been exported from a PACS to Q-Suite™ is stored in the Q-Suite™ inbox. The data found in the box will be presented to the user for selection of import

6.5 TREATMENT PLANNING

The treatment planning activity provides functionality to determine required QuiremSpheres activity for a treatment and predict the dose to the lungs for such a treatment according to instructions for use recommended by the manufacturer of QuiremSpheres. Besides that, a dose simulation can be performed to assess the expected dose to healthy liver tissue and tumor tissue based on a user defined projected holmium-166 activity. This information can add in getting better insights in the effect a certain activity of QuiremSpheres has in terms of expected dose to the tissue.

! Note ! QuiremSpheres dose simulation is meant for simulation purposes only and should not be used to plan treatment other than recommended by the IFU of QuiremSpheres

6.5.1 Activity Planning and lung dose prediction

The QuiremSpheres activity planning consist of three steps:

- 1) Activity planning
- 2) Lung dose prediction (not mandatory)
- 3) Reporting

Activity planning

In this step, first a CT or MR dataset must be selected to be used for the planning. Only 3D (multi slice) data can be selected. After selection the images are shown and can be inspected. Next, 'Target Liver' and 'Whole Liver' Volumes of Interest (VOIs) have to be delineated by the user. See section 6.8 for a description of contouring functionality. Q-Suite™ 2.1 will automatically create a VOI of type Non-target Liver based on subtraction of Whole Liver and Target Liver VOIs or a VOI of type Whole Liver based on summation of Target Liver and Non-target Liver VOIs. If neither Whole Liver nor Non-target Liver VOIs are defined by the user, the user can choose to let Q-Suite™ 2.1 automatically create a Whole Liver VOI based on summation of the Target Liver VOIs.

For the 'Target Liver' VOIs, the user has to define the aimed dose. Based on this user defined aimed dose and the volume enclosed by the corresponding user defined contour for that VOI, Q-Suite™ 2.1 calculates the QuiremSpheres activity that is needed to obtain the aimed dose in that specific segment. Inaccurate contouring of the 'Target Liver VOIs' will therefore lead to inaccurate treatment activity planning. Excluding liver volume from the contour that should be included in the VOI will lead to a lower activity than needed and including liver volume in the contour that should be excluded from the VOI will lead to a higher activity then needed.

Q-Suite™ limits the aimed dose that can be defined on the 'Target Liver VOIs' to not exceed a mean dose of 60 Gy on the Whole Liver VOI. To accurately apply this limitation, the whole liver needs to be accurately contoured. Excluding liver volume from the contour that should be included in the Whole Liver VOI will lead to a lower dose limit, including liver volume in the contour that should be excluded from the Whole Liver VOI will lead to a higher dose limit.

! Precaution ! Q-Suite™ calculates the required activity to achieve a dose on a Target Liver VOI based on the volume enclosed by the user defined contour for that VOI. Inaccurate contouring of the Target Liver VOI can therefore lead to an aimed dose that is higher or lower than intended

! Precaution ! A planned dose of 60 Gy on the whole liver volume is recommended by the QuiremSpheres IFU. Q-Suite™ limits the dose that can be planned on the Target Liver VOIs to not exceed 60 Gy on the Whole Liver VOI. Inaccurate contouring of the Whole Liver VOI can therefore lead to aimed doses that are too high

Lung Dose Prediction

For predicting the lung dose for the activity planned in step 1 the user can use one of the following options:

- a) Volumetric SPECT-CT based. For this method a QuiremScout or ^{99m}Tc-MAA SPECT-CT dataset must be selected on which 'Liver' and 'Lung' VOIs must be created using the CT image for contouring. For accurate lung dose prediction, the SPECT image should cover the entire VOIs. Based on the SPECT counts measured in these VOIs, Q-Suite™ will calculate the

lung shunt and subsequently the predicted mean dose to the lung tissue. This dose value will be displayed in the corresponding table. The user can exclude a caudal part of the lungs from analyses to decrease influence of photon scatter and/or liver movement on lung dose prediction. The part to exclude can be defined by choosing a distance from the liver VOI border.

- b) 2D Planar Scintigraphy based. For this method a planar QuiremScout or ^{99m}Tc-MAA image must be selected on which 'Liver' and 'Lung' Regions of Interest (ROIs) must be created. Based on the counts measured in these two ROIs, Q-Suite™ will calculate the lung shunt. To convert the lung shunt into mean dose to the lungs the user can either provide a pre-measured lung volume or choose to use a standard lung mass of 1 kg.
- c) Lung shunt fraction based. For this method the user must provide a lung shunt fraction obtained from previous measurements. To convert this lung shunt into mean dose to the lungs the user can either provide a pre-measured lung volume or choose to use a standard lung mass of 1 kg.

! Precaution ! Q-Suite™ predict the lung dose based on the volumes enclosed by the user defined contours for the 'Liver' and 'Lung' VOIs. Inaccurate contouring of these VOIs can therefore lead to incorrect lung dose prediction

! Precaution ! If the SPECT image does not cover the entire 'Liver' and/or 'Lung' VOIs, the predicted lung dose value may not be representative for the lung VOIs

! Note ! The standard lung mass of 1 kg is recommended by the QuiremSpheres IFU. When using lung volume for dose calculations a density of 0.3Kg/L is assumed.

Reporting

Results from activity planning and lung dose prediction can be stored in a structured report. Q-Suite™ automatically copies all the calculated activity and predicted lung dose values to the report together with all relevant details of input images and user defined input.

Screenshots of image data for both activity planning and lung dose prediction can be added to the report via their specific screens by clicking on the camera icon in the image viewers.

The report contains several areas to record user input for example the name of the user that created the report or general comments.

When finished, the dose report can be saved as a PDF file to a local disk or, if configured, as DICOM encapsulated PDF to a PACS system. See section 6.11.4 for configuration of a PACS connection.

6.5.2 Liver, lung and tumor dose simulation

The dose simulation activity consists of three steps:

- 1) Creating dose simulation sets
- 2) Simulating dose for a dose simulation set
- 3) Creating a dose simulation report

Dose simulation sets

In this step, dose simulation sets can be created. A dose simulation set is a combination of a QuiremScout or ^{99m}Tc-MAA SPECT or SPECT/CT set and a 3D image dataset that is used to identify VOIs based on morphology or function (VOI image series). The purpose of a dose simulation set is to simulate the radiation absorbed dose in tissues of interest that can be identified on the images.

Dose simulation set properties

Dose simulation sets have the following properties:

- Creation date: Date and time the dose review set was created
- Name: A user editable name
- SPECT Series: Description of the SPECT data set
- VOI Image series: Description of the VOI image data set
- Status: Indication of the relation between the SPECT/CT and VOI image series. The status can have the following values:
 - o Standalone; the SPECT/CT is combined with its own CT as VOI image series. The spatial relation is fixed and cannot be edited.
 - o Fused by data; the spatial relation between the SPECT/CT and the VOI image series is defined by the data itself (datasets share their frame of reference) but can be changed by the user

- Fused by user; the spatial relation between the SPECT/CT and the VOI image series has been set and confirmed by the user.
- Not fused; the spatial relation between the SPECT/CT and VOI image series is unknown and has to be set and confirmed by the user. Simulation sets with status 'Not fused' cannot be used for dose simulation.

Create set

New dose simulation sets can be created by the user via a dedicated selection window in which first a SPECT-CT set can be selected from the database followed by selection of a VOI image series. Supported image modalities for the VOI image series are CT, MRI, SPECT/CT and PET/CT. If the CT of the SPECT/CT is selected as VOI image series the dose simulation set will get the status 'Standalone'. After creating the set the user can edit the name of the set.

Edit fusion

If the spatial relation between the SPECT-CT and the VOI image series is not defined, i.e. the data is not fused, the user can set the relationship by applying a rigid registration. By selecting the dose simulation set to be fused and clicking 'edit fusion' the fusion screen opens. (See section 6.9). After confirming the fusion, the status of the dose review set will change to 'Fused by user'. Fusion editing can also be performed for dose simulation sets that have status "fused by data" or "fused by user" to change the existing spatial relation between the SPECT-CT and the VOI image series. When this is done for a dose simulation set with status "fused by data", the status will change to "fused by user" after confirming the fusion.

Dose simulation

In the dose simulation step, dose simulation sets can be visually inspected and dose simulations can be performed on tissues of interest based on user defined treatment activity.

For a selected dose simulation set, the SPECT image is displayed as a color map overlaid on the VOI image series in three orthogonal orientations. For controlling the image appearance see section 6.7.

To perform a dose simulation, first VOIs have to be defined for which the simulation must be performed. VOIs can be added or edited via a dedicated contouring screen (see section 6.8). Q-Suite™ 2.1 will automatically create a VOI of type Non-target Liver based on subtraction of Whole Liver and Target Liver VOIs or a VOI of type Whole Liver based on summation of Target Liver and Non-target Liver VOIs. If neither Whole Liver nor Non-target Liver VOIs are defined by the user, the user can choose to let Q-Suite™ 2.1 automatically create a Whole Liver VOI based on summation of the Target Liver VOIs. For Target Liver and Whole Liver VOIs that include a (part of) Healthy tissue VOI, Q-Suite™ 2.1 will automatically create an additional VOI of type Healthy Liver based on intersection of those VOIs with the Healthy tissue VOI. If a Healthy tissue VOI is not defined but (a) Tumor VOI(s) are defined, the user can choose to let Q-Suite™ 2.1 automatically create this additional of type 'Healthy Liver' based on subtraction of the Tumor and Liver VOIs.)

For each Target Liver VOI, the user can define the treatment activity to be used for the simulation. When at least one Target Liver VOI with a corresponding treatment activity has been defined for each VOI the following data is shown:

- Volume
- Mean dose
- Dose Volume Table (DVT)
- Dose Volume Histogram (DVH)

The maximum dose range of the DVH and the number of bins and the dose range per bin of the DVT can be set by the user with a maximum of 5 bins.

When a healthy liver VOI is present, the option becomes available to limit the maximum treatment activity that can be defined based on a user defined mean dose limit to that healthy liver tissue. When a Whole Liver VOI is present, the option becomes available to limit the maximum treatment activity that can be defined based on a user defined mean dose limit to that Whole Liver. When a lung VOI is present, the user can activate the option to account for lung shunt and the treatment activity can be limited by defining a mean dose limit for the lung tissue. Limiting the mean dose can still lead to a local dose that is higher than the set limit.

! Note ! If the SPECT image does not cover the entire VOIs used for the simulation, the simulated dose values may be incorrect or not representative for the entire VOI

! Note ! QuiremSpheres dose simulation is meant for simulation purposes only and should not be used to plan treatment other than recommended by the IFU of QuiremSpheres

Reporting

Results from dose simulation can be stored in a structured report. Q-Suite™ automatically copies the dose volume histogram and the values from the dose volume table to the report together with all relevant details of input images and user defined input.

Screenshots of the dose simulation set can be added to the report by clicking on the camera icon in the image viewers in the respective screen.

The report contains several areas to record user input for example the name of the user that created the report or general comments.

When finished, the dose report can be saved as a PDF file to a local disk or, if configured, as DICOM encapsulated PDF to a PACS system. See section 6.11.4 for configuration of a PACS connection.

6.6 TREATMENT EVALUATION

The treatment evaluation activity provides functionality to reconstruct dose from SPECT data or MR data and to perform dose evaluation based on such a dose reconstruction.

On starting a dose reconstruction for a selected patient, the user has three options: starting a new MRI based dose reconstruction, starting a new SPECT based dose reconstruction or loading a previously stored dose reconstruction composition from the database. The available options depend on the data that is stored in the database.

Each dose reconstruction workflow consists of specific steps that depend on the data that will be processed.

6.6.1 SPECT based dose reconstruction

On starting a SPECT based dose reconstruction workflow, Q-Suite™ 2.1 presents a workflow that guides the user through all required steps to convert SPECT(-CT) data into absorbed dose. For each step there will be an indication whether the step is completed or not. Workflow steps that do depend on activities in a previous step will only be available if that previous step is completed.

Image data selection

In the first step a SPECT dataset acquired after administration of Holmium-166 microspheres must be selected (see section 6.2.1 for recommendations on the type of data). Only data suitable for dose reconstruction can be selected. The data can be a SPECT only image set or a SPECT-CT set.

Reconstruction details

To convert the SPECT images (counts) into quantitative activity images (MBq) two methods are available:

- Patient specific calibration (default).
- SPECT system specific calibration using a pre-measured calibration factor.

Patient specific calibration

For the Patient specific calibration method the user must create a VOI/VOIs for calibration (see section 6.8) and provide the radiation activity that is assumed to be within that VOI/those VOIs at time of injection.

Q-Suite™ 2.1 will determine the calibration factor at treatment time based on the total SPECT counts in the VOI and the assumed radiation activity. The VOI itself shall include the entire volume where the radiation activity is assumed to be. This can for example be the liver only, the liver plus lungs or the entire body of the patient within the field of view.

! Note ! If the SPECT image does not cover the entire VOIs used for the calibration factor, the calibration factor may be incorrect which may lead to incorrect dose reconstruction

Use pre-measured factor

For this method the user must enter the date and time of injection of the microspheres used for treatment and select a pre-measured calibration factor. The factor can be selected from a list of factors that have previously been stored in the Q-Suite™ 2.1 database (see section 6.11.3).

Q-Suite™ 2.1 will determine the calibration factor at treatment time based on the decay time between injection and acquisition. Recommendations on how to measure a system specific calibration factor can be found in scientific literature [Elschot M, et al. Quantitative Monte Carlo-based holmium-166 SPECT reconstruction. Med Phys 2013;40:112502]. Prior to processing data, at least one calibration factor must be defined in the Q-Suite™ 2.1 database.

Check and Confirm

In this last step an overview is provided of the data that will be used to reconstruct dose. After inspection of this overview the user can start the reconstruction.

During processing the user is informed about the progress. After successful processing Q-Suite™ 2.1 automatically starts the Dose Evaluation workflow.

6.6.2 MRI based dose reconstruction

On starting an MRI based dose reconstruction workflow, Q-Suite™ 2.1 presents a workflow that guides the user through all required steps to convert MRI data into absorbed dose. For each step an indication is shown whether the step is completed or not. Workflow steps that depend on activities in a previous step will only be available if that previous step is completed.

Image data selection pre-treatment

In this step an MGE dataset acquired before administration of Holmium-166 microspheres must be selected (see section 6.2.1 for recommendations on the type of data). Only data suitable for dose reconstruction can be selected. After selection of the dataset the images are shown for inspection.

Creating contour sets pre-treatment

MRI-based dosimetry requires segmentation of VOIs prior to processing of the images. Two VOIs must be created: 1) a VOI of type 'Quantification', including the entire liver and 2) a VOI of type 'Noise Measure' including a region in the images that is suitable for image noise measurements.

See section 6.8 for a description of contouring functionality.

The accuracy of the output dose results does depend on the contour data. Factors that influence the accuracy of the dose calculations include:

- Including non-liver tissue in the 'Quantification' contour or excluding liver tissue from the 'Quantification' contour.
- Defining a 'Noise Measure' area that is not suitable or too small for measuring the image noise.

! Note ! The quality of the input contour data directly influences the accuracy of the calculated dose.

Image data selection post-treatment

In this step an MGE dataset acquired after administration of Holmium-166 microspheres must be selected (see 6.2.1 for recommendations on the type of data). Only data suitable for dose reconstruction can be selected.

Creating contour sets post-treatment

As for the pre-treatment data also for the post-treatment data, a VOI 'Quantification' and a VOI 'Noise Measure' must be contoured (see "Creating contour sets pre-treatment").

Reconstruction details

In addition to the image input data, the following treatment related input data is required:

- 1) the holmium content of the microspheres that were used for the treatment in percentage of weight and
- 2) the specific activity at the moment of injection of the microspheres that were used for treatment.

These values must be entered in the appropriate fields.

For the MR signal fit algorithm several control parameters are available:

- SNR threshold R_2^* fit: The Signal to Noise (SNR) threshold is used to exclude signal intensities that are too low from the calculation. The user can choose a default value of 3 or select a value from a list. The default value of 3 is taken from scientific literature (Eur Radiol 2013;23:827-35).
- Minimum echoes R_2^* fit: This number defines the minimum number of gradient echoes that should be available for a voxel to perform R_2^* fitting after applying the SNR threshold. The user can choose a default value of 2 or select a value from a list. The default value is the minimum to mathematically perform the calculation
- Range neighbors S_0 -fit: for S_0 -fitting, S_0 -values of neighboring voxels are used to estimate the S_0 value of a voxel for which the R_2^* fitting failed because of low signal and requires S_0 -fitting. Range neighbors is the maximum distance in nr of voxels between the neighboring voxel that can be used for S_0 -fitting and the voxel that requires S_0 -fitting. The user can choose a default value of 3 select a value from a list.
- Minimum neighbors S_0 -fit: This number defined the minimum number of neighbor S_0 values required to estimate the S_0 -value of the voxel that requires S_0 -fitting. The user can choose a default value of 9 or select a value from a list.

See section 8.1 for a detailed description of the fit algorithm.

To convert voxel activity values into absorbed dose the user can choose between applying a Dose Point Kernel to convert radiation activity into absorbed dose or applying the method Local Dose Deposition. The latter method assumes that all radiation energy is absorbed in a single voxel that contains the radiation source.

Check and Confirm

In this last step an overview is provided of the data that will be used to reconstruct dose. After inspection of this overview the user can start the reconstruction.

During processing the user is informed about the progress. After successful processing Q-Suite™ 2.1 automatically starts the Dose Evaluation workflow.

6.6.3 DOSE EVALUATION

The dose evaluation activity consists of three steps:

- 1) Creating dose review sets
- 2) Reviewing the dose of a dose simulation set
- 3) Creating a dose review report

Dose review sets

In this step dose review sets can be created. A dose review set is a combination of a 3D dose distribution, reconstructed by Q-Suite™ 2.1, and a 3D image dataset that is used to identify VOIs based on morphological (VOI image series). The purpose of a dose review set is to measure the radiation absorbed dose in tissues of interest that can be identified on the images.

Dose review set properties

Dose review sets have the following properties:

- Creation date: Date and time the dose review set was created
- Name: A user editable name
- Dose Series: Description of the dose data set including the type (SPECT or MRI based) and reconstruction details
- VOI Image series: Description of the VOI image data set
- Status: Indication of the relation between the Dose and VOI image series. The status can have the following values:
 - o Standalone; the Dose Series is combined with the CT of the original SPECT-CT or the original MR data. The spatial relation is defined by the data itself (datasets share their frame of reference) and cannot be edited.
 - o Fused by data; the spatial relation between the Dose and the VOI image series is defined by the data itself (datasets share their frame of reference) and cannot be edited.
 - o Fused by user; the spatial relation between the Dose and VOI image series has been set and confirmed by the user.
 - o Not fused; The spatial relation between the Dose and VOI image series is unknown and has to be set and confirmed by the user. Review sets with status 'Not fused' cannot be used for dose review.

Auto-generated sets

Q-Suite™ 2.1 automatically generates the following (standalone) dose review sets after a successful dose reconstruction:

- Dose combined with the CT of the input SPECT-CT for a SPECT-based dose reconstruction.
- Dose combined with the input MGE MRI for an MRI-based dose reconstruction.

User defined sets

New dose review sets can be created by the user via a dedicated data selection window in which a dose series and a VOI image series can be selected from the database. Supported image modalities for the VOI image series are CT, MRI, SPECT/CT and PET/CT. If the input CT or MRI of the dose reconstruction is selected as VOI image series the dose simulation set will get the status 'Standalone'. After creating the set the user can edit the name of the set.

Edit fusion

If the spatial relation between the dose series and the VOI image series is not defined, i.e., the data is not fused, the user can set the relationship by applying a rigid registration. For editing the fusion of a dose review set, the dedicated fusion screen opens (See section 6.9). After confirming the fusion, the status of the dose review set will change to 'Fused by user'.

Dose review

In the dose review step the dose can be visually inspected and dose metrics can be performed on tissues of interest.

Dose visualization

For a selected dose review set the dose is displayed as a heat map overlaid on its VOI image series in three orthogonal orientations. When moving the mouse over the dose display, the dose value of the voxel the mouse is located is displayed near the mouse cursor.

For controlling the image appearance see section 6.7

Dose Volume calculations

VOIs for dose volume calculations can be added or edited via a dedicated contouring screen (see section 6.8). Q-Suite™ 2.1 will automatically create a VOI of type Non-target Liver based on subtraction of Whole Liver and Target Liver VOIs or a VOI of type Whole Liver based on summation of Target Liver and Non-target Liver VOIs. If neither Whole Liver nor Non-target Liver VOIs are defined by the user, the user can choose to let Q-Suite™ 2.1 automatically create a Whole Liver VOI based on summation of the Target Liver VOIs. For Target Liver and Whole Liver VOIs that include a (part of) Healthy tissue VOI, Q-Suite™ 2.1 will automatically create an additional VOI of type Healthy Liver based on intersection of those VOIs with the Healthy tissue VOI. If a Healthy tissue VOI is not defined but (a) Tumor VOI(s) are defined, the user can choose to let Q-Suite™ 2.1 automatically create this additional of type 'Healthy Liver' based on subtraction of the Tumor and Liver VOIs.)

For each VOI the following data is shown:

- Volume
- Mean dose
- Dose Volume Table (DVT)
- Dose Volume Histogram (DVH)

The maximum dose range of the DVH and the number of bins and the dose range per bin of the DVT can be set by the user with a maximum of 5 bins.

! Note ! If the dose series does not cover the entire VOI used for the evaluation, the resulting dose value may be incorrect or not representative for the entire VOI

Reporting

Results from dose review can be stored in a structured report. Q-Suite™ automatically copies the dose volume histogram and the values from the dose volume table to the report together with all relevant details of input images and dose.

Screenshots of the dose review set can be added to the report by clicking on the camera icon in the image viewers in the respective screen. Screenshots of the dose reconstruction input data can be created via the dose reconstruction workflow that contains a read-only version of the dose reconstruction once a dose review set is selected.

The report contains several areas to record user input for example the name of the user that created the report or general comments.

When finished, the dose report can be saved as a PDF file to a local disk or, if configured, as DICOM encapsulated PDF to a PACS system. See section 6.11.4 for configuration of a PACS connection.

6.7 IMAGE VIEWERS

Several screens of Q-Suite™ 2.1 contain areas in which image data is shown. These image viewers all have standardized behavior and controls.

Navigation

Slice navigation: Navigation through slices is controlled by the mouse scroll wheel, a scroll bar and keyboard up/down arrow keys

Zoom: Zooming of images is controlled via right mouse click and dragging the mouse

Panning: Panning of images is controlled via left mouse click and dragging the mouse

Re-center: The zoom and pan status of the images can be restored to the original state via the re-center icon (see section 9)

Appearance

Each image viewer contains functionality to change the brightness and contrast of the (background) image and, if applicable functionality to change the color scaling and opacity of the overlay image, to show and choose a color for the border of the overlay image and to change the thickness of the contour lines.

Brightness/Contrast:	Brightness and contrast levels of the background image can be changed via sliders that can be accessed by moving the mouse over the 'brightness/contrast' icon (see section 9)
Overlay Opacity:	The opacity of the overlay image can be changed via sliders that can be accessed by moving the mouse over the 'brightness/contrast' icon (see section 9)
Color scale:	Scaling of the overlay image can be controlled by left mouse click and dragging the mouse up or down on the corresponding color bar
Show overlay border:	The part of the background image that is outside the Field of View of the overlay image can be indicated by a color that can be chosen by the user to check if the Field of View of the overlay image covers the Field of View of the background image
Reset:	All appearance settings can be reset to the default values via the corresponding button
Stroke Thickness:	The stroke thickness of the displayed contour lines can be changed via a slider that can be accessed by moving the mouse over the 'brightness/contrast' icon (see section 9)

6.8 CREATING VOLUMES OF INTEREST

Volumes of interest (VOIs) can be created for several purposes. For all these purposes Q-Suite™ 2.1 has a dedicated contouring screen in which VOIs can be created and deleted and contours can be drawn for these VOIs.

Image Views

The contouring screen contains two image viewers, one that shows the images in their original orientation (in general transverse) and one that can show the orthogonal views of the contour data or helper images (see *helper images*). Contouring is only possible in the original orientation. In case the image is a SPECT-CT or PET-CT the user can choose whether the SPECT or PET or the CT image is displayed and/or whether the SPECT or PET is shown as an overlay image.

Creating and editing VOIs

The contouring screens contain a list of VOIs that belong to the image data set that has been loaded. It depends on the purpose of creating contours what type of VOIs and how many VOIs can be defined. The table below list the available VOI types per purpose.

Purpose	VOI types	Amount
Activity planning	Target Liver	Maximum 3
	Whole Liver	1
	Non-target Liver	No limit
Lung dose prediction	Lung	Maximum 2
	Liver	1
Dose Simulation	Target Liver	Maximum 3
	Whole Liver	1
	Non-target Liver	No limit
	Healthy Tissue	1
	Tumor	No limit
MR based dose reconstruction	Lung	Maximum 2
	Quantification	1
SPECT based dose reconstruction	Noise Measure	1
	Calibration	No limit
Dose review	Target Liver	Maximum 3
	Whole Liver	1
	Healthy Tissue	1
	Tumor	No limit

For each VOI:

- a name can be defined (double click on the name)
- a color can be chosen from a list

VOIs can be deleted via the 'waste bin' icon. (see section 9)

VOI relations and dependencies

The following relational dependencies do apply to VOIs

- VOIs of the same type cannot overlap each other
- Target Liver type cannot overlap Non-target Liver type
- Target Liver type cannot exist outside Whole Liver type
- Non-target Liver type cannot exist outside Whole Liver type
- Lung cannot overlap Liver
- Target Liver type cannot overlap Non-target Liver type
- Target Liver type cannot overlap Lung
- Target Liver type cannot exist outside Whole Liver type
- Non-target Liver type cannot exist outside Whole Liver type
- Non-target Liver type cannot overlap Lung
- Lung cannot overlap Whole Liver
- Tumor cannot overlap Healthy Tissue

Screen lay-out

Two different screen lay-outs are available for contour editing, a lay-out with equally sized image viewers and a lay-out with a big contour image viewer and a small side viewer. When the lay-out with equally sized image viewers is active, the list of VOIs can be expanded and collapsed to reserve maximum space for the image viewers.

Helper image

To help in identifying tissue, additional images series can be loaded and selected for display in the side image viewer. Multiple 3D image series (CT, MR, SPECT-CT or PET-CT) can be added to the selection list. Slice navigation through the contour image series and helper image can be synchronized by navigating to corresponding image slices in both viewers and linking those slices via the corresponding icon (see section 9). When linked, on navigating in one of the viewers the other viewer will automatically display the corresponding slice of the other set based on slice position information.

Manual Contouring

On selection of a VOI, a ball-contour tool appears at the mouse cursor when the mouse is positioned in the original image view. Contouring can be controlled in the following way:

- Create initial contour on a slice:
Click and hold left mouse button while moving the mouse.
- Edit existing contour:
Click and hold left mouse button while pushing against a contour from the inside or the outside.
- Add a contour besides existing contour in the same image plan
Switch to 'add/cut' mode via corresponding switch and click outside the existing contour to create an additional contour. Temporarily switching from 'edit' to 'add/cut' can be done by pressing and holding the Shift key
- Create a cut out in an existing contour
Switch to 'add/cut' mode via corresponding switch and click inside an existing contour to create a cut out contour. Temporarily switching from 'edit' to 'add/cut' can be done by pressing and holding the Shift key
- Delete contour on a single image plane via corresponding button
- Undo last contour action via corresponding button
- Change size of ball contour tool
Click right mouse button and drag mouse

Navigating through image slices during contouring can be done by the mouse scroll wheel, keyboard up/down arrow keys or scroll bar. To pan or zoom the image, switch from 'Contour' to 'Pan/Zoom' mode via corresponding switch. Temporarily switching from 'Contour' to 'Pan/Zoom' can be done by pressing and holding the Ctrl key.

Auto contouring

Auto-contouring functionality is available for several VOI types depending on the contour image modality and the contour purpose:

- Auto-contouring of the whole liver is available for CT image data for the purpose activity planning, dose simulation and dose review to create a VOI of type Whole Liver, for lung dose prediction to create a VOI of type Liver, and for SPECT based dose reconstruction to create a VOI of type Calibration. This is a semi-automatic process that requires user input. The user needs to define contours in the liver that will be used as input seeds for the auto-contour algorithm. This process is guided by a dedicated workflow for input definition and result review and acceptance.

- Auto-contouring of the lungs is available for CT image data for the purpose lung dose prediction and dose simulation to create a VOI of type Lung and for SPECT based dose reconstruction to create a VOI of type Calibration. This is a fully automatic process that does not require user input.
- Auto-contouring of tumors is available for PET image data and SPECT image data that is not used as the surrogate marker image for the purpose of dose simulation and dose review to create VOIs of type Tumor. The process is based on thresholding for which the user needs to define the threshold values.
- Auto-contouring of healthy tissue is available for SPECT image data for the purpose of dose simulation and dose review to create VOIs of type Healthy Tissue. The process is based on thresholding for which the user needs to define the threshold values.

Load contours from database

Contours that have previously been created by the user on an image series or that have been imported into the Q-Suite™ database by the user can be loaded for a VOI. Depending on the purpose, contours can only be used for the specific image series they originally belong to, for image series that share their frame of reference or for all image series. Loading contours is guided by a dedicated workflow for contour selection and result review and acceptance. If rules apply for the VOI (e.g. no overlap) the loaded contour will be adjusted to comply if needed.

6.9 FUSING IMAGE AND CONTOUR DATA

Fusion (aligning of images or contours with respect to each other in space) can be performed for several purposes. For all these purposes Q-Suite™ 2.1 has a dedicated screen in which two datasets can be displaced manually by the user with respect to each other.

Datasets for fusion

For the purpose of dose simulation the following datasets are used

- The VOI image set is used as background series
- The SPECT-CT of the surrogate marker is used as the overlay series

The user can choose whether the SPECT or the CT image of the surrogate marker is displayed. If the VOI image is a SPECT-CT or PET-CT image, the user can also choose whether the SPECT/PET or the CT image is displayed.

For the purpose of dose review the following datasets are used

- The VOI image set is used as background series
- As the overlay series the following data is used:
 - The SPECT-CT used as input for dose reconstruction in case of SPECT-CT based dose map
 - The SPECT used as input for dose reconstruction in case of SPECT based dose map
 - The post-treatment MGE MRI data used as input for dose reconstruction in case of MRI based dose map

The user can choose whether the SPECT or the CT image of overlay is displayed in case of SPECT-CT based dose map. If the VOI image is a SPECT-CT or PET-CT image, the user can also choose whether the SPECT/PET or the CT image is displayed.

For the purpose of fusing contours loaded from database the following datasets are used

- The contour image set is used as background series
- The contour data is used as the overlay series

Visualization

The data sets are displayed in three orthogonal orientations. The following options are available to change the appearance of these sets:

Background series:

Color Mapping: changing the color appearance of the background image

Overlay series:

Opacity: changing the opacity (transparency) of the overlay from 0 to 100%

Coverage: changing how far the overlay image covers the background image (horizontal and vertical)

Color Mapping: changing the color appearance of the overlay image

If the image data sets contain contours or for the purpose of fusing contours, the following options are available:

VOIs to show: choosing the VOIs of which the contours are displayed

Thickness: changing the thickness of the displayed contours

Displacement

The overlay image can be displaced with respect to the background either by translation in a plane or rotation in a plane. This translation and rotation can be done in all three orientations when the Mouse control mode 'Translate/Rotate' is active. For translation click the left mouse button and drag the image around, for rotation click the right mouse button, move the cursor outside the circle and rotate around the circle to rotate the image.

6.10 DATABASE MAINTENANCE

All imported image data and structure sets, created dose simulation sets, created dose reconstructions and created dose review sets are stored in the Q-Suite™ 2.1 database. Deleting data can be done in the following ways:

- In the database screen all data of a selected patient can be deleted at once
- In the database screen the following data of a selected patient can be selectively deleted
 - o unused image data
 - o unused dose maps
 - o structure sets (VOIs with contours)
- In 'Select dose reconstruction composition', incomplete (interrupted) compositions can be deleted
- In 'Create and fuse dose review sets' dose review sets can be deleted.
- In 'Create and fuse dose simulation sets' dose simulation sets can be deleted

6.11 SETTINGS

To enter the settings page an additional login code must be provided. The 'Settings' page contains four screens in tabulated form: 'General', 'DICOM', 'MRI calibration' and 'SPECT calibration'.

6.11.1 General

Interface Language

The user can change the Q-Suite™ 2.1 user interface language by selecting one of the available languages. Changing the language requires a restart of Q-Suite™ 2.1.

Database limit

A limit can be set to the maximum size of the Q-Suite™ 2.1 database. Setting a limit that is lower than the data size already stored in the database is prevented.

Liver segmentation algorithm

Q-Suite™ 2.1 is shipped with two different algorithms for auto contouring of the Liver based on CT. Version 1 was introduced in Q-Suite™ 2.0. Version 2 is an updated version that might improve auto-contouring results. The user can choose which version to use.

6.11.2 MRI calibration

For MRI the Holmium relaxivity can be adjusted. The user can choose between a predefined value, based on measurements described in scientific literature [van de Maat GH et al. Eur Radiol 2013;23:827-35] or a user defined value based on own relaxivity measurements. Recommendations to measure the holmium relaxivity can be found in the referenced paper. The holmium relaxivity value that is used directly influences the MRI based dose reconstruction. Using a higher relaxivity value than valid for the used holmium microspheres will result in reconstructed dose values that are too low. Using a lower relaxivity value than valid for the used holmium microspheres will result in reconstructed dose values that are too high.

! Note ! changing the calculation parameters will have impact on the accuracy of the dose reconstructions

6.11.3 SPECT calibration

The user can store multiple SPECT calibration factors in Q-Suite™ 2.1 which can be selected during preparation of a dose reconstruction. Factors can be added, edited and removed. Additional information that can be stored for each factor includes: used system, used collimator and used reconstruction method. Only the calibration factor is mandatory, the other fields are for information only.

6.11.4 DICOM

Local Import folder

A default folder can be set for importing image data in the Q-Suite™ 2.1 database. This folder is automatically used for import from a local folder on starting a new Q-Suite™ 2.1 session. This folder can be set by clicking 'Change' and selecting a specific folder of choice in the Explorer window that appears.

PACS communication settings

For PACS communication for data import and export, the local application entity and remote application entities can be configured.

For storing reports to a PACS the following must be defined:

- Local AE title: application entity title of the local Q-Suite™ system
- Remote AE title: application entity title of the PACS server to connect to
- Remote Host name or IP address: host name or IP address of the PACS server
- Remote Port Number: port number used by the PACS for handling storage requests
- Remote AE is archive: enabled

For retrieving data from a PACS the following must be defined

- Local AE title: application entity title of the local Q-Suite™ system
- Local IP address: IP address of the local Q-Suite™ system (C-MOVE only)
- Local Port number: port number to be used by Q-Suite™ for handling storage requests (C-MOVE only)
- Remote AE title: application entity title of the PACS server to connect to
- Remote Host name or IP address: host name or IP address of the PACS server
- Remote Port Number: port number used on the PACS server for query/retrieve requests
- Query/Retrieve protocol: protocol to be used for retrieving data (C-GET or C-MOVE)

For pushing data from a PACS to the Q-Suite™ inbox the following must be defined

- Local AE title: application entity title of the local Q-Suite™ system
- Local IP address: IP address of the local Q-Suite™ system
- Local Port number: port number to be used by Q-Suite™ system for handling storage requests
- Remote AE title: application entity title of the PACS server from which to receive data
- Remote AE pushes to inbox: enabled

After configuration of a remote application entity, the connection can be tested via the corresponding button

7. TROUBLESHOOTING AND MAINTENANCE

7.1 REPORTING INCIDENTS

If an serious incident related to Q-Suite™ occurs, report the incident to Quirem Medical B.V.: info.quirem@terumo-europe.com. Seious incidents should also be reported to the competent authorities.

7.2 TROUBLESHOOTING

For reporting errors in the software or help on using Q-Suite™ functionality, please contact your local Q-Suite™ representative.

7.3 MAINTENANCE

The Q-Suite™ 2.1 software does not require any preventative maintenance, such as log files to be checked or to be cleared, database maintenance, etc.

7.4 SUPPORT

The lifetime of Q-Suite™ 2.1 is 5 years. During this time, Quirem Medical B.V. will periodically provide patches to fix non-critical issues on the newest release. Critical updates (safety bugs) will be provided for all releases within their lifetime. Quirem Medical B.V. may decide to release a critical update to more or all releases in use in the field. In case of a safety incident that relates to the use of Q-Suite™ 2.1, Quirem Medical B.V. will provide adequate follow-up actions during the lifetime as required by law.

8. TECHNICAL DATA

8.1 CALCULATIONS

! Note ! *The accuracy of the dose values calculated by Q-Suite™ 2.1 always depends on the quality and the content of the input images and the way the input data is prepared. Recommendations for image acquisition and preparation of data are provided in section 6.2.*

8.1.1 Volume calculations

The total volume of the VOI is determined by the total number of image voxels included by the contours multiplied by the voxel volume where the voxel volume is determined based on the DICOM data in the following way:

$$\text{pixel spacing}_1 \times \text{pixel spacing}_2 \times \text{spacing between slices} \quad [1]$$

With pixel spacing_1 and pixel spacing_2 the first and second element of the DICOM tag *Pixel Spacing* (Tag ID: 0028, 0030) and spacing between slices the value of DICOM Tag *Spacing Between Slices* (Tag ID: 0018,0088).

8.1.2 Voxel inclusion VOIs

Since contour lines created by the user are continuous and smooth whereas the underlying image data is a discrete voxel grid, visually a part of a certain voxel may be outside the contour whereas the other part is inside the contour. Whether a voxel is included in a VOI or not is determined based on weighting. If the center point of a voxel is inside the contour it is included in the VOI, otherwise it is excluded. Since the ratio between the perimeter and the surface of the contour increases for decreasing contour size, the smaller the drawn contour, the bigger the discrepancy between the displayed contour line and actual included volume.

8.1.3 Treatment activity calculations

Treatment activity calculation per user defined Target liver VOI is performed based on the following formula:

$$A_i [\text{MBq}] = D_i [\text{Gy}] \times 63 [\text{MBq/J}] \times M_i [\text{kg}] \quad [2]$$

Where A_i is the activity for Target Liver VOI i , M_i is the mass of Target Liver VOI i , D_i is the user defined aimed dose for Target Liver VOI i and 63 a holmium specific factor. The mass of a Target Liver VOI is calculated based on the total volume included in the user defined contours for that VOI multiplied by a density of 1.06 kg/L for liver tissue. See 8.1.1 for the volume calculation method.

The resulting Whole Liver dose is calculated based on the following formula:

$$D_{\text{whole liver}} [\text{Gy}] = \frac{\sum A_i [\text{MBq}]}{M_{\text{whole liver}} [\text{kg}] \times 63 [\text{MBq/J}]} \quad [3]$$

Where $\sum A_i$ is the total activity of all Target Liver VOIs i , $M_{\text{whole liver}}$ is the mass of the Whole Liver. The mass of a Target Liver VOI is calculated based on the total volume included in the user defined contours for that VOI multiplied by a density of 1.06 kg/L for liver tissue. See 8.1.1 for the volume calculation method.

8.1.4 Image based lung dose prediction

Lung dose prediction is performed in two steps. In the first step, the fraction of SPECT counts in the lung VOIs is calculated based on the following formulas:

$$L_i = \frac{\text{counts of lung } i}{\text{counts of total lung} + \text{counts of liver}} \times 100 \quad [4]$$

$$L = \frac{\text{counts of total lung}}{\text{counts of total lung} + \text{counts of liver}} \times 100 \quad [5]$$

With L_i the SPECT count fraction per lung VOI i and L the SPECT count fraction for the total Lung VOI, $\text{counts of lung } i$ the sum of voxel values (counts) of all voxels included in the user defined lung VOI i , $\text{counts of total lung}$ the sum of voxel values (counts) of all voxels included in all user defined Lung VOIs and counts of liver the sum of all voxels values (counts) included in the user defined VOI liver.

In the second step, the mean lung dose for the user defined lung VOIs and the mean dose for the total lung VOI is calculated based on the following formula:

$$D_{lung\ i} \text{ [Gy]} = \frac{A[\text{MBq}] \times L_i/100}{M_{lung\ i}[\text{kg}] \times 63 [\text{MBq / l}]} \quad [6]$$

$$D_{lung} \text{ [Gy]} = \frac{A[\text{MBq}] \times L/100}{M_{lung} [\text{kg}] \times 63 [\text{MBq / l}]} \quad [7]$$

With $D_{lung\ i}$ the mean dose per of lung VOI i , D_{lung} the mean dose of the total lung VOI, A the planned treatment activity, L_i the SPECT count fraction in lung VOI i , L the SPECT count fraction of the total lung VOI, $M_{lung\ i}$ the mass of the lung VOI i , M_{lung} the mass of the total lung VOI and 63 a holmium specific factor. For planar scintigraphy data, only equation [x] is used.

$M_{lung\ i}$ and M_{lung} are calculated based on the total volume included in the user defined contours for that VOI multiplied by a density of 0.3 kg/L for lung tissue in case of SPECT-CT data.

In case of planar scintigraphy data, M_{lung} is calculated by the volume defined by the user multiplied by a density of 0.3 kg/L for lung tissue

When excluding the caudal part of the lungs from lung dose analyses, lung VOI voxels that overlap with the liver VOI + defined margin are excluded from the lung VOI both for calculation of *counts of total lung/counts of lung_i* and calculation of $M_{lung}/M_{lung\ i}$.

8.1.5 Dose simulation

Dose simulation is performed as follows. In the first step, per Target Liver VOI, SPECT counts are converted into activity using the following formulas:

$$CF_{Target\ Liver\ VOI} \text{ [MBq/count]} = \frac{A_{Target\ Liver\ VOI} \text{ [GBq]} \times 1000}{Counts_{Target\ Liver\ VOI}} \quad [8]$$

$$A_{Target\ Liver\ voxel} = Counts_{voxel} \times CF_{Target\ liver\ VOI} \quad [9]$$

With $Counts_{Target\ Liver}$ the sum of voxel values (counts) of all voxels included in the user defined Target Liver VOI, $A_{Target\ Liver}$ is the treatment activity defined by the user for the dose simulation for that Target Liver VOI, $Counts_{voxel}$ is the voxel value (counts) of a certain voxel that belongs to the Target Liver VOI and A_{voxel} is the activity calculated for that voxel.

When accounting for lung shunt is activated, the calibration factor per Target Liver VOI is corrected for the lung shunt fraction according to:

$$LSF = \frac{Counts_{Total\ Lung}}{Counts_{Total\ Lung} + Counts_{Total\ Target\ Liver}} \quad [10]$$

$$CF_{Target\ Liver\ VOI} \text{ [MBq/count]} = \frac{A_{Target\ Liver\ VOI} \text{ [GBq]} \times (1 - LSF) \times 1000}{Counts_{Target\ Liver\ VOI}} \quad [11]$$

Where LSF is the total lung shunt fraction, $Counts_{Total\ Lung}$ is the sum of voxel values (counts) of all voxels included in all user defined Lung VOIs, $counts_{Total\ Target\ Liver}$ is the sum of voxel values (counts) of all voxels included in all user defined Target Liver VOIs.

Per lung VOI, SPECT counts are converted into activity using the following formulas:

$$CF_{Lung} \text{ [MBq/count]} = \frac{A_{Target\ Liver} \text{ [GBq]} \times 1000}{Counts_{Total\ Lung} + Counts_{Total\ Target\ Liver}} \quad [12]$$

$$A_{Lung\ voxel} = Counts_{voxel} \times CF_{Lung} \quad [13]$$

As a second step, the dose per voxel D_{voxel} for voxels included in user defined VOIs is calculated based on the following formula:

$$D_{voxel} \text{ [Gy]} = \frac{A_{voxel} \text{ [MBq]}}{M_{voxel} [\text{kg}] \times 63 [\text{MBq / l}]} \quad [14]$$

With M_{voxel} the mass of the specific voxel. M_{voxel} is calculated based on the voxel volume (see 8.1.1) and the density of the tissue of the VOI the voxel belongs to.

The following densities are used:

- VOI type Whole Liver: 1.06 kg/L
- VOI type Target Liver: 1.06 kg/L
- VOI type Healthy Liver: 1.06 kg/L

- VOI type Lung: 0.3 kg/L
- VOI type Tumor: 1.06 kg/L

8.1.6 SPECT based dose reconstruction

Dose reconstruction based on SPECT as input is performed as follows. First, SPECT counts are converted into activity per voxel A_{voxel} using the following formula:

$$A_{voxel} = counts_{voxel} \times CF \quad [15]$$

Here $counts_{voxel}$ is the voxel value (counts) of a certain voxel and CF the calibration factor which either a value selected by the user (in case of pre-measured method) or calculated as follows (in case of patient specific calibration method)

$$CF [MBq/count] = \frac{A_{assumed} [MBq]}{Counts_{Calibration\ VOI}} \quad [16]$$

Here $Counts_{calibration\ VOI}$ is the sum of voxel values (counts) of all voxels included in all user defined calibration VOIs and $A_{assumed}$ is the user defined radiation activity that is assumed to be within that VOI/those VOIs at time of injection.

Second, the dose per voxel D_{voxel} is calculated based on convolution of the activity with a holmium Dose Point Kernel for liver tissue taking into account the voxel mass or based on multiplication of the activity with a holmium specific energy factor taking into account the voxel mass. For both methods a density of 1.06 kg/L is used.

8.1.7 MRI based dose reconstruction

Dose calculations using input MR images include a sequence of processing steps to convert image intensities into absorbed dose. The effect that the presence of the paramagnetic holmium-166 microspheres induces changes in the MR parameter R_2^* is utilized to determine the local concentration of microspheres. Measuring this change requires pre-treatment and post-treatment MGE datasets as described in section 6.2.

First the MRI images are converted into radiation activity per voxel via the following method:

- 1) R_2^* values are determined by applying a mono-exponential fit to the signal intensities of both the pre- and a post-treatment MGE data sets using

$$S(t) = S(0) \cdot e^{-R_2^* t} \quad [17]$$

With t the echo time of a gradient echo, $S(t)$ the gradient echo voxel signal at echo time t and $S(0)$ the signal at $t = 0$. Signal values $S(t)$ below the user chosen SNR threshold are excluded from the fit where the noise is determined from the standard deviation of the voxel values included in the user defined 'Noise measure' VOI

- 2) For voxels for which R_2^* cannot not be determined because the number of valid gradient echo signals after SNR threshold exclusion is lower than the number of required echoes, S_0 -fitting is applied. First $S(0)$ for such a voxel is estimated by the mean $S(0)$ value of neighboring voxels using region growing. Starting with a selection region of 3 by 3 voxels surrounding the undetermined voxel, voxels are selected for the fit was successful. The selection region is extended until at least the *minimum neighbours S_0 -fit* are included or *range neighbours S_0 -fit* is reached. Second, $S(0)$ is included in the fit to determine R_2^* for the voxel. In case no $S(t)$ meets the SNR-threshold, R_2^* is estimated based on noise level and first echotime according to

$$R_{2\ lim}^* = \frac{\ln(S(0)/SNR_{threshold})}{TE_1} \quad [18]$$

- 3) ΔR_2^* values are determined by

$$\Delta R_{2\ voxel}^* = R_{2\ voxel\ post}^* - R_{2\ mean\ pre}^* \quad [19]$$

With $R_{2\ mean\ pre}^*$ the mean R_2^* value of voxels included in the user defined liver VOI on the pre-treatment data and $R_{2\ voxel\ post}^*$ the R_2^* value of voxels of the post-treatment data included in the user defined 'Liver' VOI.

- 4) ΔR_2^* values are converted into amount of holmium $amount_{Ho}$ by

$$amount_{Ho} = \frac{\Delta R_2^*}{r_2^*} \times voxel\ volume \quad [20]$$

With r_2^* the chosen holmium relaxivity factor corrected for the used Magnetic Field Strength for image acquisition and the user defined holmium content of the microspheres

- 5) The holmium amount is converted into units of radiation activity by applying the user defined specific activity of the microspheres.

Second, the dose per voxel D_{voxel} for voxels included in the user defined liver VOIs is calculated based on convolution of the activity with a holmium Dose Point Kernel for liver tissue taking into account the voxel mass or based on multiplication of the activity with a holmium specific energy factor taking into account the voxel mass. For both methods a density of 1.06 kg/L is used.

8.2 ACCURACY

8.2.1 SPECT based dose reconstruction

Since dose calculations using input SPECT images are performed based on a direct translation of voxel intensities (counts) into absorbed dose, the SPECT voxel intensities should represent the actual activity concentration in that specific volume. All factors that lead to inaccurate reconstruction of the activity concentration in the final SPECT image will therefore introduce deviations in the dose calculations. Factors known for their influence of image quality include patient movement during acquisition, detector and collimator response and used reconstruction techniques.

Dose deviations introduced by Q-Suite™ 2.1 in the calculation process of transforming voxel intensity information into absorbed dose have been measured based on computer generated digital phantom data representing body, liver and tumor volume with specified intensities for each of these volumes correlating with a certain activity concentration. Measurements have been performed for a range of activity concentrations in the tumor volume.

A deviation of less than 1% has been determined for a dose range of 0 to 1000 Gy for this digital phantom data. Any of the above mentioned factors that influence the accuracy can lead to a deviation that is larger than 1%.

8.2.2 MRI Based dose reconstruction

Factors that directly influence the accuracy of MRI based dose calculations include the presence of magnetic field distortions not caused by the microspheres, the signal-to-noise level in the images, the presence of image artifacts induced by for example patient/organ movement and used reconstruction techniques.

The implemented R_2^* based calculation model is sensitive to three other effects:

- deviations in voxel R_2^* values from the mean value of the entire liver volume (including all other tissue) in the pre-treatment image set will induce dose deviations.
- changes in voxel R_2^* values between the pre and post-treatment case induced by for example changes in tissue characteristics will induce dose deviations.
- very high concentrations of microspheres will lead to signal decay too fast to be able to measure, leading to underestimation of the dose.

The accuracy of Q-Suite™ 2.1 for MR based dose calculations has been validated on computer generated datasets simulating real image data in respect to expected image intensities but free of noise and artifacts. For these datasets the following observations have been made:

- in case the voxels of interest do have pre-treatment R_2^* values equal to the mean value of the entire volume (including liver and tumor) and there are no tissue induced changes, errors are limited to less than 1% for a dose range of 0 to 1000 Gy when appropriate echo times are used.
- In case the voxels of interest do have pre-treatment R_2^* values that are lower or higher than the mean value of the entire volume (including liver and tumor), the dose is respectively under- or overestimated where the actual deviation depends on the specific activity of the microspheres. Dose errors of 9 Gy were observed for R_2^* deviations of 10 s^{-1} in case of a specific activity of 6.7 MBq/mg but can increase in case of larger R_2^* deviations and higher specific activity

Any of the factors that were described before that influence the MR based dose calculations can lead to local deviations larger than the values mentioned in this section.

8.3 SYSTEM CONFIGURATIONS

Q-Suite™ 2.1 configurations

There is only a single configuration of Q-Suite™ 2.1. Customization of Q-Suite™ 2.1 is not supported.

Required and recommended computer systems

Q-Suite™ 2.1 requires a x64 based computer running Microsoft Windows 7, 8.1 or 10. In addition, the Microsoft .NET framework version 4.7.2 framework (or higher) should be installed.

The following hardware is recommended to ensure a good user experience:

- Dual-core CPU, 2.0GHz or higher
- 16 GB RAM
- SSD with 20 GB free disk space
- Screen resolution 1920 x 1080 or higher

Data size and performance

The responsiveness of Q-Suite™ 2.1 and the time needed to process data depends on the size of the image data to be processed and the computer system configuration like processor speed and available RAM. Q-Suite™ 2.1 accepts large image datasets, but the user shall be aware that performing dose reconstructions and evaluation on high resolution/large volume datasets may cause slow or delayed system response.

8.4 UNITS

The following quantities and units are used in Q-Suite™ 2.1:

Quantity	Unit
Material content	% (weight/weight)
Specific activity	MBq/mg
Dose	Gy
Volume	mL
Activity	GBq, MBq
Magnetic field strength	T
Energy	keV
Time duration	H, s, ms,
Relaxivity	$s^{-1} \cdot ml \cdot mg^{-1} \cdot T^{-1}$
Distance	mm
Weight	kg
Height	cm

9. SYMBOLS AND NAMES

The following symbols are used in Q-Suite™ 2.1

	Manufacturer		Contour was drawn manually
	Batch Number		Settings
	Catalogue number	 eIFU on DVD www.quirem.com/IFU	Consult instructions for use (electronic format on DVD or refer to: www.quirem.com/IFU)
	Copyright		CE-mark that show European Conformity
	Brightness/contrast		Add image snapshot to report
	Patient		Delete
	Complete		Close Remove from list
	Attention		Expand Down
	Search		Collapse Up
	Increase size		Re-center
	Add		Decrease size
	Export		Link
	DICOM inbox is empty, not empty		Clean up
	Dose plan is locked		Show note
	Edit		Local DICOM Storage SCP is running Test connectivity with remote DICOM AE
	Horizontal		Vertical
	Unlock dose plan		Undo changes
	Previous		Next
	Toggle visibility		Quantity

Q-Suite™ is a trademark of Quirem Medical

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