|    |          |       |
|----|----------|-------|
| en | Needle   | 2-3   |
| fr | Aiguille | 4-5   |
| de | Kanüle   | 6-7   |
| es | Aguja    | 8-9   |
| pt | Agulha   | 10-11 |
| it | Ago      | 12-13 |
| nl | Naald    | 14-15 |
| sv | Kanyl    | 16-17 |
| da | Kanyle   | 18-19 |
| no | Kanyle   | 20-21 |
| fi | Neula    | 22-23 |
| el | Βελόνα   | 24-25 |
| ru | Игла     | 26-27 |
| pl | Igła     | 28-29 |
| hu | Tű       | 30-31 |
| cs | Jehla    | 32-33 |
| sk | Ilha     | 34-35 |
| tr | İğne     | 36-37 |
| et | Nõel     | 38-39 |
| lv | Adata    | 40-41 |
| lt | Adata    | 42-43 |
| sl | Igla     | 44-45 |
| sr | Igla     | 46-47 |
| ro | Ac       | 48-49 |
| bg | Игла     | 50-51 |
| uk | Голка    | 52-53 |
| hr | Igla     | 54-55 |

## Instructions for Use

**ENGLISH**

### SYMBOL EXPLANATION

|                |                                                                                            |                                                                   |                                 |
|----------------|--------------------------------------------------------------------------------------------|-------------------------------------------------------------------|---------------------------------|
| <b>LOT</b>     | Batch code                                                                                 | Do not re-use                                                     | Non pyrogenic                   |
| <b>REF</b>     | Catalogue number                                                                           | Do not use if package is damaged and consult instructions for use | Single Sterile Barrier System   |
| <b>Rx only</b> | Caution: Federal law (USA) restricts this device to sale by or on the order of a physician | Fragile, handle with care                                         | Sterilized using ethylene oxide |
|                | Consult Instructions for Use                                                               | Keep away from sunlight                                           | Temperature limit               |
|                | Contents                                                                                   | Keep dry                                                          | Use-by date                     |
|                | Date of Manufacture                                                                        | Manufacturer                                                      |                                 |
|                | Do not re-sterilize                                                                        | Medical Device                                                    |                                 |

### PRODUCT DESCRIPTION

Non-Toxic

| Needle Gauge | Needle Length | Needle Outer Diameter | Needle Length | Wall Thickness | Bevel Angle |
|--------------|---------------|-----------------------|---------------|----------------|-------------|
| 27G          | 1/2"          | 0.40 mm               | 12 mm         | Thin Wall (TW) | 12°         |
| 29G          | 1/2"          | 0.33 mm               | 12 mm         | Thin Wall (TW) | 12°         |

### INTENDED PURPOSE

The K-Pack Enhance Needle is a sterile, hypodermic needle, for single use, intended to inject fluids into, or withdraw fluids from parts of the body below the surface of the skin.

### INDICATIONS

The K-Pack Enhance Needle is for general application - for injection of fluids or withdrawal of fluids.

### CONTRA-INDICATIONS

No contra-indications.

### PATIENT TARGET GROUP

Intended for general application.

### INTENDED USERS

Healthcare professional or lay person.

### CLINICAL BENEFIT

The K-Pack Enhance Needle has an indirect clinical benefit (indirect medical effect) since it is used for injection or withdrawal of fluids.

### WARNINGS

- Do not use if unit package (cap and case) is damaged.
- Use immediately after opening the unit package (cap and case).
- Do not re-cap.
- For single use only. Do not re-use. Do not re-sterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- The needle is made of stainless steel containing nickel and cobalt. Cobalt is classified as CMR\* 1B and is present in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

\*CMR = Carcinogenic, mutagenic or toxic to reproduction (CLP Regulation EU 1272/2008)

### PRECAUTIONS


- Do not use for ophthalmic injections.
- If the needle is bent or damaged, no attempt should be made to straighten the needle or use the product.
- After use, dispose of safely as medical waste in a sharps disposal container and/or according to health institution policies. The product is biohazardous and is physically hazardous due to its sharp edge.

### DIRECTIONS FOR USE


Aseptic technique, proper preparation and continued protection of the injection site are essential. Observe Universal Precautions on **ALL** patients.

**Precaution** Handle with care to avoid needle sticks. Keep hands behind the needle at all times during use and disposal.

#### COMPONENTS OVERVIEW


### 1. PREPARING FOR USE – OPENING

**Caution:** Avoid touching the hub after removing the cap.


#### COMPONENTS OVERVIEW

### 2. PREPARING FOR USE – ATTACHING


Firmly attach the hub onto the (pre-)filled syringe.

**Caution:** Do not overtighten.

**Luer slip:** pushing and twisting the needle

**Luer lock:** twisting the needle


Pull the case straight off.


### 3. INJECTION AND/OR WITHDRAWAL OF FLUIDS

Perform injection and/or withdrawal of fluids according to established technique.

### 4. DISPOSAL

**Caution:** Do not cover the needle with the case after use.

After use, dispose needle and syringe safely as medical waste in a sharps disposal container and/or according to health institution policies.

### PRECAUTIONS FOR STORAGE

Store between 2°C and 30°C. Keep dry. Keep away from sunlight. Fragile, handle with care.

### REPORT OF INCIDENT

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

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