en	SYRINGE WITH FIXED NEEDLE	2-3
fr	SERINGUE AVEC AIGUILLE SERTIE	4-5
de	SPRITZE MIT EINGESCHWEISSTER KANÜLE	6-7
es	JERINGA CON AGUJA PREMONTADA	8-9
pt	SERINGA COM AGULHA PREMONTADA	10-11
it	SIRINGA CON AGO FISSO	12-13
nl	SPUIT MET VASTE INJECTIENAALD	14-15
sv	INJEKTIONSSPRUTA MED INTEGRERAD KANYL	16-17
da	SPRØJTE MED INTEGRERET KANYLE	18-19
no	SPRØYTE MED INTEGRERT KANYLE	20-21
fi	RUISKU KIINTEÄLLÄ KANYYLILLÄ	22-23
el	ΣΥΡΙΓΓΑ ΜΕ ΜΗ ΑΠΟΣΠΩΜΕΝΗ ΒΕΛΟΝΑ	24-25
ru	ШПРИЦ СО ВСТРОЕННОЙ ИГЛОЙ	26-27
pl	STRZYKAWKA POŁĄCZONA Z IGŁĄ	28-29
hu	TŰVEL SZERELT FECSKENDŐ	30-31
cs	STŘÍKAČKA S PERMANENTNĚ FIXOVANOU JEHLOU	32-33
sk	STRIEKAČKA SO ZABUDOVANOU IHLOU	34-35
tr	SABIT İĞNELİ ŞİRINGA	36-37
et	SÜSTAL FIKSEERITUD NÕELAGA	38-39
lv	ŠĻIRCE AR FIKSĒTU ADATU	40-41
lt	ŠVIRKŠTAS SU NENUIMAMA ADATA	42-43
sl	BRIZGA S FIKSNO IGLO	44-45
sr	ŠPRIC SA FIKSIRANOM IGLOM	46-47
ro	SERINGA CU AC FIX	48-49
bg	СПРИНЦОВКА С ВГРАДЕНА ИГЛА	50-51
uk	ШПРИЦ ІЗ ФІКСОВАНОЮ ГОЛКОЮ	52-53
hr	ŠPRICA S FIKSNOM IGLOM	54-55

Instructions for Use

ENGLISH

SYMBOL EXPLANATION

	Batch code
	Catalogue number
	Consult Instructions for Use
	Contents
	Date of Manufacture
	Do not re-sterilize

	Do not re-use
	Do not use if package is damaged and consult instructions for use
	Fragile, handle with care
	Keep away from sunlight
	Keep dry
	Manufacturer

	Medical Device
	Non pyrogenic
	Single Sterile Barrier System
	Sterilized using ethylene oxide
	Temperature limit
	Use-by date

Rx only **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

PRODUCT DESCRIPTION

Non-Toxic

Needle Gauge	Needle Length	Needle Outer Diameter	Needle Length	Wall Thickness	Syringe Nominal Capacity/ Nominal Volume (v)	Dead Space Volume
27G	1/2"	0.4 mm	12 mm	Regular Wall (RW)	1 ml	0.01 ml
Tolerance on Graduated Capacity						
Syringe Nominal Capacity/ Nominal Volume (v)		When syringe filled to a volume less than half of the nominal capacity			When syringe filled to a volume equal or greater than half of the nominal capacity	
1 ml		± (1.50% of v + 2% of expelled volume)			2% of expelled volume	

INTENDED PURPOSE

TERUMO SYRINGE with Fixed Needle for 1 mg/0.2 ml leuprolide acetate is a hypodermic syringe intended for manual aspiration of leuprolide acetate and for the injection of this solution into parts of the body below the surface of the skin.

INDICATIONS

The device is indicated for subcutaneous injection of leuprolide acetate in the palliative treatment of advanced prostatic cancer.

CONTRAINDICATIONS

No contraindications.

PATIENT TARGET GROUP

Men with advanced prostatic cancer.

INTENDED USERS

Healthcare professional or lay person.

CLINICAL BENEFIT

The device has an indirect clinical benefit (indirect medical effect) since it is used for subcutaneous injection of leuprolide acetate in the palliative treatment of advanced prostatic cancer.

WARNINGS

- Do not use if unit package is damaged.
- Use immediately after opening the unit package.
- For single use only. Do not reuse. Do not re-sterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- The needle is made of stainless steel containing nickel and cobalt. Cobalt is classified as CMR* 1B and is present in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.

*CMR = Carcinogenic, mutagenic or toxic to reproduction (CLP Regulation EU 1272/2008)










PRECAUTIONS

- The syringe is intended for use immediately after filling, as it is not suitable for containing fluids over extended periods of time.
- If the needle is bent or damaged, no attempt should be made to straighten the needle or use the product.
- After use, dispose of safely as medical waste in a sharps disposal container and/or according to health institution policies. The product is biohazardous and is physically hazardous due to its sharp edge.

Instructions for Use
ENGLISH
DIRECTIONS FOR USE

Aseptic technique, proper skin preparation and continued protection of the site are essential. Observe Universal Precautions on ALL patients. **PRECAUTION** Handle with care to avoid needle sticks. Keep hands behind the needle at all times during use and disposal.

1. Open the blister package.
2. Pull the syringe plunger back to the graduation line for the prescribed dose. This fills the syringe with air.
3. Remove the cap. Do not touch the needle.
4. Place the bottle on a clean, flat surface and push the needle through the center of the rubber stopper on the bottle. Push the plunger all the way in to inject air into the bottle.
5. Keep the needle in the bottle. Lift the bottle and turn it straight upside down. Check to see that the needle tip is in the solution.
6. With the needle tip in the solution, slowly pull back the plunger until the syringe fills to the graduation line for the prescribed dose. If any air bubbles appear in the syringe, remove them by pushing the plunger up slowly. With the needle tip still in the solution, pull the plunger back until it is once more at the graduation line for the prescribed dose.
7. Perform injection procedure according to established technique.
8. Do not re-cap.
9. Dispose of safely.

PRECAUTIONS FOR STORAGE

Store between 2°C and 30°C. Keep dry. Keep away from sunlight. Fragile, handle with care.

REPORT OF INCIDENT

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

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