en	SYRINGE WITH FIXED NEEDLE	2-3
fr	SERINGUE AVEC AIGUILLE SERTIE	4-5
de	SPRITZE MIT EINGESCHWEISSTER KANÜLE	6-7
es	JERINGA CON AGUJA PREMONTADA	8-9
pt	SERINGA COM AGULHA PREMONTADA	10-11
it	SIRINGA CON AGO FISSO	12-13
nl	SPIJT MET VASTE INJECTIENAALD	14-15
sv	INJEKTIONSSPRUTA MED INTEGRERAD KANYL	16-17
da	SPRØJTE MED INTEGRERET KANYLE	18-19
no	SPRØYTE MED INTEGRERT KANYLE	20-21
fi	RUISKU KIINTEÄLLÄ KANYYLILLÄ	22-23
el	ΣΥΡΙΓΓΑ ΜΕ ΜΗ ΑΠΟΣΠΩΜΕΝΗ ΒΕΛΟΝΑ	24-25
ru	ШПРИЦ СО ВСТРОЕННОЙ ИГЛОЙ	26-27
pl	STRZYKAWKA POŁĄCZONA Z IGŁĄ	28-29
hu	TŰVEL SZERELT FECSKENDŐ	30-31
cs	STŘÍKAČKA S PERMANENTNĚ FIXOVANOU JEHLOU	32-33
sk	STRIEKAČKA SO ZABUDOVANOU IHLOU	34-35
tr	SABIT İĞNELİ ŞIRINGA	36-37
et	SÚSTAL FIKSEERITUD NÕELAGA	38-39
lv	ŠĻIRCE AR FIKSĒTU ADATU	40-41
lt	ŠVIRKŠTAS SU NENUIMAMA ADATA	42-43
sl	BRIZGA S FIKSNO IGLO	44-45
sr	ŠPRIC SA FIKSIRANOM IGLOM	46-47
ro	SERINGA CU AC FIX	48-49
bg	СПРИНЦОВКА С ВГРАДЕНА ИГЛА	50-51
uk	ШПРИЦ ІЗ ФІКСОВАНОЮ ГОЛКОЮ	52-53
hr	ŠPRICA S FIKSNOM IGLOM	54-55

**Instructions for Use**
**ENGLISH**
**SYMBOL EXPLANATION**

<b>LOT</b>	Batch code	Do not re-use	<b>MD</b>	Medical Device
<b>REF</b>	Catalogue number	Do not use if package is damaged and consult instructions for use		Non pyrogenic
<b>i</b>	Consult Instructions for Use	Fragile, handle with care		Single Sterile Barrier System
<b>#</b>	Contents	Keep away from sunlight		Sterilized using ethylene oxide
	Date of Manufacture	Keep dry	$2^{\circ}\text{C}$ $-30^{\circ}\text{C}$	Temperature limit
	Do not re-sterilize	Manufacturer		Use-by date

**PRODUCT DESCRIPTION**

Non-Toxic

Needle Gauge	Needle Length	Needle Outer Diameter	Needle Length	Wall Thickness	Syringe Nominal Capacity/ Nominal Volume (v)	Dead Space Volume
27G	1/2"	0.4 mm	12 mm	Regular Wall (RW)	1 ml	0.01 ml
Tolerance on Graduated Capacity						
Syringe Nominal Capacity/ Nominal Volume (v)		When syringe filled to a volume less than half of the nominal capacity		When syringe filled to a volume equal or greater than half of the nominal capacity		
1 ml		$\pm (1.50\% \text{ of } v + 2\% \text{ of expelled volume})$		2% of expelled volume		

**INTENDED PURPOSE**

TERUMO SYRINGE with Fixed Needle for 600 IU/ml FSH (Follicle Stimulating Hormone) is a hypodermic syringe intended for manual aspiration of GONAL-f® Multi-Dose 600 IU/ml (follitropin alfa for injection) and for the injection of this solution into parts of the body below the surface of the skin.

**INDICATIONS**

The device is indicated for subcutaneous injection of fertility hormones to treat women suffering from hormonal deficiencies.

**CONTRAINDICATIONS**

No contraindications.

**PATIENT TARGET GROUP**

Women suffering from hormonal deficiencies.

**INTENDED USERS**

Healthcare professional or lay person.

**CLINICAL BENEFIT**

The device has an indirect clinical benefit (indirect medical effect) since it is used for subcutaneous injection of fertility hormones to treat women suffering from hormonal deficiencies.

**WARNINGS**


- Do not use if unit package is damaged.
  - Use immediately after opening the unit package.
  - For single use only. Do not reuse. Do not re-sterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
  - The needle is made of stainless steel containing nickel and cobalt. Cobalt is classified as CMR\* 1B and is present in a concentration above 0.1% weight by weight. Current scientific evidence supports that medical devices manufactured from stainless steel alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects.
- \*CMR = Carcinogenic, mutagenic or toxic to reproduction (CLP Regulation EU 1272/2008)

**PRECAUTIONS**

- The syringe is intended for use immediately after filling, as it is not suitable for containing fluids over extended periods of time.
- If the needle is bent or damaged, no attempt should be made to straighten the needle or use the product.
- After use, dispose of safely as medical waste in a sharps disposal container and/or according to health institution policies. The product is biohazardous and is physically hazardous due to its sharp edge.

**Instructions for Use****ENGLISH****DIRECTIONS FOR USE**

1. Peel open blister
2. Remove cap
3. Push plunger down
4. Fill syringe
5. Remove air bubbles
6. Inject
7. Do not re-cap
8. Dispose of safely

**PRECAUTIONS FOR STORAGE**

Store between 2°C and 30°C. Keep dry. Keep away from sunlight.  
Fragile, handle with care.

**REPORT OF INCIDENT**

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and to your national authority.

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