Ultimaster Tansei - Patient Information

| | |
|---|-----|
| Patient Information – ENGLISH (English) | 2 |
| Informations pour le patient – FRANÇAIS (French) | 6 |
| Informationen für den Patienten – DEUTSCH (German) | 10 |
| Información para el paciente – ESPAÑOL (Spanish) | 14 |
| Informações para o paciente – PORTUGUÊS (Portuguese)..... | 18 |
| Informazioni per il Paziente – ITALIANO (Italian) | 22 |
| Informatie voor de patiënt – NEDERLANDS (Dutch)..... | 26 |
| Patientinformation – SVENSKA (Swedish)..... | 30 |
| Oplysninger til patient – DANSK (Danish) | 34 |
| Informasjon til pasienten – NORSK (Norwegian)..... | 38 |
| Tietoja potilaalle – SUOMI (Finnish) | 42 |
| Πληροφορίες για τον ασθενή – ΕΛΛΗΝΙΚΆ (Greek) | 46 |
| Информация для пациента – РУССКИЙ (Russian) | 50 |
| Informacja dla pacjenta – POLSKI (Polish)..... | 54 |
| Információ betegeknek – MAGYAR (Hungarian) | 58 |
| Informace pro pacienta – ČESKÝ (Czech) | 62 |
| Informácie pre pacienta – SLOVENSKÝ (Slovak) | 66 |
| Hasta için bilgiler – TÜRKÇE (Turkish) | 70 |
| Teave patsiendile - EESTI (Estonian)..... | 74 |
| Informācija pacientam – LATVISKI (Latvian)..... | 78 |
| Informacija apie pacientą – LIETUVIŠKAI (Lithuanian)..... | 82 |
| Podatki za bolnika – SLOVENŠČINA (Slovenian) | 86 |
| Informacije za pacijenta – SRPSKI (Serbian)..... | 90 |
| Informații pentru pacient – ROMÂNĂ (Romanian)..... | 94 |
| Информация за пациента – БЪЛГАРСКИ (Bulgarian)..... | 98 |
| Інформація для пацієнта – УКРАЇНСЬКА (Ukrainian) | 102 |
| Informacije za pacijenta – HRVATSKI (Croatian)..... | 106 |

Patient Information – ENGLISH (English)

Ultimaster Tansei Sirolimus eluting coronary stent system

You recently received or will receive a sirolimus eluting coronary stent implanted in a coronary artery of your heart. After the procedure your physician will give you an implant card. **Please carry this patient implant card at all times and show it to any medical personnel who may be treating you.**

The following information about your stent is important for you to know:

1. When is Ultimaster Tansei stent used?

Ultimaster Tansei stent is used to improve blood supply to the heart muscle, alleviating your symptoms, such as chest pain and discomfort or shortness of breath.

2. Product description

The Ultimaster Tansei stent is a mesh-like tube made of a metal alloy, coated with a drug and a polymer. The stent is placed and expanded against the wall of the coronary artery to increase blood supply to the heart muscle. The polymer allows the drug to be released in a controlled way into the wall of the artery to reduce the chance of a reblockage.

Table 1: Information on the materials and substances used in the Ultimaster Tansei stent:

| Stent component | Material | Substance | Quantity |
|-----------------|-----------------------|------------|--|
| Stent platform | Cobalt Chromium alloy | Cobalt* | 46.38-56.95 (w/w%) |
| | | Chromium | 19.00-21.00 (w/w%) |
| | | Nickel | 9.00-11.00 (w/w%) |
| | | Carbon | 0.05-0.15(w/w%) |
| | | Manganese | 1.00-2.00 (w/w%) |
| | | Silicon | 0-0.40 (w/w%) |
| | | Phosphorus | 0-0.04 (w/w%) |
| | | Sulphur | 0-0.03 (w/w%) |
| | | Tungsten | 14.00-16.00 (w/w%) |
| | | Iron | 0-3.00 (w/w%) |
| Drug | Sirolimus | | 36-148 µg (stents with diameter 2.25; 2.5; 2.75 or 3.0 mm, depending on stent length) 35-151 µg (stents with diameter 3.5 or 4.0 mm, depending on stent length) |

| | | |
|---|---|--|
| Polymer | Poly(D,L-lactide-co-caprolactone) copolymer | 49-202 µg (stents with diameter 2.25; 2.5; 2.75 or 3.0 mm, depending on stent length) 51-220 µg (stents with diameter 3.5 or 4.0 mm, depending on stent length) |
| * Caution: This device contains Cobalt (CAS N°7440-48-4), classified as CMR 1B, in a concentration above 0.1 % weight per weight. Current scientific evidence supports that medical devices manufactured from alloys containing cobalt do not cause an increased risk of cancer or adverse reproductive effects. (CMR 1B: carcinogenic and toxic to reproduction (CLP regulation EU 1272/2008)) | | |

Caution:

Inform your physician in case of:

- possible allergy to Cobalt-Chromium alloy and nickel
- possible allergy to sirolimus or related compounds, to lactide polymers and caprolactone polymers
- pregnancy/breastfeeding

3. Potential adverse events

Adverse events are listed in alphabetical order. **Please consult your physician for further explanation on the potential adverse events.**

Potential adverse events associated with the procedure of coronary stent placement include but are not limited to:

| | |
|--|---------------------------------------|
| Abrupt vessel closure | Infection and pain at insertion site |
| Acute myocardial infarction | Ischemia, myocardial |
| Allergic reaction to anti-coagulation and/or anti-thrombotic therapy, contrast material, or stent and/or delivery system materials or any other PCI mandatory medication | Myocardial infarction |
| | Nausea and vomiting |
| | Prolonged angina |
| | No reflow |
| | Pseudoaneurysm |
| Aneurysm | Renal failure |
| Arrhythmias, including ventricular fibrillation and ventricular tachycardia | Restenosis of stented segment |
| | Respiratory failure |
| Arteriovenous fistula | Rupture of native and bypass graft |
| Cardiac tamponade | Stent compression |
| Cardiogenic shock | Stent embolization |
| Death | Stent migration |
| Emboli, distal (air, tissue or thrombotic emboli) | Stent thrombosis / occlusion |
| Emergent Coronary Artery Bypass Surgery | Stroke / cerebrovascular accident |
| Failure to deliver the stent to the intended site | Thrombosis (acute, subacute, or late) |
| Fever | Total occlusion of coronary artery |
| Heart failure | Unstable or stable angina pectoris |
| Hematoma | Vessel dissection |
| Hemorrhage, requiring transfusion | Vessel perforation |
| Hypotension / Hypertension | Vessel spasm |

Potential adverse events that may be associated with sirolimus drug are listed below and are consistent for oral intake of sirolimus. In case of stent implantation however, sirolimus is only locally delivered into the vessel wall and amounts of the drug released into the blood are low. Therefore, it is very unlikely that any of the below adverse events (apart from an allergic reaction) associated with oral intake of sirolimus will occur:

| |
|---|
| Abnormal liver function tests |
| Anemia |
| Arthralgias |
| Changes in lipid metabolism which may include hypertriglyceridemia or hypercholesterolemia |
| Diarrhea |
| Hypersensitivity including anaphylactic / anaphylactoid type of reactions |
| Hypokalemia |
| Immune suppression, especially in patients with hepatic insufficiency or who are taking medications that inhibit CYP3A4 or P-glycoprotein |
| Infections |
| Interstitial lung disease |
| Leukopenia |
| Lymphoma and other malignancies |
| Myalgia |
| Thrombocytopenia |

If you experience any of the above listed adverse events, please report it to your treating physician.

4. Symptoms that could indicate that the device is malfunctioning

If you experience any symptom other than what your physician told you before the procedure, immediately seek medical attention. Symptoms of potential device malfunction after procedure may include sudden shortness of breath and chest pain.

5. Medication and examinations after stent placement

After stent placement it is extremely important that you follow the specific medication regimen your physician prescribed. In case you need to discontinue your medication prematurely, please consult your physician.

Follow the physician's instructions for any further follow-up examinations after stent placement.

6. Lifetime of the implanted stent

This device is a permanent implant and is not intended to be removed.

7. MRI (Magnetic Resonance Imaging) information

The Ultimaster Tansei stent has been proven to be MRI conditional. This means that you can safely undergo an MRI scan under certain conditions. When you inform your physician on the type of stent you have (e.g. by showing him/her your implant card), your physician will be able to determine these safe scanning conditions.

8. Passing through metal detectors

The Ultimaster Tansei stent will not set off metal detectors, you can safely pass through metal detectors during security checks (e.g. airport, stores, ...).

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