

# Accessories

## TR Band™

Code number	Description
XX*RF06	TR BAND™ Radial artery compression device

Unit per box: 5

## Haemostatic Valve

Code number	Description	Compatible Guide Wire
RF*VA713M	Haemostatic Valve, Compatible until 7 Fr catheter, with side tube and 3-way stopcock	0.010" / 0.26 mm and more

Unit per box: 5

## Torque Device

Code number	Description	Compatible Guide Wire
XX*RF02M	TORQUE DEVICE for guide wires 0.010" / 0.26 mm to 0.038" / 0.97 mm	0.010" / 0.26 mm to 0.038" / 0.97 mm

Unit per box: 5

## Obturator

- In polypropylene with a rounded tip
- Hub with system of connection to the sheath's hub
- High flexibility and high resistance to kinking
- Prevents clotting reaction inside the sheath
- Color coded

Code number	Description	Compatible sheath size	Length	Color code
XX*RF050410M	OBTURATOR for 4 Fr x 10 cm Sheath	4 Fr	10 cm	Red
XX*RF050510M	OBTURATOR for 5 Fr x10 cm Sheath	5 Fr	10 cm	Grey
XX*RF050525M	OBTURATOR for 5 Fr x 25 cm Sheath	5 Fr	25 cm	Grey
XX*RF050610M	OBTURATOR for 6 Fr x 10 cm Sheath	6 Fr	10 cm	Light green
XX*RF050625M	OBTURATOR for 6 Fr x 25 cm Sheath	6 Fr	25 cm	Light green
XX*RF050710M	OBTURATOR for 7 Fr x 10 cm Sheath	7 Fr	10 cm	Orange
XX*RF050725M	OBTURATOR for 7 Fr x 25 cm Sheath	7 Fr	25 cm	Orange
XX*RF050810M	OBTURATOR for 8 Fr x 10 cm Sheath	8 Fr	10 cm	Light blue
XX*RF050825M	OBTURATOR for 8 Fr x 25 cm Sheath	8 Fr	25 cm	Light blue
XX*RF050910M	OBTURATOR for 9 Fr x 10 cm Sheath	9 Fr	10 cm	Black

Unit per box: 10

## Vessel Dilators

- For predilation of puncture site with Seldinger technique
- Color coded hub

Code number	Description	Length	Outer diameter	Compatible Guide wire	Color code
RF*VD40K10M	4 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.035" guide wire	15.5 cm	4 Fr	0.035" / 0.89 mm	Red
RF*VD50K10M	5 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.035" guide wire	15.5 cm	5 Fr	0.035" / 0.89 mm	Grey
RF*VD50N10M	5 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.038" guide wire	15.5 cm	5 Fr	0.038" / 0.97 mm	Grey
RF*VD60K10M	6 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.035" guide wire	15.5 cm	6 Fr	0.035" / 0.89 mm	Light green
RF*VD60N10M	6 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.038" guide wire	15.5 cm	6 Fr	0.038" / 0.97 mm	Light green
RF*VD70K10M	7 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.035" guide wire	15.5 cm	7 Fr	0.035" / 0.89 mm	Orange
RF*VD70N10M	7 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.038" guide wire	15.5 cm	7 Fr	0.038" / 0.97 mm	Orange
RF*VD80K10M	8 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.035" guide wire	15.5 cm	8 Fr	0.035" / 0.89 mm	Light blue
RF*VD80N10M	8 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.038" guide wire	15.5 cm	8 Fr	0.038" / 0.97 mm	Light blue
RF*VD90K10M	9 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.035" guide wire	15.5 cm	9 Fr	0.035" / 0.89 mm	Black
RF*VD90N10M	9 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.038" guide wire	15.5 cm	9 Fr	0.038" / 0.97 mm	Black
RF*VD10K10M	10 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.035" guide wire	15.5 cm	10 Fr	0.035" / 0.89 mm	Purple
RF*VD10N10M	10 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.038" guide wire	15.5 cm	10 Fr	0.038" / 0.97 mm	Purple
RF*VD11K10M	11 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.035" guide wire	15.5 cm	11 Fr	0.035" / 0.89 mm	Yellow
RF*VD11N10M	11 Fr x 15.5 cm VESSEL DILATOR, compatible with 0.038" guide wire	15.5 cm	11 Fr	0.038" / 0.97 mm	Yellow

Unit per box: 10