

# Radifocus® Guide wire M

## Standard type



### Guide wires

#### Length < 180 cm

**Nitinol hydrophilic guide wires covered with polyurethane and hydrophilic coating for an absolute maneuverability and absence of friction.**

- Extra flexible and non traumatic tapered tip
- Extra hydrophilic thanks to the exclusive “M” polymer coating
- Super elastic alloy core in Nitinol
- One-piece construction for a better push and torque transmission

### Specifications

| Code number | Diameter         | Length | Flexible length | Distal curve              |
|-------------|------------------|--------|-----------------|---------------------------|
| RF-GS18153M | 0.018" / 0.46 mm | 150 cm | 30 mm           | Straight                  |
| RF-GA18153M | 0.018" / 0.46 mm | 150 cm | 30 mm           | Angled                    |
| RF-GS18183M | 0.018" / 0.46 mm | 180 cm | 30 mm           | Straight                  |
| RF-GA18183M | 0.018" / 0.46 mm | 180 cm | 30 mm           | Angled                    |
| RF-GS25123M | 0.025" / 0.64 mm | 120 cm | 30 mm           | Straight                  |
| RF-GA25123M | 0.025" / 0.64 mm | 120 cm | 30 mm           | Angled                    |
| RF-GS25153M | 0.025" / 0.64 mm | 150 cm | 30 mm           | Straight                  |
| RF-GA25153M | 0.025" / 0.64 mm | 150 cm | 30 mm           | Angled                    |
| RF-GS25183M | 0.025" / 0.64 mm | 180 cm | 30 mm           | Straight                  |
| RF-GA25183M | 0.025" / 0.64 mm | 180 cm | 30 mm           | Angled                    |
| RF-GS32123M | 0.032" / 0.81 mm | 120 cm | 30 mm           | Straight                  |
| RF-GA32123M | 0.032" / 0.81 mm | 120 cm | 30 mm           | Angled                    |
| RF-GS32153M | 0.032" / 0.81 mm | 150 cm | 30 mm           | Straight                  |
| RF-GA32153M | 0.032" / 0.81 mm | 150 cm | 30 mm           | Angled                    |
| RF-GS32183M | 0.032" / 0.81 mm | 180 cm | 30 mm           | Straight                  |
| RF-GA32183M | 0.032" / 0.81 mm | 180 cm | 30 mm           | Angled                    |
| RF-GA35053M | 0.035" / 0.89 mm | 50 cm  | 30 mm           | Angled                    |
| RF-GA35083M | 0.035" / 0.89 mm | 80 cm  | 30 mm           | Angled                    |
| RF-GS35123M | 0.035" / 0.89 mm | 120 cm | 30 mm           | Straight                  |
| RF-GA35123M | 0.035" / 0.89 mm | 120 cm | 30 mm           | Angled                    |
| RF-GS35151M | 0.035" / 0.89 mm | 150 cm | 10 mm           | Straight                  |
| RF*GS35153M | 0.035" / 0.89 mm | 150 cm | 30 mm           | Straight                  |
| RF-GS35155M | 0.035" / 0.89 mm | 150 cm | 50 mm           | Straight                  |
| RF-GS35158M | 0.035" / 0.89 mm | 150 cm | 80 mm           | Straight                  |
| RF-GA35151M | 0.035" / 0.89 mm | 150 cm | 10 mm           | Angled                    |
| RF*GA35153M | 0.035" / 0.89 mm | 150 cm | 30 mm           | Angled                    |
| RF-GA35155M | 0.035" / 0.89 mm | 150 cm | 50 mm           | Angled                    |
| RF-GA35158M | 0.035" / 0.89 mm | 150 cm | 80 mm           | Angled                    |
| RF*GR35183M | 0.035" / 0.89 mm | 180 cm | 30 mm           | J Angled (1.5 mm)         |
| RF*GB35153M | 0.035" / 0.89 mm | 150 cm | 30 mm           | J Angled (3 mm)           |
| RF-GS35181M | 0.035" / 0.89 mm | 180 cm | 10 mm           | Straight                  |
| RF-GS35183M | 0.035" / 0.89 mm | 180 cm | 30 mm           | Straight                  |
| RF-GA35181M | 0.035" / 0.89 mm | 180 cm | 10 mm           | Angled                    |
| RF*GA35183M | 0.035" / 0.89 mm | 180 cm | 30 mm           | Angled                    |
| RF*GE35185M | 0.035" / 0.89 mm | 180 cm | 50 mm           | Double Angulation (Bolia) |
| RF*GE35188M | 0.035" / 0.89 mm | 180 cm | 80 mm           | Double Angulation (Bolia) |
| RF-GS38151M | 0.038" / 0.97 mm | 150 cm | 10 mm           | Straight                  |
| RF-GS38153M | 0.038" / 0.97 mm | 150 cm | 30 mm           | Straight                  |
| RF-GA38151M | 0.038" / 0.97 mm | 150 cm | 10 mm           | Angled                    |
| RF*GA38153M | 0.038" / 0.97 mm | 150 cm | 30 mm           | Angled                    |
| RF-GS38183M | 0.038" / 0.97 mm | 180 cm | 30 mm           | Straight                  |
| RF-GA38183M | 0.038" / 0.97 mm | 180 cm | 30 mm           | Angled                    |

Other code numbers are available on special demand. For any further information, please contact your Terumo local representative

**Units per box: 5**

# Radifocus® Guide wire M

## Standard type (For catheter exchange)



### Guide wires

**Length > 220 cm**

**Nitinol hydrophilic guide wires covered with polyurethane and hydrophilic coating for an absolute maneuverability and absence of friction.**

- Extra flexible and non traumatic tapered tip
- Extra hydrophilic thanks to the exclusive “M” polymer coating
- Super elastic alloy core in Nitinol
- One-piece construction for a better push and torque transmission

### Specifications

| Code number | Outer diameter   | Length | Flexible length | Distal curve |
|-------------|------------------|--------|-----------------|--------------|
| RF-GA18223M | 0.018" / 0.46 mm | 220 cm | 30 mm           | Angled       |
| RF-GA18263M | 0.018" / 0.46 mm | 260 cm | 30 mm           | Angled       |
| RF-GS18263M | 0.018" / 0.46 mm | 260 cm | 30 mm           | Straight     |
| RF-GS25263M | 0.025" / 0.64mm  | 260 cm | 30 mm           | Straight     |
| RF-GA25263M | 0.025" / 0.64mm  | 260 cm | 30 mm           | Angled       |
| RF-GS32263M | 0.032" / 0.81 mm | 260 cm | 30 mm           | Straight     |
| RF-GA32263M | 0.032" / 0.81 mm | 260 cm | 30 mm           | Angled       |
| RF-GS35223M | 0.035" / 0.89 mm | 220 cm | 30 mm           | Straight     |
| RF-GA35223M | 0.035" / 0.89 mm | 220 cm | 30 mm           | Angled       |
| RF-GS35261M | 0.035" / 0.89 mm | 260 cm | 10 mm           | Straight     |
| RF-GS35263M | 0.035" / 0.89 mm | 260 cm | 30 mm           | Straight     |
| RF-GA35261M | 0.035" / 0.89 mm | 260 cm | 10 mm           | Angled       |
| RF*GA35263M | 0.035" / 0.89 mm | 260 cm | 30 mm           | Angled       |
| RF-GA35303M | 0.035" / 0.89 mm | 300 cm | 30 mm           | Angled       |
| RF-GA38223M | 0.038" / 0.97 mm | 220 cm | 30 mm           | Angled       |
| RF-GS38263M | 0.038" / 0.97 mm | 260 cm | 30 mm           | Straight     |
| RF-GA38263M | 0.038" / 0.97 mm | 260 cm | 30 mm           | Angled       |

Other code numbers are available on special demand. For any further information, please contact your Terumo local representative

**Units per box: 5**

# Radifocus® Guide wire M

## Stiff type



### Guide wires

**Stiff Nitinol hydrophilic guide wires with excellent torque control and shape memory as well as higher stiffness and resistance.**

- Excellent torque control and superior stiffness with reinforced super elastic metallic core
- Superior smoothness and safety with “M” coating and one-piece construction
- Easy to insert due to its round and flexible tip
- Smooth surface to prevent thrombus formation
- Internal radiopaque structure

### Specifications

| Code number | Type       | Diameter         | Length | Flexible length | Distal curve      |
|-------------|------------|------------------|--------|-----------------|-------------------|
| RF-PA18153M | Stiff      | 0.020" / 0.51 mm | 150 cm | 30 mm           | Angled            |
| RF-PA18263M | Stiff      | 0.020" / 0.51 mm | 260 cm | 30 mm           | Angled            |
| RF-PS25153M | Stiff      | 0.025" / 0.64 mm | 150 cm | 30 mm           | Straight          |
| RF-PA25153M | Stiff      | 0.025" / 0.64 mm | 150 cm | 30 mm           | Angled            |
| RF-PA25183M | Stiff      | 0.025" / 0.64 mm | 180 cm | 30 mm           | Angled            |
| RF-PS25263M | Stiff      | 0.025" / 0.64 mm | 260 cm | 30 mm           | Straight          |
| RF-PA25263M | Stiff      | 0.025" / 0.64 mm | 260 cm | 30 mm           | Angled            |
| RF-PS35083M | Stiff      | 0.035" / 0.89 mm | 80 cm  | 30 mm           | Straight          |
| RF-PA35083M | Stiff      | 0.035" / 0.89 mm | 80 cm  | 30 mm           | Angled            |
| RF-PS35153M | Stiff      | 0.035" / 0.89 mm | 150 cm | 30 mm           | Straight          |
| RF-PA35151M | Stiff      | 0.035" / 0.89 mm | 150 cm | 10 mm           | Angled            |
| RF*PA35153M | Stiff      | 0.035" / 0.89 mm | 150 cm | 30 mm           | Angled            |
| RF-PS35183M | Stiff      | 0.035" / 0.89 mm | 180 cm | 30 mm           | Straight          |
| RF*PA35183M | Stiff      | 0.035" / 0.89 mm | 180 cm | 30 mm           | Angled            |
| RF*PA35188M | Stiff      | 0.035" / 0.89 mm | 180 cm | 80 mm           | Angled            |
| RF*HR35183M | Half stiff | 0.035" / 0.89 mm | 180 cm | 30 mm           | J Angled (1.5 mm) |
| RF*HR35303M | Half stiff | 0.035" / 0.89 mm | 300 cm | 30 mm           | J Angled (1.5 mm) |
| RF*HN35183M | Half stiff | 0.035" / 0.89 mm | 180 cm | 30 mm           | J curve           |
| RF-PS35263M | Stiff      | 0.035" / 0.89 mm | 260 cm | 30 mm           | Straight          |
| RF-PA35263M | Stiff      | 0.035" / 0.89 mm | 260 cm | 30 mm           | Angled            |
| RF-PS38083M | Stiff      | 0.038" / 0.97 mm | 80 cm  | 30 mm           | Straight          |
| RF-PS38153M | Stiff      | 0.038" / 0.97 mm | 150 cm | 30 mm           | Straight          |
| RF-PA38153M | Stiff      | 0.038" / 0.97 mm | 150 cm | 30 mm           | Angled            |
| RF-PS38183M | Stiff      | 0.038" / 0.97 mm | 180 cm | 30 mm           | Straight          |
| RF-PA38183M | Stiff      | 0.038" / 0.97 mm | 180 cm | 30 mm           | Angled            |
| RF-PA38263M | Stiff      | 0.038" / 0.97 mm | 260 cm | 30 mm           | Angled            |

Other code numbers are available on special demand. For any further information, please contact your Terumo local representative

**Unit per box: 5**

# Radifocus® Guide wire M

## Shapeable type

### Guide wires



Nitinol hydrophilic guide wires with high torque control and shapeable tip for special sites particularly difficult to reach.

- Tapered pre-formable tip to reach difficult access sites
- Hydrophilic “M” Coating and one-piece construction for superior sliding and easy control
- Super elastic metallic core made of Nitinol for an excellent torque control
- Smooth surface to prevent thrombus formation
- Internal radiopaque structure

### Specifications

| Code number | Diameter         | Length | Flexible Length | Distal curve |
|-------------|------------------|--------|-----------------|--------------|
| RE*GL35153M | 0.035" / 0.89 mm | 150 cm | 30 mm           | Straight     |
| RE*GK35153M | 0.035" / 0.89 mm | 150 cm | 30 mm           | Angled       |
| RE*GL35183M | 0.035" / 0.89 mm | 180 cm | 30 mm           | Straight     |
| RE*GK35183M | 0.035" / 0.89 mm | 180 cm | 30 mm           | Angled       |

Other code numbers are available on special demand. For any further information, please contact your Terumo local representative

**Units per box: 5**

# Radifocus® Guide wire M with Radiopaque tip marker



## Guide wires

**Nitinol hydrophilic guide wires with small diameter, high torque control and highly visible distal gold marker tip for micro vessels access.**

- Small diameters (0.011" / 0.25 mm, 0.016" / 0.41 mm and 0.018" / 0.46 mm) specially adapted for micro vessels access
- Hydrophilic "M" Coating for superior sliding
- One-piece construction for navigability
- High visibility of the tip due to the distal gold marker
- Highly flexible for an easy access to small and distal micro vessels

## Specifications

| Code number  | Diameter         | Length | Flexible length | Distal curve | Distal marker |
|--------------|------------------|--------|-----------------|--------------|---------------|
| RF*RA10188TM | 0.011" / 0.25 mm | 180 cm | 80 mm           | Angled 45°   | Gold marker   |
| RF*RB10188TM | 0.011" / 0.25 mm | 180 cm | 80 mm           | Angled 70°   | Gold marker   |
| RF*RA16188TM | 0.016" / 0.41 mm | 180 cm | 80 mm           | Angled 45°   | Gold marker   |
| RF*RB16188TM | 0.016" / 0.41 mm | 180 cm | 80 mm           | Angled 70°   | Gold marker   |
| RF*RA18188TM | 0.018" / 0.46 mm | 180 cm | 80 mm           | Angled 45°   | Gold marker   |
| RF*RB18188TM | 0.018" / 0.46 mm | 180 cm | 80 mm           | Angled 70°   | Gold marker   |

Other code numbers are available on special demand. For any further information, please contact your Terumo local representative

**Unit per box: 1**