



Winged Blood Sampling Set with pre-attached Holder

User-friendliness in a small package



**The ready-to-use solution
for reliable phlebotomy**

- Solid, hard blister packaging to secure sterility and integrity
- No kinking of tubing, tight and proof connection of luer adapter
- Easy opening for accurate preparation before use
- Broad range for all major applications and demands

 **TERUMO**[®]

We keep life flowing

Available in a tight-fitting, hard blister packaging

The Terumo's winged blood sampling set with pre-attached holder is individually wrapped in a solid hard blister. Its design allows you to check the integrity of the product before use.

The strong fixation of the needle protector secures the well known sharpness of the needle tip and the functionality of the safety shield - during transport and while opening.



No kinking of tubing, or disconnection of winged set from holder

The hard blister contains a large chamber in which the tubing of the winged set is placed - thus keeping it as straight as possible, in order to avoid any kinking. Disconnection - even partially - of the luer adapter and luer cone is avoided by a special fixation feature inside of the blister.

Easy opening, accurate preparation

The large knobs on the blister, allow for an easy opening by peeling the blister paper off in a soft, smooth motion.

Once opened, you can easily grasp the unit by the holder as additional space is foreseen to increase your comfort.

The blister design will prevent the product from falling out of its packaging while opening it.



Winged Blood Sampling Sets with pre-attached Holder

CE 0197

Surshield™ Safety Winged Blood Sampling Set with Holder					
Code		Diameter of Needle	Length Tubing	Safety	Packaging
MN-SVS21H18		21G	180 mm	Yes	12 x 40
MN-SVS21H30		21G	300 mm	Yes	12 x 40
MN-SVS23H18		23G	180 mm	Yes	12 x 40
MN-SVS23H30		23G	300 mm	Yes	12 x 40

CE 0197

Conventional Winged Blood Sampling Set with Holder					
Code		Diameter of Needle	Length Tubing	Safety	Packaging
MN-SV21H18		21G	180 mm	No	12 x 40
MN-SV21H30		21G	300 mm	No	12 x 40
MN-SV23H18		23G	180 mm	No	12 x 40
MN-SV23H30		23G	300 mm	No	12 x 40

Every day, as healthcare professionals, you are working to ensure life and to contribute to the well-being of your patients. Every day, Terumo Europe provides you with the clinical systems, means and solutions to keep life flowing.



We keep life flowing

TERUMO CORPORATION
44-1, 2-chome Hatagaya Shibuya-ku Tokyo 151-0072 JAPAN
Tel. +81/3 3374.8111 Fax +81/3 3374.8196

TERUMO EUROPE N.V.
Researchpark Zone 2 Haasrode Interleuvenlaan 40 B-3001 Leuven BELGIUM
Tel. +32/16.38.12.11 Fax +32/16.40.02.49

TERUMO EUROPE N.V. BENELUX SALES DIVISION
Researchpark Zone 2 Haasrode Interleuvenlaan 40 B-3001 Leuven BELGIUM
Tel. +32/16.39.25.80 Fax +32/16.39.25.99 THE NETHERLANDS
Tel. 0800/022.03.96 Fax 0800/022.04.14

LABORATOIRES TERUMO FRANCE S.A.
1, Parc Ariane Bâtiment Uranus Rue Hélène Boucher 78284 Guyancourt Cedex FRANCE
Tel. +33/1 30.96.13.00 Fax +33/1 30.43.60.85

TERUMO DEUTSCHLAND GmbH
Hauptstrasse 87 65760 Eschborn GERMANY
Tel. +49/61.96.80.230 Fax +49/61.96.80.23.200

TERUMO DEUTSCHLAND GmbH Zweigniederlassung
Wettingen Bahnhofstrasse 88 CH-5430 Wettingen SWITZERLAND
Tel. +41/56.419.10.10 Fax +41/56.419.10.11

TERUMO EUROPE ESPAÑA SL
Avda. Juan Carlos I, N° 13 - 7ª Planta Edificio Torre La Garena 28806 Alcalá de Henares (Madrid) SPAIN
Tel. general: +34/902.10.12.98 Tel. cust. serv.: +34/902.20.31.40 Fax general: +34/902.10.13.58 Fax cust. serv.: +34/902.10.17.57

TERUMO ITALIA SRL
Via Simone Martini 143/145 I-00142 Roma ITALY
Tel. +39/06.519.61.420 Fax +39/06.503.04.07

TERUMO SWEDEN AB
Nya Varvet, Byggnad 90 Sven Källfelts Gata 18 SE-426 71 Västra Frölunda SWEDEN
Tel. +46/31.748.58.80 Fax +46/31.748.58.90

www.terumo.com
www.terumo-europe.com

TERUMO DENMARK Filial af TERUMO SWEDEN AB
Tel. +45/7020.93.80 Fax +45/7020.94.80

TERUMO UK Ltd.
Tamesis The Causeway Egham Surrey TW20 9AW UNITED KINGDOM
Tel. +44/1784.476.200 Fax +44/1784.476.201

TERUMO EUROPE N.V. Moscow Representative office
4, Ivana Franko St 121108 Moscow RUSSIA
Tel. +7/495.234.3303 Fax. +7/495.234.3303

®: Registered Trademark Europe N.V.

Addresses are subject to change, please consult our website.