


1961 Type Aortic Perfusion Cannulae




1963 Type Aortic Perfusion Cannulae

1961 Type Aortic Perfusion Cannulae

▶ Ordering Information

PRESSURE DROP VS. FLOW RATE 1961 TYPE AORTIC PERFUSION CANNULAE
Terumo in vitro evaluation. Data obtained with water.



DESCRIPTION	CATALOG NO.	UNITS/CASE
Wire reinforced, multi-notched suture ring		
24 Fr, 25 cm (9") long	L8267	10

✓ Features


1961 Type Aortic Perfusion Cannulae

- Smooth, angled tip
- Multi-notched suture ring
- Wire reinforced

1963 Type Aortic Perfusion Cannulae

▶ Ordering Information

PRESSURE DROP VS. FLOW RATE 1963 TYPE AORTIC PERFUSION CANNULAE
Terumo in vitro evaluation. Data obtained with water.



DESCRIPTION	CATALOG NO.	UNITS/CASE
Wire reinforced, short angled tip		
22 Fr, 27 cm (11") long	L7710	10
24 Fr, 27 cm (11") long	L7711	10

✓ Features

1963 Type Aortic Perfusion Cannulae

- Short, smooth tip
- Wire reinforced
- Retention ring




1965 Type Aortic Perfusion Cannulae

1965 Type Aortic Perfusion Cannulae

PRESSURE DROP VS. FLOW RATE 1965 TYPE AORTIC PERFUSION CANNULAE

Terumo in vitro evaluation. Data obtained with water.



▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
Wire reinforced, angled tip, 3 side holes		
22 Fr, 25 cm (9") long	L7185	10
24 Fr, 25 cm (9") long	L7186	10
26 Fr, 25 cm (9") long	L7187	10

✓ Features

1965 Type Aortic Perfusion Cannulae


- Anti-kink wire reinforced tubing
- Smooth blood flow path
- Suture ring
- Side hole exit ports



Aortic Perfusion Cannulae

PRESSURE DROP VS. FLOW RATE AORTIC PERFUSION CANNULAE

Terumo in vitro evaluation. Data obtained with water.



Aortic Perfusion Cannulae

▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
10 Fr, straight taper, 5.7" (14.5 cm) long	L211010	25
14 Fr, straight taper, 5.5" (14.0 cm) long	L211014	25
18 Fr, straight taper, 5.3" (13.4 cm) long	L211018	25
24 Fr, straight taper, 15.0" (38.1 cm) long	L211024	25
24 Fr, straight taper, blue stripe, 9.0" (22.9 cm) long	L211224	25
24 Fr, straight taper, blue stripe, 3/8" connector, 6.9" (17.5 cm) long	L211324	25

✓ Features

Aortic Perfusion Cannulae

- Smooth tapered tip and tube
- Retention ring
- Compact length