




Sarns™ Flexible Arterial Cannulae


To perfuse the ascending aorta or femoral artery

PRESSURE DROP VS. FLOW RATE SARNS™ FLEXIBLE ARTERIAL CANNULAE

Terumo in vitro evaluation. Data obtained with water.



Sarns™ Flexible Arterial Cannulae



▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
6.5 mm (20 Fr) OD with 3/8" connector, suture ring, 9.5" (24 cm) long	13010	20
7.5 mm (22 Fr) OD with 3/8" connector, suture ring, 9.5" (24 cm) long	13020	20
8.0 mm (24 Fr) OD with 3/8" connector, suture ring, 9.5" (24 cm) long	13030	20
8.0 mm (24 Fr) OD with cm markings, 10.0" (25 cm) long	144776	20

✓ Features

Sarns™ Flexible Arterial Cannulae

- Smooth, beveled tip
- Movable suture ring
- Facilitates aortic or femoral cannulation