



Terumo® AL8X Arterial Blood Filter

Pall's proven filtration technology with Terumo's X Coating™

SPECIFICATIONS	AL8X
Housing material	Polycarbonate
Filter material	Polyester
Filter	40 microns
Blood ports	3/8" (10 mm)
Prime volume	170 mL
Max. blood flow rate	8.0 L/min
Pressure drop (blood flow)	22 mmHg (6.0 L/min)

▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
Arterial blood filter, 40 microns, sterile	AL8X	10
 Holders		
Holder for AL8, AL8X	FH8J	1



Terumo® AL8X
Arterial Blood Filter