



## CAPIOX<sup>®</sup> Bubble Traps

Unique design for effective air removal for adult and pediatric use

SPECIFICATIONS	BT15	BT05
Intended use	In adult extracorporeal circuits (QB < 6.5 L/min)	In neonatal extracorporeal circuits (QB < 1.6 L/min)
Blood port	3/8"	1/4"
Vent port	Female luer	Female luer
Filter material	Polyester	Polyester
Filter pore size	170 microns	170 microns
Filter pore area	60%	60%
Housing	Polycarbonate	Polycarbonate
Prime volume	150 mL	50 mL



CAPIOX<sup>®</sup> BT15 Bubble Trap, adult



CAPIOX<sup>®</sup> BT05 Bubble Trap, pediatric

### ▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
Pediatric bubble trap .....	CX*BT05 .....	10
Adult bubble trap .....	CX*BT15 .....	10
Bubble trap holder for BT05 and BT15 .....	XX*XH041 .....	1

### ✓ Features

#### CAPIOX<sup>®</sup> Bubble Traps

- Easy priming
- Efficient air removal
- Minimal blood trauma
- Low prime volume