



CAPIOX[®] Bubble Traps

Unique design for effective air removal for adult and pediatric use

SPECIFICATIONS

BT05

Intended use	In neonatal extracorporeal circuits (QB < 1.6 L/min)
Blood port	1/4"
Vent port	Female luer
Filter material	Polyester
Filter pore size	170 microns
Filter pore area	60%
Housing	Polycarbonate
Prime volume	50 mL



CAPIOX[®] BT05 Bubble Trap

▶ Ordering Information

DESCRIPTION

CATALOG NO.

UNITS/CASE

Pediatric	CX*BT05	10
Bubble trap holder for BT05	XX*XH041	1



Features

CAPIOX[®] Bubble Traps

- Easy priming
- Efficient air removal
- Minimal blood trauma
- Low prime volume