



# CAPIOX<sup>®</sup> Flexible Venous Reservoirs

Flexible reservoirs with unique venous inlet tube design for maximum air removal with minimal blood stagnation for adult and pediatric use

### SPECIFICATIONS

Reservoir size	Inlet port size	CR port size	Outlet port size
100 mL*	1/4"	1/4"	1/4"
400 mL*	3/8"	3/8"	1/4"

## ▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
100 mL* .....	CX*VRA0101T2 .....	10
400 mL* .....	CX*VRA0401T2 .....	4
Holder for 100 mL* reservoir .....	XX*XH011 .....	1

\*Double luer port with detachable temperature probe



### Features

CAPIOX<sup>®</sup> Flexible Venous Reservoirs

- Allows maximum air removal
- Minimizes blood stagnation
- Minimizes blood damage with simple design