



## CAPIOX® Bubble Traps

Unique design for effective air removal for pediatric use

### SPECIFICATIONS

### BT05

Intended use	In neonatal extracorporeal circuits (QB < 1.6 L/min)
Blood port	1/4"
Vent port	Female luer
Filter material	Polyester
Filter pore size	170 microns
Filter pore area	60%
Housing	Polycarbonate
Prime volume	50 mL



CAPIOX® BT05 Bubble Trap

### ▶ Ordering Information

#### DESCRIPTION

#### CATALOG NO. UNITS/CASE

Pediatric bubble trap .....	CX*BT05 .....	10
Bubble trap holder for BT05 .....	XX*XH041 .....	1

### ✓ Features

#### CAPIOX® Bubble Traps

- Easy priming
- Efficient air removal
- Minimal blood trauma
- Low prime volume