



CAPIOX® Cardiotomy Reservoirs

Cardiotomy reservoirs for high-efficiency filtration, defoaming, and removal of foreign particles

SPECIFICATIONS	CR10 NX
Housing material	Polycarbonate
Reservoir capacity	1000 mL
Maximum blood flow	1.5 L/min
Filter material	Polyester
Filter efficiency	Removes 90% or more of particles with ≥ 20 micron diameter
Defoamer material	Polyurethane foam
Blood inlet ports	Five suction ports: 3/16" – 1/4" rotatable
Prime port	1/4" filtered
Luer ports	Four filtered One non-filtered
Vent port	1/4"
Blood outlet port	1/4"

▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
Ⓞ Cardiotomy reservoir with filter and X Coating™, 1000 mL	CX*CR10NX	10
Holders		
Ⓞ Holder for CR10NX	XX*CXH05R	1
Adapter for SX holder for CR10NX	XX*CXH05AD	1



CAPIOX® CR10NX Cardiotomy Reservoir — 1000 mL, X Coating™



Holder for CAPIOX® CR10NX

*X*coating™ is Terumo's own biocompatible hydrophilic polymer surface coating

Ⓞ standard item – delivery from stock