



CAPIOX[®] Flexible Venous Reservoirs

Flexible reservoirs with unique venous inlet tube design for maximum air removal with minimal blood stagnation for pediatric use

SPECIFICATIONS

Reservoir size	Inlet port size	CR port size	Outlet port size
100 mL*	1/4"	1/4"	1/4"
400 mL*	3/8"	3/8"	1/4"

▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
100 mL*	CX*VRA0101T2	10
400 mL*	CX*VRA0401T2	4
Holder for 100 mL reservoir	XX*XH011	1

*Double luer port with detachable temperature probe

✓ Features

CAPIOX[®] Flexible Venous Reservoirs

- Allows maximum air removal
- Minimizes blood stagnation
- Minimizes blood damage with simple design