


Tender Flow™ Pediatric Venous Return Cannulae

Developed for pediatric surgeons and their smallest patients

PRESSURE DROP VS. FLOW RATE TENDER FLOW™ PEDIATRIC VENOUS RETURN CANNULAE
Terumo in vitro evaluation. Data obtained with water.



Tender Flow™ Pediatric Venous Return Cannulae

Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
Straight: 12" (30.5 cm) long, bent: 11.5" (29.2 cm) long		
2.7 mm (8 Fr) with 1/4" flare, WR	813541	10
2.7 mm (8 Fr) with 1/4" flare, WR, bent tip	813542	10
3.3 mm (10 Fr) with 1/4" flare, WR	813543	10
3.3 mm (10 Fr) with 1/4" flare, WR, bent tip	813544	10
4.0 mm (12 Fr) with 1/4" flare, WR	813545	10
4.0 mm (12 Fr) with 1/4" flare, WR, bent tip	813546	10
4.7 mm (14 Fr) with 1/4" flare, WR	813547	10
4.7 mm (14 Fr) with 1/4" flare, WR, bent tip	813548	10
5.3 mm (16 Fr) with 1/4" flare, WR	813549	10
5.3 mm (16 Fr) with 1/4" flare, WR, bent tip	813550	10
6.0 mm (18 Fr) with 1/4" flare, WR	813551	10
6.0 mm (18 Fr) with 1/4" flare, WR, bent tip	813552	10
6.0 mm (18 Fr) with 3/8" flare, WR	813553	10
6.0 mm (18 Fr) with 3/8" flare, WR, bent tip	813554	10
6.7 mm (20 Fr) with 1/4" flare, WR	813555	10
6.7 mm (20 Fr) with 1/4" flare, WR, bent tip	813556	10
6.7 mm (20 Fr) with 3/8" flare, WR	813557	10
6.7 mm (20 Fr) with 3/8" flare, WR, bent tip	813558	10
7.3 mm (22 Fr) with 1/4" flare, WR	813559	10
7.3 mm (22 Fr) with 1/4" flare, WR, bent tip	813560	10
7.3 mm (22 Fr) with 3/8" flare, WR	813561	10
7.3 mm (22 Fr) with 3/8" flare, WR, bent tip	813562	10
8.0 mm (24 Fr) with 3/8" flare, WR	813565	10
8.0 mm (24 Fr) with 3/8" flare, WR, bent tip	813566	10

Features

- Tender Flow™ Pediatric Venous Return Cannulae
- Wide range of sizes
- Step-down design
- Straight and bent tips

WR: wire reinforced

Ⓢ standard item – delivery from stock