




Terumo® Modular Perfusion Circuit


ROCsafeRX – Increased performance, reduced prime volume, ultimate safety

▶ Ordering Information

DESCRIPTION	CATALOG NO.	UNITS/CASE
Holders for ROCsafe and ROCsafe RX		
Holder for ROCsafe modular perfusion circuit	XX*CXH18.....	1
Extension pole mount bracket for ROCsafe holder	814858.....	1
Holders for ROCsafe accessory kits		
Holder for accessory pack with flexible reservoir CX-ROCA	816181.....	1
Holder for accessory pack with flexible reservoir CX-ROCA, short arm	816173.....	1
Holder for accessory pack with hardshell reservoir CX-ROCB	801804.....	1
Holder for accessory pack with hardshell reservoir CX-ROCB, short arm	801139.....	1



Order # 816173



Order # 816181



Order # XX*CXH18




Order # 814858



Order # 801804




Order # 801139



✓ Features

Accessory kit with flexible reservoir

- Can be connected to ROCsafe during operation using quick connectors
- For storage and de-airing of large amounts of vented blood
- Intended to be used during complex procedures like valve replacements or combined cases



✓ Features

Accessory kit with hardshell reservoir

- Can be connected to ROCsafe during operation using quick connectors
- Converts ROCsafe to an open circuit
- Intended to be used in emergency situations