


Specifications



Ordering Information

Model number	Insertion profile O.D. (actual)	Expanded sheath O.D.	Expanded sheath I.D.	Re-collapsed sheath O.D.	Working length	Expandable length	Tip outer diameter
SR-1925	13.5 Fr	23 Fr 7.67 mm	19 Fr 6.33 mm	~30% reduction in profile ²	25 cm	20 cm	5.3 Fr
SR-1935	13.5 Fr	23 Fr 7.67 mm	19 Fr 6.33 mm	~30% reduction in profile ²	35 cm	30 cm	5.3 Fr
SR-2225	15 Fr	26 Fr 8.67 mm	22 Fr 7.33 mm	~30% reduction in profile ²	25 cm	20 cm	5.3 Fr
SR-2235	15 Fr	26 Fr 8.67 mm	22 Fr 7.33 mm	~30% reduction in profile ²	35 cm	30 cm	5.3 Fr
SR-2425	15 Fr	28.5 Fr 9.50 mm	24 Fr 8.00 mm	~30% reduction in profile ²	25 cm	20 cm	5.3 Fr
SR-2435	15 Fr	28.5 Fr 9.50 mm	24 Fr 8.00 mm	~30% reduction in profile ²	35 cm	30 cm	5.3 Fr

2. Data on file.

Solopath™ is a trademark of Onset Medical Corporation



Terumo Corporation
44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-0072,
Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8196

Terumo Europe NV
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

Terumo Europe NV
Interventional Systems
Liaison Office
Immeuble Inside
7-9 Rue Eugène et Armand Peugeot
92500 Rueil Malmaison,
France
Tel.: +33 147 16 09 30
Fax: +33 147 16 09 49

Terumo Europe NV
Benelux Sales Division
Belgium:
Tel.: +32 16 39 25 80
The Netherlands:
Tel.: 0800 0220396

Terumo Europe NV
Emerging Markets Division
Tel.: +32 16 38 13 08

Terumo Deutschland GmbH
Tel.: +49 6196 80 230

Terumo Deutschland GmbH
Zweigniederlassung Spreitenbach
Tel.: +41 56 419 10 10

Terumo Deutschland GmbH
Zweigniederlassung Österreich
Interventional Systems
Tel.: +43 2236 378020

Laboratoires Terumo France S.A.S.
Tel.: +33 130 96 13 00

Terumo Italia S.r.l.
Tel.: +39 0651 96 14 20

Terumo Europe España SL
Tel.: +34 9021 01 298

Terumo UK Ltd
Tel.: +44 1276 480440

Terumo Russia LLC
Tel.: +7 495 988 4740

Terumo Sweden AB
Tel.: +46 3174 85 880

Terumo Denmark
Filial of Terumo Sweden AB
Tel.: +45 7020 93 80

Terumo Corporation Dubai Branch
(Middle East and Africa)
Tel.: +971 4 221 2220

Terumo BCT Tıbbi Cihazlar Dağıtım
ve Hizmetleri A.Ş.
Tel.: +90 216 590 0074






© Registered Trademark
Published by Terumo Europe NV

15465GB0314PFRU1214E

SoloPath™ RE-COLLAPSIBLE

Re-collapsible Balloon Access System

Fits right in...
...and right out




The first and only **expandable**
& **re-collapsible** femoral access sheath.


Finally a sheath that goes in *easy* & comes out *easier* than ever before.¹

Introducing the first and only femoral access sheath for large-bore access that lets you:

- Insert at a low profile
- Expand to a predictable operating profile
- Collapse and remove at a low profile



Fits right in




Smallest Insertion Profile

Insertion Profile
(19 Fr Introducer)




Optimal operating profile




Large Operating Profile

Operating Profile
(19 Fr Introducer)




Fits right out




Smallest Withdrawal Profile

Withdrawal Profile
(19 Fr Introducer)




Lowest Insertion Force¹

24 Fr Insertion Force Comparisons



Lowest Withdrawal Force¹

24 Fr Total Retraction Force Comparisons



¹ Data on file. IDTE Force Testing TIS2013-005.