




Fits right in...



...and right out.

The first and only **expandable & re-collapsible** femoral access sheath.

SoloPath[®] RE-COLLAPSIBLE

Re-collapsible Balloon Access System

Finally—a sheath that goes in **easy**
and comes out **easier** than ever before.¹

**Introducing the first and only femoral access sheath
for large-bore access that lets you:**


- Insert at a low profile
- Expand to a predictable operating profile
- Collapse and remove at a low profile



SoloPath® RE-COLLAPSIBLE


Re-collapsible Balloon Access System

Fits right in...




Smallest Insertion Profile

Insertion Profile
(19 Fr Introducer)



Lowest Insertion Force¹


24 Fr Insertion Force Comparisons



SoloPath® RE-COLLAPSIBLE

Re-collapsible Balloon Access System

Optimal operating profile...



Large Operating Profile



SoloPath® RE-COLLAPSIBLE

Re-collapsible Balloon Access System


Fits right out.

External balloon enables controlled collapsation for sheath removal at lower profile

Hydrophilic coating to minimize force of withdrawal resistance
– Lower-force withdrawal reduces potential for vessel avulsion


Smallest Withdrawal Profile

Withdrawal Profile
(19 Fr Introducer)



Lowest Withdrawal Force¹

24 Fr Total Retraction Force Comparisons



SoloPath® RE-COLLAPSIBLE

Re-collapsible Balloon Access System

SOLOPATH® Re-collapsible Access System

Ordering Information

MODEL NUMBER	INSERTION PROFILE O.D. (ACTUAL)	EXPANDED SHEATH O.D.	EXPANDED SHEATH I.D.	RE-COLLAPSED SHEATH O.D.	WORKING LENGTH	EXPANDABLE LENGTH	TIP DIAMETER
SR-1925	13.5 Fr	23 Fr 7.67mm	19 Fr 6.33mm	~30% reduction in profile ²	25 cm	20 cm	5.3 Fr
SR-1935	13.5 Fr	23 Fr 7.67mm	19 Fr 6.33mm	~30% reduction in profile ²	35 cm	30 cm	5.3 Fr
SR-2225	15 Fr	26 Fr 8.67mm	22 Fr 7.33mm	~30% reduction in profile ²	25 cm	20 cm	5.3 Fr
SR-2235	15 Fr	26 Fr 8.67mm	22 Fr 7.33mm	~30% reduction in profile ²	35 cm	30 cm	5.3 Fr
SR-2425	15 Fr	28.5 Fr 9.50mm	24 Fr 8.00mm	~30% reduction in profile ²	25 cm	20 cm	5.3 Fr
SR-2435	15 Fr	28.5 Fr 9.50mm	24 Fr 8.00mm	~30% reduction in profile ²	35 cm	30 cm	5.3 Fr

SoloPath® RE-COLLAPSIBLE

Re-collapsible Balloon Access System

For Rx only.

Before using refer to **Instructions for Use** for indications, contraindications as well as warnings and precautions @ www.terumo.com

*All other brand names are registered trademarks of their respective owners.

1. Data on file. IDTE Force Testing TIS2013-005. **2.** Data on file.