

Tercross™

0.014" OTW PTA Balloon Dilatation Catheter

A workhorse PTA Balloon
dedicated for Btk* vessels

Outstanding crossability

- ▶ Lowest balloon crossing profile, hydrophilic coated, provides reliable crossability
- ▶ Extreme support, OTW**, seamless, polymer shaft offers enhanced pushability, adequate kink resistance and fast deflation times

Extensive pressure working range

- ▶ The balloon can be inflated up to a Rated Burst Pressure (RBP) of 20 atm (2026 kPa) for efficient and successful recanalization

Wide size mix

- ▶ A wide choice of balloon diameters (Ø 1.25 - 4.0 mm) and lengths (↔20 - 200 mm) cover all clinical needs




Peripheral Vascular Solutions

*Btk= Below the Knee
**OTW= Over the Wire

TERUMO®

Tercross™



Specifications

Catheter design:	Over the Wire	Coating:	hydrophilic
Guide wire compatibility:	0.014" / 0.36 mm	Shaft diameter: proximal:	Ø1.25 mm, 1.5 mm: 3.2 Fr (1.07 mm)
Usable shaft length:	100 cm, 148 cm		Ø2.0-4.0 mm: 3.6 Fr (1.21 mm)
Balloon folding:	2 folded (Ø1.25 mm, 1.5 mm)	distal:	Ø1.25 mm, 1.5 mm: 2.5 Fr (0.83 mm)
	3 folded (Ø2.0-4.0 mm)		Ø2.0-4.0 mm: 3.0 Fr (1.00 mm)

Ordering information

Shaft length: 100 cm

Balloon diameter	Balloon length / Product code					
	20 mm	40 mm	80 mm	120 mm	150 mm	200 mm
1.25	BD-T1220P4E	-	-	-	-	-
1.50	BD-T1520P4E	-	-	-	-	-
2.00	-	BD-T2040P4E	BD-T2080P4E	BD-T20120P4E	BD-T20150P4E	BD-T20200P4E
2.50	-	BD-T2540P4E	BD-T2580P4E	BD-T25120P4E	BD-T25150P4E	BD-T25200P4E
3.00	-	BD-T3040P4E	BD-T3080P4E	BD-T30120P4E	BD-T30150P4E	BD-T30200P4E
3.50	-	BD-T3540P4E	BD-T3580P4E	BD-T35120P4E	BD-T35150P4E	BD-T35200P4E
4.00	-	BD-T4040P4E	BD-T4080P4E	BD-T40120P4E	BD-T40150P4E	BD-T40200P4E

Shaft length: 148 cm

Balloon diameter	Balloon length / Product code					
	20 mm	40 mm	80 mm	120 mm	150 mm	200 mm
1.25	BD-T1220Q4E	-	-	-	-	-
1.50	BD-T1520Q4E	-	-	-	-	-
2.00	-	BD-T2040Q4E	BD-T2080Q4E	BD-T20120Q4E	BD-T20150Q4E	BD-T20200Q4E
2.50	-	BD-T2540Q4E	BD-T2580Q4E	BD-T25120Q4E	BD-T25150Q4E	BD-T25200Q4E
3.00	-	BD-T3040Q4E	BD-T3080Q4E	BD-T30120Q4E	BD-T30150Q4E	BD-T30200Q4E
3.50	-	BD-T3540Q4E	BD-T3580Q4E	BD-T35120Q4E	BD-T35150Q4E	BD-T35200Q4E
4.00	-	BD-T4040Q4E	BD-T4080Q4E	BD-T40120Q4E	BD-T40150Q4E	BD-T40200Q4E



Terumo Corporation
44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-0072, Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8196

Terumo Europe N.V.
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

**Terumo Europe
Interventional Systems
Liaison Office**
Immeuble Krebs
8, rue Jean-Pierre Timbaud
CS 10720
78067 Saint Quentin
en Yvelines Cedex, France
Tel.: +33 139 30 95 75
Fax: +33 139 30 95 70

**Terumo Europe N.V.
Benelux Sales Division**
Belgium:
Tel.: +32 16 39 25 80
The Netherlands:
Tel.: 0800 0220396

**Terumo Europe N.V.
Emerging Markets Division**
Tel.: +32 16 38 13 08

Terumo Deutschland GmbH
Tel.: +49 6196 80 230

**Terumo Deutschland GmbH
Zweigniederlassung Spreitenbach**
Tel.: +41 56 419 10 10

**Terumo Deutschland GmbH
Zweigniederlassung Österreich
Interventional Systems**
Tel.: +43 2236 378020

Laboratoires Terumo France S.A.
Tel.: +33 130 96 13 00

Terumo Italia S.r.l.
Tel.: +39 0651 96 14 20

Terumo Europe España SL
Tel.: +34 9021 01 298

Terumo UK Ltd
Tel.: +44 1784 476 200

Terumo Russia LLC
Tel.: +7 495 988 4740

Terumo Sweden AB
Tel.: +46 3174 85 880

**Terumo Denmark
Filial of Terumo Sweden AB**
Tel.: +45 7020 93 80

**Terumo Corporation Dubai Branch
(Middle East and Africa)**
Tel.: +971 4 221 2220



* Registered Trademark
Published by Terumo Europe N.V.