








Terumo[®] Surshield[®] Versatus[™]

**Straight I.V. Catheters designed
to prevent sharps injuries**



Surshield[®] sharps protection. Our point? Your safety!

-  Passive activation of the protection system
-  Full and immediate encapsulation of the needle point at removal
-  No exposure to contaminated needle point
-  Permanent protection after use
-  Unchanged insertion technique



 **TERUMO[®]**

We keep life flowing



TERUMO's Surshield®, an integrated sharps protection feature¹, on TERUMO's I.V. catheters, represents an important step ahead in reducing risks² for hospital personnel. Healthcare workers will appreciate the easy to use protection that equips Surshield® Versatus™ I.V. catheters. The Versatus™ I.V. catheters with Surshield® sharps protection, set new safety standards in sharps injury prevention. The comfort and efficiency, for which TERUMO products are widely appreciated by both medical staff and patients, are fully preserved.

Polyurethane I.V. catheters

The radio-opaque catheter is made of polyurethane (PUR) which is firm at room temperature and becomes softer at body temperature. The comfort of your patients will be enhanced as the catheter is kink resistant and biocompatible. The colour coding of the hub allows you to easily recognize the catheter size.



Conical catheter tip

The perfectly tapered catheter tip allows the catheter to enter the vein minimizing venous trauma and allows for an easy placement. The bevel/catheter tip distance reduces the risk of penetrating the vein.

Flashback

The flashback chamber is transparent for easy detection of the flashback. A venting plug, having an hydrophobic filter, allows the passage of air and ensures closure to avoid any risk of contamination. The luer lock cap allows the extremity of the catheter to be closed after withdrawal of the internal needle.



Surshield® Versatus™ - Straight I.V. Catheters

Product code	Colour code	Diameter external catheter		Diameter internal catheter	Catheter length	Flow rate	Obturator
		Gauge	mm	mm	mm	ml/min	Product code
SR+DS2419SX	Yellow	24	0.74	0.55	19	29	SO+DM2419W
SR+DS2225SX	Blue	22	0.90	0.65	25	42	SO+DM2225W
SR+DS2032SX	Pink	20	1.00	0.75	32	59	SO+DM2032W
SR+DS1832SX	Green	18	1.30	0.95	32	103	SO+DM1832W
SR+DS1845SX	Green	18	1.30	0.95	45	96	SO+DM1845W
SR+DS1745SX	White	17	1.50	1.15	45	155	SO+DM1745W
SR+DS1645SX	Grey	16	1.75	1.35	45	225	SO+DM1645W
SR+DS1445SX	Orange	14	2.00	1.55	45	290	SO+DM1445W

Unit box : 50 pieces
Shipping carton : 200 pieces

Shelf life: 3 years



Sterilised by ethylene oxide

¹ Eucomed position paper - Preventing 'sharps' injury. P5 .3

² Ippolito and Puro : Device Specific Risk of NSI in Italian HCW. Jama 1994 ; 272 : 607-10.

Every day, as healthcare professionals, you are working to ensure life and to contribute to the well-being of your patients. Every day, Terumo Europe provides you with the clinical systems, means and solutions to keep life flowing.



TERUMO®

We keep life flowing

TERUMO CORPORATION
44-1, 2-chome Hatagaya Shibuya-ku Tokyo 151-0072 JAPAN
Tel. +81/3 3374.8111 Fax +81/3 3374.8196

TERUMO EUROPE N.V.
Researchpark Haasrode 1520 Interleuvenlaan 40 B-3001 Leuven BELGIUM
Tel. +32/16.38.12.11 Fax +32/16.40.02.49

TERUMO EUROPE N.V. BENELUX SALES DIVISION
Researchpark Haasrode 1515 Interleuvenlaan 40 B-3001 Leuven BELGIUM
Tel. +32/16.39.25.80 Fax +32/16.39.25.99 THE NETHERLANDS
Tel. 0800/022.03.96 Fax 0800/022.04.14

LABORATOIRES TERUMO FRANCE S.A.
1, Parc Ariane Bâtiment Uranus Rue Hélène Boucher 78284 Guyancourt Cedex FRANCE
Tel. +33/1 30.96.13.00 Fax +33/1 30.43.60.85

TERUMO DEUTSCHLAND GmbH
Hauptstrasse 87 65760 Eschborn GERMANY
Tel. +49/61.96.80.230 Fax +49/61.96.80.23.200

TERUMO DEUTSCHLAND GmbH
Zweigniederlassung Wettingen Bahnhofstrasse 88 CH-5430 Wettingen SWITZERLAND
Tel. +41/56.419.10.10 Fax +41/56.419.10.11

TERUMO EUROPE ESPAÑA SL
Avda. Juan Carlos I, Nº 13 - 7ª Planta Edificio Torre La Garena 28806 Alcalá de Henares (Madrid) SPAIN
Tel. general: +34/902.10.12.98 Tel. cust. serv.: +34/902.20.31.40 Fax general: +34/902.10.13.58 Fax cust. serv.: +34/902.10.17.57

TERUMO ITALIA SRL
Via Simone Martini 143/145 I-00142 Roma ITALY
Tel. +39/06.519.61.420 Fax +39/06.503.04.07

TERUMO SWEDEN AB
Nya Varvet, Byggnad 90 Sven Källfelts Gata 18 SE-426 71 Västra Frölunda SWEDEN
Tel. +46/31.748.58.80 Fax +46/31.748.58.90

TERUMO DENMARK Filial af TERUMO SWEDEN AB
Tel. +45/7020.93.80 Fax +45/7020.94.80

TERUMO UK Ltd.
Tamesis The Causeway Egham Surrey TW20 9AW UNITED KINGDOM
Tel. +44/1784.476.200 Fax +44/1784.476.201

TERUMO EUROPE N.V.
Moscow Representative office 4, Ivana Franko St 121108 Moscow RUSSIA
Tel. +7/495.234.3303 Fax. +7/495.234.3303

©: Registered Trademark Published by Terumo Europe N.V.

Addresses are subject to change, please consult our website.

www.terumo.com
www.terumo-europe.com