


Terusafe

Sharpscontainer



*Optimal disposal
with Terumo QUICK FIT*



Increased protection for healthcare workers

- *Safe disposal of sharp devices*
- *Various apertures for different options to disconnect*
- *Secure closure to prevent needle spill*
- *QUICK FIT - one hand release*
- *Penetration proof*
- *Latex-free*


Terusafe

Reduce the risk of percutaneous injuries

The TERUSAFE boxes are equipped with a variety of apertures enabling different options of sharp disposal. TERUSAFE boxes cannot be pierced or reached into. A secure closure prevents needle spill while transported.

All round boxes (0.5 L - 2.5 L) can be easily secured using an adhesive adapter. Therefore the boxes can swiftly be mounted near the patient's bedside.

- Aperture to remove **QUICK FIT** needles by pushing
- Aperture to remove **QUICK FIT** needles by twisting
- General aperture** for the disposal of Quick Safe-TE holder, single-use holder, scalpels, winged sampling needles, tubes, ...
- V-aperture** to remove needles by pulling
- U-aperture** to remove needles by pulling
- Aperture for insulin pen needles**
- Venoject® aperture** to unscrew Venoject® needles from their holder



Product References

Product Code	Volume	Shape	Pieces/Carton	Quick Fit Slot	Fixation
XX+Q6000	6.0 L	Round	50	No	No
XX+Q2500	2.5 L	Rectangular	45	Yes	No
XX+Q2510	2.5 L	Round	45	Yes	Yes
XX+Q2000	2.0 L	Round	65	No	Yes
XX+Q2010	2.0 L	Round	65	Yes	Yes
XX+Q1500	1.5 L	Round	80	No	Yes
XX+Q1510	1.5 L	Round	80	Yes	Yes
XX+Q0510LB	0.5 L	Round / Long	200	Yes	Yes*
XX+Q0510	0.5 L	Round / Long	200	Yes	Yes*
XX+Q0500	0.5 L	Round	120	No	Yes
XX+Q0300	0.3 L	Rectangular	125	No	No
XX+Q0310LB	0.25 L	Oval / Pocket-Shape	200	Yes	No
XX+Q0310	0.25 L	Oval / Pocket-Shape	200	Yes	No
XX+Q0001	Adhesive adapter (50 pcs/bag)				
XX+Q0001SB	SafeBOX table mount (10 pcs/unit) - for items indicated with *				



Terumo Corporation
44-1, 2-chome, Hatagaya
Shibuya-ku, Tokyo 151-0072,
Japan
Tel.: +81 3 3374 8111
Fax: +81 3 3374 8196

Terumo Europe N.V.
Interleuvenlaan 40
3001 Leuven, Belgium
Tel.: +32 16 38 12 11
Fax: +32 16 40 02 49

Terumo Europe N.V.
Benelux Sales Division
Belgium:
Tel.: +32 16 39 25 80
The Netherlands:
Tel.: 0800 0220396

Terumo Europe N.V.
Emerging Markets Division
Tel.: +32 16 38 13 08

Terumo Deutschland GmbH
Tel.: +49 6196 80 230

Terumo Deutschland GmbH
Zweigniederlassung Spreitenbach
Tel.: +41 56 419 10 10

Laboratoires Terumo France S.A.
Tel.: +33 130 96 13 00

Terumo Italia S.r.l.
Tel.: +39 0651 96 14 20

Terumo Europe España SL
Tel.: +34 902 10 12 98

Terumo UK Ltd
Tel.: +44 1784 476 200

Terumo Russia LLC
Tel.: + 7 495 988 4740

Terumo Sweden AB
Tel.: +46 3174 85 880

Terumo Denmark
Tel.: +45 7020 93 80

Terumo Corporation Dubai Branch
(Middle East and Africa)
Tel.: + 971 4 221 2220



LA 33GB-1106TE-III(12.13)E

Published by Terumo Europe N.V.