

### Tubing Options

Type	Size	Colour	SHORE Hardness	Xcoated	DEHP Free
PVC	1/2x3/32	Clear	70sh	x	Yes
PVC	1/2x3/32	Blue Stripe	70sh	x	Yes
PVC	1/2x3/32	Clear	63sh	x	Yes
PVC	1/2x3/32	Blue Stripe	67sh	x	Yes
PVC	1/2x3/32	Clear	65sh	x	Yes
PVC	1/2x3/32	Clear	65sh	Yes	No
PVC	1/2x3/32	Clear	70sh	Yes	No
PVC	3/8x1/16	Clear	70sh	x	Yes
PVC	3/8x3/32	Clear	70sh	x	Yes
PVC	3/8x3/32	Red Stripe	70sh	x	Yes
PVC	3/8x3/32	Blue Stripe	70sh	x	Yes
PVC	3/8x3/32	Clear	63sh	x	Yes
PVC	3/8x3/32	Clear	65sh	x	Yes
PVC	3/8x3/32	Red Stripe	67sh	x	Yes
PVC	3/8x3/32	Blue Stripe	67sh	x	Yes
PVC	3/8x3/32	Red Stripe	70sh	Yes	No
PVC	3/8x3/32	Blue Stripe	70sh	Yes	No
PVC	1/4x1/16	White Stripe	70sh	Yes	No
PVC	1/4x1/16	Yellow Stripe	70sh	Yes	No
PVC	1/4x1/16	Red Stripe	70sh	Yes	No
PVC	1/4x1/16	Blue Stripe	70sh	Yes	No
PVC	1/4x1/16	Green Stripe	70sh	Yes	No
PVC	1/4x1/16	Clear	65sh	Yes	No
PVC	1/4x1/16	Clear	70sh	Yes	No
PVC	1/4x3/32	Clear	65sh	Yes	No
PVC	1/4x1/16	Green Stripe	70sh	x	Yes
PVC	1/4x1/16	Blue Stripe	70sh	x	Yes
PVC	1/4x1/16	Red Stripe	70sh	x	Yes
PVC	1/4x1/16	Yellow Stripe	70sh	x	Yes
PVC	1/4x1/16	Clear	70sh	x	Yes
PVC	1/4x3/32	Clear	70sh	x	Yes
PVC	1/4x1/16	Red Stripe	67sh	x	Yes
PVC	1/4x1/16	Clear	63sh	x	Yes
PVC	1/4x1/16	White Stripe	67sh	x	Yes
PVC	1/4x1/16	Blue Stripe	67sh	x	Yes
PVC	1/4x1/16	Red Stripe	67sh	x	Yes
PVC	1/4x1/16	Yellow Stripe	67sh	x	Yes
PVC	1/4x1/16	Green Stripe	67sh	x	Yes
PVC	3/16x16	Clear	68sh	Yes	No
PVC	3/16x16	Clear	70sh	x	Yes
PVC	4.75x1.0	Clear	70sh	x	Yes
PVC	4.0x0.8	Clear	70sh	x	Yes
PVC	1/8x1/16	Clear	70sh	x	Yes
PVC	5/16x1/16	Clear	70sh	x	Yes
PVC	2.0x1.05	Clear	70sh	x	Yes
PVC	3/8x1/16	Clear	70sh	Yes	No
PVC	3/8x3/32	Clear	65sh	Yes	No
PVC	3/8x3/32	Clear	70sh	Yes	No
PVC	1/8x1/32	Clear	70sh	x	Yes
PVC	1.5x2.65	Clear	70sh	x	Yes
PVC	1/8x16	Clear	63sh	x	Yes
PVC	3/16x1/16	Clear	63sh	x	Yes
PVC	3.0x0.75	Clear	69sh	x	Yes
Silicone	1/2x3/32	Clear	N/A	x	Yes
Silicone	1/2x1/8	Clear	N/A	x	Yes
Silicone	3/8x1/16	Clear	N/A	x	No
Silicone	3/8x3/32	Clear	N/A	x	No
Silicone	3/8x1/8	Clear	N/A	x	No
Silicone	3/16x16	Clear	N/A	x	Yes
Silicone	1/4x1/16	Clear	N/A	x	Yes
Silicone	1/4x3/32	Clear	N/A	x	Yes
Silicone	5/16x1/16	Clear	N/A	x	Yes
Silicone	5/16x3/32	Clear	N/A	x	Yes
Silicone	1/8x1/16	Clear	N/A	x	Yes
Silicone	5/8x1/8	Clear	N/A	x	Yes
Tygon	1/2x3/32	Clear	N/A	x	Yes
Tygon	3/8x3/32	Clear	N/A	x	Yes
Tygon	3/8x1/16	Clear	N/A	x	Yes
Tygon	1/4x1/16	Clear	N/A	x	Yes

### Standard Packs

Material	Material Description
CX-FAC02	TUBING SET ANC DEVICE 1/4 SHUNT BYPASS
CX-FAC02X	1/4 SHUNT BYPASS LINE
CX-FAC03X	BYPASS LINE, 3/8 XCOAT
CX-FAC04X	1/2 SHUNT BYPASS LINE
CX-FAC05	GAS LINE
CX-FAC07	QUICK - PRIME LINE
CX-FAC109	1/2 RESERVOIR ADAPTER
CX-FAC113	3/8 RESERVOIR ADAPTER
CX-FAC116	3/16 PERFUSION ADAPTER
CX-FAC119	SUCKER WITH 10CM PVC 1/4X1/16
CX-FAC127	STANDARD SPARE GAS LINE
CX-FAC13	CONNECTOR, MALE-MALE LUER ADAPTER
CX-FAC16X	1/4 x 1/4 CONNECTOR (X COATED)
CX-FAC17	3/8 x 3/8 CONNECTOR
CX-FAC17X	3/8 x 3/8 CONNECTOR X COATED
CX-FAC20	1/4 X 3/8 CONNECTOR
CX-FAC20X	CONNECTOR 1/4x3/8 XCOATED
CX-FAC21	1/4 x 1/2 CONNECTOR
CX-FAC22	3/8 x 1/2 CONNECTOR
CX-FAC22X	3/8 x 1/2 CONNECTOR (X COATED)
CX-FAC25X	1/4 X 1/4 LL CONNECTOR X COATED
CX-FAC26	3/8X3/8 LL CONECTOR
CX-FAC28	3/16 x 1/4 LL CONNECTOR
CX-FAC29	1/4 x 3/8 LL CONNECTOR
CX-FAC29X	1/4 x 3/8 LL CONNECTOR (X COATED)
CX-FAC30	3/8 x 1/2 LL CONNECTOR
CX-FAC31	CONNECTOR, 3/16" Y-CONNECTOR
CX-FAC32	1/4x1/4x1/4 CONNECTOR
CX-FAC32X	1/4x1/4x1/4 CONNECTOR (X-COATED)
CX-FAC33	3/8X3/8X3/8 Y CONNECTOR
CX-FAC34	1/2X1/2X1/2 Y CONNECTOR
CX-FAC35	1/4 X 3/16 X 3/16 Y-CONNECTOR
CX-FAC36	3/8X1/4X1/4 Y-CONNECTOR
CX-FAC36X	3/8 x 1/4 x 1/4 Y CONNECTOR (X COATED)
CX-FAC37	3/8X3/8X1/4 Y-CONNECTOR
CX-FAC37X	3/8 X 3/8 X 1/4 CONNECTOR
CX-FAC38	1/2X1/2X3/8 CONNECTOR
CX-FAC38X	1/2 x 1/2 x 3/8 CONNECTOR X COATED
CX-FAC39	1/2X1/2X1/4 Y CONNECTOR
CX-FAC40	1/2x3/8x3/8 CONNECTOR

CX-FAC42	1/4X1/4X1/4 LL Y CONNECTOR
CX-FAC42X	CONNECTOR 1/4 Y LL X COATED
CX-FAC43	3/8X3/8X3/8 LL Y CONNECTOR
CX-FAC43X	3/8" Y-CONNECTOR LUER XCOAT
CX-FAC44	3/8X3/8 LL Y CONNECTOR
CX-FAC45	3/8X1/4X1/4 LL CONNECTOR
CX-FAC74	GAS FILTER
CX-FAC75	TERUMO STANDARD
CX-FAC97	CARDIOPLEGIA TABLE LINE
CX-FCP42	(CP50) 4:1 RATIO
CX-FHC05	PRECONNECTED HAEMO HC DP03HC W D/BAG
CX-FHC17	PRECONNECTED HAEMO HC DP09HC W D/BAG

### Major Components

Category	Description	Size
Blood Monitoring	Terumo CDI Cuvette	1/4x1/4
Blood Monitoring	Terumo CDI Cuvette	3/8x3/8
Blood Monitoring	Terumo CDI Cuvette	1/2x1/2
Cardioplegia	Terumo CP50	N/A
Cardioplegia	Stainless Steel Cooling Coil	N/A
Cardioplegia	Aluminium Cooling Coil	N/A
Hemoconcentrator	Medica DP03HC	Paediatric
Hemoconcentrator	Medica DP07HC	Adult
Hemoconcentrator	Medica DP09HC	Adult
Filter	Prebypass Filter, BT15X	N/A
Sucker	Sucker	N/A
Valve	Relief Check Valves, One Way Valves, Pillow Valves	N/A
Reservoir	Hardshell Reservoir, Cardiomy Reservoir	N/A
Oxygenator	FX05RW	Infant
Oxygenator	FX05RE	Infant
Oxygenator	FX15RE30C (Advance)	Paediatric
Oxygenator	FX15RW30C (Advance)	Paediatric
Oxygenator	FX15RE40C (Advance)	Small Adult
Oxygenator	FX15RW40C (Advance)	Small Adult
Oxygenator	FX25REC (Advance)	Adult
Oxygenator	FX25RWC (Advance)	Adult
Oxygenator	FX25RE (Standard)	Adult
Oxygenator	FX25RW (Standard)	Adult
Connectors	Various	Many Sizes
MLL's and FLL's	Various	Many Sizes
Stopcocks	Various	Many Sizes
Caps	Various	

### Packaging Type

Packaging Options	Tubing Pack Type	Oxygenator Preconnections Possible
Masterbox Packaging	Tubing Pack with Oxygenator	No
Tyvek Lidded Tray	Tubing Pack with Oxygenator	Yes
Waste Bin Packaging	Tubing Pack with Oxygenator	Some preconnections
Eco Packaging	Standalone Tubing Pack	N/A
Tyvek Lidded Tray	Standalone Tubing Pack	N/A
HC Box	Standalone Hemoconcentrator set	N/A
CP Box	Standalone Cardioplegia Set	N/A
CP Box	Accessory Lines	N/A