

Quality time for better care


"Quality time for better care" is the brand promise of Terumo Medical Care Solutions.



HOSPITAL CARE

AGANI™

Hypodermic Needle



Your needles of choice
for a wide range of applications

The safety and
comfort of your patients
is in your hands



Secure and easy handling increases safety and user comfort

The steeper point angle of the short bevel² shortens the length of the needle bevel, thereby reducing the risk of vascular perforation.

The outside diameter of the needle can be recognized by the ISO colour-coded hub.

- Easy identification of the required size
- Easy recognition of flashback through the semi-transparent hub

All needles come sterile, individually packed in colour-coded blister packaging

- Supplied in boxes of 100 needles
- Comprehensive information on the packaging



² Available on AN*1838S1 and AN*1938S1

Choose from a wide range from 18G up to 30G



Agani™ Hypodermic Needles

The range includes a broad variety of needle lengths and diameters



18G Cannula (ø 1.2 mm)
Length 38 mm
Colour Pink



19G Cannula (ø 1.1 mm)
Length 25, 38, 50 mm
Colour Cream



20G Cannula (ø 0.9 mm)
Length 25, 38 mm
Colour Yellow



21G Cannula (ø 0.8 mm)
Length 16, 25, 38, 50 mm
Colour Green



22G Cannula (ø 0.7 mm)
Length 32, 38 mm
Colour Black



23G Cannula (ø 0.6 mm)
Length 16, 25, 32 mm
Colour Blue



24G Cannula (ø 0.55 mm)
Length 25 mm
Colour Purple



25G Cannula (ø 0.5 mm)
Length 16, 25 mm
Colour Orange



26G Cannula (ø 0.45 mm)
Length 13, 23 mm
Colour Brown



27G Cannula (ø 0.4 mm)
Length 16, 19 mm
Colour Grey



30G Cannula (ø 0.3 mm)
Length 13 mm
Colour Yellow

CE 0123 Does not contain components made of PVC or natural rubber latex.

For further ordering information, please refer to the Agani™ Safety Hypodermic needle product information card.

Smooth and easy puncture increase patient comfort

Make the administration of medication as painless as possible with the special geometry of the sharp lancet point and silicone-coated cannula that allows a smooth penetration and prevents any knocking effect¹.



Precisely shaped needle point with silicone coating.


Agani™ Hypodermic needles have a sharp lancet tip and an outside layer of silicone lubricant.

The sharp point

- Facilitates the penetration of the skin
- Helps to minimize the trauma

The silicone coating

- Provides good gliding properties
- Allows a smooth insertion



1 Data on file

Life is easier with needles from Terumo

Agani hypodermic needles provide you and your patients with the best possible levels of care and comfort.

- Sharp lancet point – easy puncture
- Silicone coating – smooth needle movement
- Thin walled needles³ – increased flow rate
- Short bevel options – reduced risk of vascular perforation
- Large portfolio – suitable for a wide range of applications



Injection & Access

**Contributing to Society
through Healthcare**



TERUMO EUROPE N.V.

Terumo Corporation Terumo Europe N.V.
+81 3 6742 8500 +32 16 38 12 11

www.terumo-europe.com

Your Terumo partner

All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties. Published by Terumo Europe N.V.

³ Available on specific item codes