

Quality time for better care

"Quality time for better care" is the brand promise of Terumo Medical Care Solutions.



HOSPITAL CARE

Administration Sets



- EASY HANDLING
- ENHANCED SAFETY
- Terufusion™
- OPTIMIZED WORKFLOW

Supporting your infusion therapy

A clear view on infusion with a comprehensive portfolio of solutions



Easy handling

- Easy setting following the instruction screen of the pump
- Sharp 2-way spike to penetrate bags and bottles easily
- Adjustable roller clamp with notch to mount tubing during priming

Increased patient safety

- AFF clip for prevention of accidental bolus infusion⁺
- Does not contain DEHP
- Y-site with NLCD* for prevention of needle stick injuries⁺

Highly accurate infusion

- Reliable high accuracy combination of Terufusion™ pump & sets[#]
- Durable tubes for continuous accurate infusion[§]

* Needleless connecting device

⁺ Available on specific sets

[#] Flow rate accuracy: within $\pm 5\%$

[§] Always consult the product label and IFU for a complete overview of warnings, cautions and/or precautions prior to actual use.



General

Product code	TI*PU300WY	TI*PJ300WY01	PT+PU400WY	PT+PU402WY
AFF clip	Yes	Yes	Yes	Yes
Filter (µm)	15	15	15	15
Pump tube (mm)	2000	2000	2000	2000
Injector port	None	Needle port	Needle port	T-Port ⁷ x 1
Drip rate (drops/ml)	20	20	20	20
Material - tubing	PVC ¹	Polybutadiene ⁴	PVC ¹	PVC ¹
Packaging (unit/carton)	20/200	20/200	25/250	25/250
Priming volume (ml) ²	15.7	16.7	18.1	18.1
Colour ³	Transparent	Transparent	Transparent	Transparent
Pressure resistance (bar)	1.5	1.5	1.5	1.5
Total length of sets	2150	2440	2380	2380
Product Illustration				

General (amber)

TI*PA300WY01E	TI*PA300WN01E
Yes	No
15	15
2000	2000
Needle port	Needle port
20	20
PVC ¹	PVC ¹
20/200	20/200
18.2	18.2
Amber	Amber
1.5	1.5
2640	2640

Transfusion

TB*PU300SYE
Yes
200
2200
None
20
PVC ¹
20/200
17.2
Transparent
1.5
2420

Specific

SA*QTF301NN	SA*QNT310MN	SA*QNT310YN	SA*QNT320MN
Yes	No	Yes	No
0.2	40	40	40
1220	1220	1220	1220
T-Port ⁶ x 1	3WSC ⁵ x 1	3WSC ⁵ x 1	3WSC ⁵ x 2
20	20	20	20
Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴
20/160	20/160	20/160	20/160
15.4	13.5	11.8	13.9
Transparent	Transparent	Transparent	Transparent
1.5	1.5	1.5	1.5
2450	2400	1900	2450

Specific

Product code	SA*QNT320YN	SA*QNT330MN	SA*QNT330YN	SA*QTF301UN
AFF clip	Yes	No	Yes	No
Filter (µm)	40	40	40	0.2
Pump tube (mm)	1220	1220	1220	1220
Injector port	3WSC ⁵ x 2	3WSC ⁵ x 3	3WSC ⁵ x 3	T-Port ⁶ x 1
Drip rate (drops/ml)	20	20	20	20
Material - tubing	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴
Packaging (unit/carton)	20/160	20/160	20/160	20/160
Priming volume (ml) ²	12.2	14.3	14.3	15.4
Colour ³	Transparent	Transparent	Transparent	Transparent
Pressure resistance (bar)	1.5	1.5	1.5	1.5
Total length of sets	1950	2500	2500	2450
Product Illustration				

SA*QTF302NN	SA*QTF302UN	SA*QTF303NN	SA*QTF303UN	SA*QTT301LN
Yes	No	Yes	No	No
0.2	0.2	0.2	0.2	15
1220	1220	1220	1220	1220
T-Port ⁶ x 2	T-Port ⁶ x 2	T-Port ⁶ x 3	T-Port ⁶ x 3	T-Port ⁶ x 1
20	20	20	20	20
Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴
20/160	20/160	20/160	20/160	20/160
15.1	15.1	15.8	15.8	11.7
Transparent	Transparent	Transparent	Transparent	Transparent
1.5	1.5	1.5	1.5	1.5
2350	2350	2500	2500	1900

SA*QTT301MN	SA*QTT301RNA	SA*QTT301YN	SA*QTT302MN	SA*QTT302YN
No	Yes	Yes	No	Yes
40	40	40	40	40
1220	2200	1220	1220	1220
T-Port ⁶ x 1	T-Port ⁷ x 1	T-Port ⁶ x 1	T-Port ⁶ x 2	T-Port ⁶ x 2
20	20	20	20	20
Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴	Polybutadiene ⁴
20/160	20/160	20/160	20/160	20/160
11.7	18.4	11.7	12.2	12.2
Transparent	Transparent	Transparent	Transparent	Transparent
1.5	2.0	1.5	1.5	1.5
1900	2600	1900	2050	2050

Terufusion - Full focus infusion system



Terumo offers a comprehensive portfolio of solutions for the complete infusion line-up. User-friendly syringe and volumetric pumps and racks, as well as a wide range of highly accurate, easy-to-use disposables, such as an extended range of infusion sets and syringes.