

Quality time for better care

"Quality time for better care" is the brand promise of Terumo Medical Care Solutions.



HOSPITAL CARE

# Administration Sets



Supporting your  
infusion therapy

# A clear view on infusion with a comprehensive portfolio of solutions



## Easy handling

- Easy setting following the instruction screen of the pump
- Sharp 2-way spike to penetrate bags and bottles easily
- Adjustable roller clamp with notch to mount tubing during priming

## Increased patient safety

- AFF clip for prevention of accidental bolus infusion\*
- Does not contain DEHP
- Y-site with NLCD\* for prevention of needle stick injuries\*

## Highly accurate infusion

- Reliable high accuracy combination of Terufusion™ pump & sets#
- Durable tubes for continuous accurate infusion§

\* Needleless connecting device

\* Available on specific sets

# Flow rate accuracy: within  $\pm 5\%$

§ Always consult the product label and IFU for a complete overview of warnings, cautions and/or precautions prior to actual use.





General

Product code	Ti*PU300WY	Ti*PJ300WY01	PT+PU400WY
AFF clip	Yes	Yes	Yes
Filter (µm)	15	15	15
Pump tube (mm)	2000	2000	2000
Injector port	None	Needle port	Needle port
Drip rate (drops/ml)	20	20	20
Material - tubing	PVC <sup>1</sup>	Polybutadiene <sup>4</sup>	PVC <sup>1</sup>
Packaging (unit/carton)	20/200	20/200	25/250
Priming volume (ml) <sup>2</sup>	15.7	16.7	18.1
Colour <sup>3</sup>	Transparent	Transparent	Transparent
Pressure resistance (bar)	1.5	1.5	1.5
Total length of sets	2150	2440	2380
Product Illustration			

General (amber)

Ti*PA300WY01E	Ti*PA300WN01E
Yes	No
15	15
2000	2000
Needle port	Needle port
20	20
PVC <sup>1</sup>	PVC <sup>1</sup>
20/200	20/200
18.2	18.2
Amber	Amber
1.5	1.5
2640	2640

Transfusion

TB*PU300SYE
Yes
200
2200
None
20
PVC <sup>1</sup>
20/200
17.2
Transparent
1.5
2420

Specific

SA*QTF301NN	SA*QNT310MN	SA*QNT310YN	SA*QNT320MN
Yes	No	Yes	No
0.2	40	40	40
1220	1220	1220	1220
T-Port <sup>5</sup> x 1	3WSC <sup>6</sup> x 1	3WSC <sup>6</sup> x 1	3WSC <sup>6</sup> x 2
20	20	20	20
Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>
20/160	20/160	20/160	20/160
15.4	13.5	11.8	13.9
Transparent	Transparent	Transparent	Transparent
1.5	1.5	1.5	1.5
2450	2400	1900	2450

Specific

Product code	SA*QNT320YN	SA*QNT330MN	SA*QNT330YN
AFF clip	Yes	No	Yes
Filter (µm)	40	40	40
Pump tube (mm)	1220	1220	1220
Injector port	3WSC <sup>6</sup> x 2	3WSC <sup>6</sup> x 3	3WSC <sup>6</sup> x 3
Drip rate (drops/ml)	20	20	20
Material - tubing	Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>
Packaging (unit/carton)	20/160	20/160	20/160
Priming volume (ml) <sup>2</sup>	12.2	14.3	14.3
Colour <sup>3</sup>	Transparent	Transparent	Transparent
Pressure resistance (bar)	1.5	1.5	1.5
Total length of sets	1950	2500	2500
Product Illustration			

SA*QTF301UN	SA*QTF302NN	SA*QTF302UN	SA*QTF303NN	SA*QTF303UN
No	Yes	No	Yes	No
0.2	0.2	0.2	0.2	0.2
1220	1220	1220	1220	1220
T-Port <sup>5</sup> x 1	T-Port <sup>5</sup> x 2	T-Port <sup>5</sup> x 2	T-Port <sup>5</sup> x 3	T-Port <sup>5</sup> x 3
20	20	20	20	20
Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>
20/160	20/160	20/160	20/160	20/160
15.4	15.1	15.1	15.8	15.8
Transparent	Transparent	Transparent	Transparent	Transparent
1.5	1.5	1.5	1.5	1.5
2450	2350	2350	2500	2500

SA*QTT301LN	SA*QTT301MN	SA*QTT301YN	SA*QTT302MN	SA*QTT302YN
No	No	Yes	No	Yes
15	40	40	40	40
1220	1220	1220	1220	1220
T-Port <sup>5</sup> x 1	T-Port <sup>5</sup> x 1	T-Port <sup>5</sup> x 1	T-Port <sup>5</sup> x 2	T-Port <sup>5</sup> x 2
20	20	20	20	20
Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>	Polybutadiene <sup>4</sup>
20/160	20/160	20/160	20/160	20/160
11.7	11.7	11.7	12.2	12.2
Transparent	Transparent	Transparent	Transparent	Transparent
1.5	1.5	1.5	1.5	1.5
1900	1900	1900	2050	2050

1 Does not contain DEHP - 2 Total volume of the infusion line (drip chamber empty) - 3 Tubing & Drip chamber - 4 Does not contain PVC - 5 SURPLUG AD T-Port Needleless connecting device - 6 SURPLUG AD 3-way stopcock

# Terufusion - Full focus infusion system



Terumo offers a comprehensive portfolio of solutions for the complete infusion line-up. User-friendly syringe and volumetric pumps and racks, as well as a wide range of highly accurate, easy-to-use disposables, such as an extended range of infusion sets and syringes.